# Reusable EKG / ECG Leadwire

English

**Indications:** EKG/ECG leadwires are used to convey electrocardiographic signals of a patient from ECG electrodes placed on the patient's skin to an electrocardiograph or patient monitor via EKG/ECG trunk cable. ECG cables are intended to be used by qualified medical personnel on infant, pediatric, and adult patients in hospitals and hospital-type facilities.

**Contraindications:** No known contraindications to the use of the product.

**Adverse Reactions:** No known adverse reactions to the use of the product.

**Clinical Benefits:** Enables ECG monitoring or recording with diagnostic EKG equipment or ECG monitoring equipment.

**Product Description:** FMT EKG/ECG leadwires have different trunk cable-side connectors and connector wiring configurations for different electrocardiograph or patient monitor manufacturers. Therefore, the same appearance does not guarantee full compatibility with the device. Always consider the information on the product's label and/or primary packaging. There are 4.0 mm banana, 3.2 mm straight pin, snap, or grabber-type patient terminations available.

**E2Y-XZ**, 1.6mm DIN Type leadwire  
**E3Y-XZ**, Siemens/Drager Type leadwire  
**E4Y-5Z**, GE Type leadwire  
**E5Y-XZ**, Philips Type leadwire  
**E6Y-XZ**, Spacelabs Type leadwire  
**E9Y-5Z**, Datascope Type lead leadwire  
**E4Y-3ZT**, GE Type leadwire  
**E4Y-5ZT**, GE Type leadwire  
**E7Y-3Z**, Nihon Kohden Type leadwire  
**E7Y-6Z**, Nihon Kohden Type leadwire  
**E8Y-3Z**, Philips Type leadwire  
**E8Y-5Z**, Philips Type leadwire

**MTX2Z**, leadwire for Drager MultiMed Plus and MultiMed Plus OR Multiparameter cable, (only with grabber).

(X shows number of the leadwires, 3, 5 or 6. Y shows patient termination type 1=Snap, 2=Grabber. Z shows color code I=IEC, A=AHA/AAMI, -OR suffix indicates ESU protection)

**E101YZ/10**, leadwire for Marquette/Hellige EKG trunk cable.

**E102YZ/10**, leadwire for Philips M3703C EKG module.

**E103YZ/10**, leadwire for Philips Trim EKG module.

**E104YZ/10**, leadwire for GE AM4 / AM5 EKG module.

**E105YZ/10**, leadwire for GE AM6 / CAM14 EKG module.

**E112YZ/10**, leadwire for Norav 1200W EKG module.

**E117YZ/10**, leadwire for Quinton PM EKG module.

**E106YZ/5 - E107YZ/5**, leadwire for GE Multi-Link EKG trunk cable.

**E108YZ/5 - E109YZ/5**, leadwire for GE Multi-Link MAC 400/1600 EKG trunk cable.

**E110YZ/5 - E111YZ/5**, leadwire for Philips M1663A, M1949A EKG trunk cable.

**E113YZ/4 - E114YZ/6**, leadwire for Medtronic/Physio Control Lifepack 12/20 EKG trunk cable.

**E118YZ/10**, leadwire for Mortara AM12 and WAM EKG module.

**E118YZ/5A F - E118YZ/5B**, leadwire for Mortara AM12 and WAM EKG module.

**E119GZ/10**, leadwire for Mortara AM12 and WAM Stress EKG module (only with grabber).

**E119GY/5A - E119GY/5B**, leadwire for Mortara AM12 and WAM Stress EKG Module (only with grabber).

**E120YZ/10**, leadwire for Welch Allyn CP100/200 module.

**E121YZ/10**, leadwire for Philips PageWriter TC 30/50/70 module.

**E121YZ/3A - E121YZ/3B**, leadwire for Philips PageWriter TC 30/50/70 module.

**E122GZ/6**, grouped leadwire for Schiller Argus Pro Lifecare2 EKG trunk cable (only with grabber).

(X shows number of the leadwires, 3, 5 or 6. Y shows patient termination type B=4.0mm Banana, P=3.2mm Straight Pin, S=Snap, G=Grabber. Z shows color code I=IEC, A=AHA/AAMI. -OR suffix indicates ESU protection)

## Directions for Use:

- 1- Before use, check that the product is intact and clean.
- 2- Connect the EKG/ECG leadwire to the compatible electrocardiograph or patient monitor by attaching the leadwire plugs to the yoke of the compatible trunk cable and check for correct operation
- 3- Before applying the ECG electrodes ensure the patient skin is clean, dry, and properly prepared. Always apply the electrodes to the intact skin of a patient.
- 4- Attach leadwire patient terminations to the electrodes. Make sure that the cables do not pull the electrodes.
- 5- Apply the electrodes firmly on prepared sites.
- 6- Ensure the leads conform to the body contours and no strain is placed on the electrodes.
- 7- Verify that the ECG signal appears on the electrocardiograph or patient monitor screen. For detailed information on ECG recording or monitoring, please refer to the electrocardiograph or patient monitor instructions for use.

## Warnings:

- 1- All EKG/ECG leadwires are designed for use with the specific yoke of the ECG trunk cable of a compatible patient monitor or electrocardiograph. Verify the compatibility of the patient monitor or electrocardiograph, trunk cable, and ECG leadwires and all functions properly before use, otherwise degraded performance may result.
- 2- Before each use, visually inspect the ECG leadwires. The ECG leadwires should be free of visible defects and damage. If there are breaks, scratches, or cracks on the cable jacket or plastics insulation of the connector or patient terminations, or there are visibly bent, damaged pins of the connector or patient terminal discontinue using the cable. Never use damaged leadwires or one with exposed electrical circuitry. Even if there is no exposed electrical circuit, the scratches or cracks in the cable jacket or plastic insulation may cause the accumulation of pathogens and reduce the effectiveness of cleaning and disinfection of the product. Dispose of damaged leadwires according to hospital disposal practices, local environmental and waste disposal laws, and regulations.
- 3- Do not use the ECG leadwires during Magnetic Resonance Imaging (MRI) scanning or in an MRI environment, as this may result in physical harm.
- 4- Keep the ECG cable connectors as far as possible from electromagnetic sources such as an electromedical device's power transformer. Also, avoid placing the EKG cable alongside the device's power cord. Strong electromagnetic fields may adversely interfere with the function of an ECG cable.
- 5- As with all medical equipment, carefully route cables to reduce the possibility of patient entanglement or strangulation.
- 6- To prevent shock hazards do not connect ECG leadwires to any power outlet.
- 7- To avoid a shock hazard and interface with nearby electrical equipment, keep electrodes and patient cables away from grounded metal and other electrical equipment.
- 8- When connecting and disconnecting the ECG leadwire and the patient terminations pull on the connector or plug itself, not on the cable.
- 9- Do not modify the ECG leadwire in any way. Modifications may affect patient safety, performance, and accuracy and void the warranty.
- 10- Mishandling of the ECG leadwires could result in damage to internal wires or loss of electrical isolation. Do not tie the ECG cable and leadwires in knots or bunch tightly. Do not wrap ECG leadwires around equipment to avoid damaging internal wires. Do not leave the ECG cable dangling from the instrument. Keep ECG cable leadwires away from moving parts that might pinch or crush them. When not in use, loosely coil the ECG leadwire to a bend diameter of at least 10 cm and store it at room temperature.
- 11- Dirty, compromised or unprepared skin may cause incorrect ECG readings. Using previously opened or out-of-date electrodes may degrade the ECG signal quality. ECG electrodes should be removed from their sealed package immediately prior to use.
- 12- Mating connectors with force without proper alignment may damage the connectors and disrupt the electrical continuity.
- 13- For further information and warnings, read the instruction manual of the ECG device before using the ECG leadwires.

### Manual Cleaning and Disinfection:

**Manual Cleaning:** The EKG/ECG leadwires must be cleaned and disinfected immediately after each use. Remove heavy debris from the EKG/ECG leadwire by wiping with an appropriate, moist, lint-free cloth, sponge, or gauze pad. Prepare the detergent solution with reverse osmosis or distilled water according to the manufacturer's instructions. Saturate a synthetic felt or sponge pad with the cleaning solution. Wipe all surfaces of the EKG/ECG leadwire until all visible contamination is removed. Wash the pad in the cleaning solution and repeat wiping the cable. Use a plastic cleaning brush for adhered particles or dried secretions. Wipe or brush all surfaces at least six strokes back and forth. At least allow a contact time of 2 minutes. If visible contamination is still present on the cable at the end of the cleaning process, repeat the cleaning process. Moist a clean cloth, sponge, or gauze pad with tap water or reverse osmosis and thoroughly wipe the EKG/ECG leadwire to remove any debris and all traces of the cleaning solution. Dry manually by using a lint-free cloth or absorbent paper and, in particular, for drying cavities, sterile compressed air.

**Manual Disinfection:** Prepare the disinfectant solution according to the manufacturer's instructions. Only disinfect clean products that have dried slightly. Saturate a synthetic felt or sponge pad with the disinfectant solution. Wipe all surfaces of the EKG/ECG leadwire. Wash the pad in the disinfectant solution and repeat wiping the cable. Wipe or brush all surfaces at least six strokes back and forth. At least allow a contact time of 6 minutes. Moist a clean cloth, sponge, or gauze pad with reverse osmosis or distilled water and thoroughly wipe the EKG/ECG leadwire to remove any residual debris and all traces of disinfectant solution. Dry manually by using a lint-free cloth or absorbent paper and, in particular, for drying cavities, sterile compressed air.

EKG/ECG leadwire, especially metallic parts of the plugs and connector should never be immersed in cleaning, disinfecting, or rinsing solutions, or placed in an ultrasonic cleaning device. Soaking will reduce the life of EKG/ECG leadwire by oxidation of metal parts and hardening of the cable jacket. Avoid using excessive force when cleaning or disinfecting by hand which may damage the internal wires of the cable and reduces the product life. Please do not clean with hard or frictional materials, abrasive cleaners, or solvents. Cleaning agents containing phenol or chloride are not suitable. The EKG/ECG leadwire should be completely dry before use, do not use wet EKG/ECG leadwire. These instructions were validated using the below detergents and disinfectants.

Cleaning agent (manual cleaning): Enzo® / Cidezyme®, enzymatic detergent, Johnson&Johnson (add 8ml enzymatic detergent per 1 liter distilled water, for cable with dried-on organic matter use 16 ml enzymatic detergent per 1 liter distilled water/ warm water)

Disinfectant (manual disinfection): Cidex OPA®, 0.55% ortho-phthalaldehyde solution, Johnson&Johnson (at a minimum temperature of 20°C/68°F). Reverse osmosis/ distilled water.

**Caution:** During the cleaning and disinfection of EKG/ECG leadwires, disposable gloves, protective eyewear, and a filtration mask should be used to reduce the risk of transmission of infectious agents by splashing or creating aerosols.

**Warning:** FMT EKG/ECG leadwires are sold NON-STERILE. Clean and disinfect the cables before the first and every use.

**Warning:** Disconnect the cable from the EKG/ECG device or the trunk cable's yoke before cleaning or disinfecting it.

**Warning:** EKG/ECG leadwires should be cleaned of overburden before disinfection to improve the effectiveness.

**Warning:** Do not sterilize ECG leadwires by autoclave, radiation, or steam.

**Warning:** Avoid contact with strong, aromatic, chlorinated, ketone, ether, or ester solvents. Do not use bleach on electrical contacts or connectors.

### Environmental Requirements:

Operating Temperature : 5°C - 45°C (41°F - 113°F)

Storage Temperature : 0°C - 50°C (32°F - 122°F)

Relative Humidity : 20 - 80% (not condensing)

During the storage, the product should be protected from sunlight. It is recommended to store the product in its original package until the first use. Used cables should be loosely coiled in 10 to 15 cm diameter (4 to 6 inches) and must be stored individually in a protective container or packaging.

**Compatibility:** To ensure compatibility, the EKG/ECG leadwire should only be used with the equipment for which it has been designed and is specified on the label on the product's primary packaging. For product models and compatible devices, please refer to the catalogue at [metkomedical.com/PDF/01DiagnosticEKGCablesandAccessories.pdf](http://metkomedical.com/PDF/01DiagnosticEKGCablesandAccessories.pdf) for diagnostic EKG leadwires and [metkomedical.com/PDF/02ECGCablesforPatientMonitoring.pdf](http://metkomedical.com/PDF/02ECGCablesforPatientMonitoring.pdf) for monitoring ECG leadwires.

### Safety:

Degree of protection from electric shocks: type CF

Degree of protection against the ingress of water: Keep dry

The FMT EKG/ECG leadwires are sold NON-STERILE.

The FMT EKG/ECG leadwires are not made with natural rubber latex and PVC.

The FMT EKG/ECG leadwires have not been evaluated for safety in the MR environment. Scanning a patient who has this device may result in patient injury.

For an explanation of symbols, refer to the separate leaflet "Description of Symbols" located within the product package.

**Reporting of Serious Incidents:** Any serious incident related to the use of this product should be reported to both the manufacturer and the health authority/competent authority where the product is installed. Either contact your local representative or report to: [metko@metkomedical.com](mailto:metko@metkomedical.com)

**Caution:** Federal Law (U.S.A.) restricts this device to sale by or on the order of a physician.

**Warranty:** FMT EKG/ECG leadwires are under six (6) months warranty against material and workmanship defects from the date of original purchase. In the warranty period, METKO will be responsible for repairing or changing the cable free of charge if the defect is proven. This warranty does not extend to any product that has been subject to misuse, wrong disinfection, neglect, or accident, or that has been damaged by causes external to the product or that has been used in violation of the operating instructions supplied with the product. The product life is 3 years from the production date.

FMT® is a registered trademark of Metko A.Ş.

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