Reusable High Frequency Bipolar and Monopolar Cable

English

Indications: Reusable high-frequency (HF) bipolar and monopolar cables are intended for transferring high-frequency electrical power/energy from the compatible electrosurgical generator (high-frequency unit) to the bipolar/monopolar electrosurgical instruments without changing. High-Frequency bipolar and monopolar cables are intended to be used with appropriate electrosurgical equipment by trained and qualified medical personnel in hospitals and hospital-type facilities.

Contraindications: While not an absolute contraindication to the use of electrosurgery, implantable electrical devices (IEDs) such as cardiac and gastric pacemakers, implantable cardiac defibrillators (ICDs), and cochlear implants, and deep brain, nerve, spinal cord, or bone stimulators create unique risks. HF current during electrosurgery can interfere with the functioning of these devices and devices can be damaged by HF current. Patients with these devices require a thorough pre-procedure evaluation and may require intraoperative monitoring and post-procedure device assessment.

Adverse Reactions / Events: Adverse events have been reported when using electrosurgical devices, such as tissue damage at the wrong site caused by unintentional activation, fires involving surgical drapes and other flammable materials, burns caused by alternating current paths where the patient or user is in contact with exposed metal, explosions caused by an electrosurgical spark in the flammable gas environment, organ perforation and sudden massive bleeding.

Clinical Benefits: Enables conducting of high-frequency current during electrosurgery.

Product Description: FMT HF cables are manufactured from Thermoplastic Elastomer (TPE) and silicone. The connectivity depends on the specific types of connectors, both on the generator side and the instrument side. Use only compatible electrosurgical generators and instruments.

ESU-BP/XX (3m), ESU-BP/XXL (5m), ESU-BP/XWE (4m), High-Frequency Bipolar Cable
ESU-MP/XX (3m), ESU-MP/XXL (5m), High-Frequency Monopolar Cable (XX: Code for Compatibility)
Maximum accessory rated voltage (U_{max}): High-frequency monopolar cables 4.0 kVp, High-frequency bipolar cables 0.8 kVp, High-frequency monopolar resectoscope cables 2.0 kVp

When combining with other electrosurgical devices, ensure that the output parameters of the electrosurgical generator do not exceed the cable's rated voltage. In any combination with another electrosurgical accessory, the maximum rated voltage of the combination corresponds to the lowest rated voltage of the accessories used.

Directions for Use:

Prior to use, HF cables and their connectors should be inspected for visual defects. Loose or damaged connector mating parts (broken or missing components), corrosion of metal contacts, damage in the insulation of cable (i.e. cuts, cracks, brittle or solid, permanent bends, burnt or melted), and wires or metal parts that have come out of the cable/connector insulation, loose or improperly attached connectors, and not completely dried and ventilated cables should not be used.

- 1- HF cables are used to transmit the HF current required to cut or coagulate tissues to high-frequency electrodes, forceps, or instruments. Therefore, a suitable HF electrode or instrument should be connected to the socket side of the HF cable. Connection to the generator is made with the plug of the HF cable.
- 2- Check whether the plug is suitable for the correct socket on the generator. Careless push or pushing with force may result in damage to the plug and/or the socket on the generator.
- 3- The instrument connector(s) should be plugged into the socket(s) of the HF cable covering all the metal contacts of the instrument connector.
- 4- When connecting the cable it is important to check whether good and continuous contact is obtained. This is obtained by pushing the HF cable plug into the mechanical borders of the socket on the generator covering all the metal contacts.
- 5- Before disconnecting the cables always switch the HF unit off. When connecting and disconnecting cables, hold only the plastic part of the cable connector. Otherwise, it may result in a burn or electric shock in patients or users. Do not pull on the cable. In preparation, transport, sterilization, and storage cables should not be bent, twisted tightly, or folded. Such movements reduce the life of the product. Do not load heavy or sharp objects on the cables. This may damage the cable insulation. Damage to the insulation may result in burns with the unwanted electric current passes.
- 6- Activation of the HF current is done by the foot or hand switch or if available by autostart mode of the unit. Directives in the user manual of the HF unit should be followed.

Special Information: The socket on the unit for bipolar handles with 2-line (bipolar) cables can be mixed with other bipolar coagulation instruments. For example, the bipolar coagulation connection cable in ERBE fits the monopolar cut and coagulation socket of high-frequency units of AESCULAP, BERCHTOLD, INTEGRA, and MARTIN. For this reason before use make sure that the product is correct and suitable for the HF unit.

Warnings:

- 1- The user of HF bipolar/monopolar cables should have the necessary knowledge, expertise, and education to use the products. Read all instruction manuals for HF cable, generator, and connected instruments thoroughly. An insufficient understanding of the dangers, warnings, cautions, and information in these manuals can result in death, serious injury, or equipment damage.
- 2- Before first and each use, an electrical continuity test must be performed and the insulation of the HF cable must be inspected for damage. If there is a break in the cable wire or the cable becomes otherwise electrically discontinuous, arcing may occur in the patient-return or active circuit and may burn the patient or user or create a fire. Visually inspect the entire HF cable. Visual inspection alone may not be sufficient to ensure that the insulation is intact. Inspection may be done visually under magnification or with a high-voltage insulation testing device. Do not use HF cables with brittle or defective insulation or exhibiting pressure points.
- 3- Reusable HF cables should be inspected, cleaned, sterilized, and dried/aerated before first and each subsequent use. Improper or incomplete reprocessing can cause infection to the patient or medical personnel. The user is recommended to implement a program to track the number of uses of the cables. The serial number on the autoclave record tag attached to the HF cable allows the user to keep track of the HF cable and record reprocessing cycles or the service life of each HF cable.
- 4- Do not perform electrosurgery in the presence of flammable anesthetics or oxidizing gases, other flammable gases, or in close proximity to volatile solvents, fluids, or objects, or in oxygen-enriched environments, as explosions and fire could result. Do not place electrosurgical instruments near or in contact with flammable materials (such as gauze or surgical drapes). Instruments that are activated or hot from use may cause a fire.
- 5- High-frequency accessories connected to the electrosurgery unit but not temporarily used should not be touched to the patient and held in a different position. The active electrode should be placed in a clean, dry, well-insulated safety holster when not in use to prevent inadvertent activation or contamination. Inadvertent contact with the patient may result in burns.
- 6- The HF cables should not touch directly to the patient's skin or other cables should be kept at a minimum length and should be plugged in without any friction to any place. Do not lay high-frequency cables in loops. HF cables should not be positioned parallel to the camera cables to avoid interference.
- 7- During the use of high-frequency equipment in patients with cardiac pacemakers, implanted heart defibrillators, or pace electrodes, there is a danger of ruining the function or the device itself. The pacemaker or electrodes should never be touched with high-frequency electrodes. Before using advice from cardiology departments must always be taken.
- 8- HF cables should only be connected to the high-frequency generator when the unit is switched off. They should not be connected while preparing the electrosurgical instrument. Failure to do so may result in injury or electrical shock to the patient or user.
- 9- Avoid unnecessary or prolonged activation. This may cause injuries due to capacitive coupling with other surgical equipment. Activate only when touching or immediately next to the intended point of contact. Keep the voltage/power settings at the lowest possible level to achieve the desired effect. Keep the active electrodes clean. When the effectiveness of the instrument is reduced, do not increase the power of the generator. Instead, clean the active electrode contacts.
- 10- If using irrigation near a HF cable, care should be taken to ensure that the liquid does not enter the connection between the electrosurgical fitting and the cable. There is a high risk of electrical shock if liquid enters the fitting.
- 11- Aspirate fluid from the area before activating the HF instrument. Conductive fluids (e.g., blood or saline) in direct contact with or in close proximity to an active electrode may carry electrical current or heat away from target tissues, which may cause unintended burns to the patient. Flammable solutions pooling under the patient or in body depressions may cause burns to the patient.
- 12- Potentially hazardous conditions may exist when accessories are not appropriate for the type of generator output used. If unusual operating characteristics develop or are observed during procedures, the use of the product should be discontinued.
- 13- Connect bipolar cables to the bipolar receptacle and monopolar cables to the monopolar receptacle only to avoid accidental accessory activation or hazards.
- 14- Make sure that the patient is electrically isolated against electrical conductive elements and does not come into direct contact with other metal parts or electric conductive instruments which are grounded or which have an appreciable capacitance to ground. Metal devices can provide an alternate pathway for the electrical current and be responsible for patient burns. Avoid skin-to-skin contact of the patient.
- 15- HF cables should not be wrapped around metal objects or bundled together. This may induce currents that could lead to shocks, fires, or injury to the patient or user.
- 16- Due to concerns about the carcinogenic and infectious potential of electrosurgical byproducts (such as tissue smoke and aerosols), protective eyewear, filtration masks, and effective smoke evacuation equipment should be used in both open and laparoscopic procedures.
- 17- Always provide a safe working environment. Avoid sharp edges and trip hazards. Avoid kinking or sharply bending the cables or damaging the insulation. Failure to do so may result in damage to the insulation and result in burns with the unwanted electric current passes.
- 18- Always plug and unplug the HF cables by holding on to the connector. Do not touch the metal parts of the connector. Never pull at the cable.
- 19- Cables whose connectors have not been thoroughly rinsed and dried/aerated may cause electrical shocks or burns to the patient or user.
- 20- The patient monitoring electrodes, such as ECG electrodes, should be placed as far apart from the surgical site as feasible. Avoid using ECG needle electrodes for monitoring due to the risk of patient burns from the alternate electrical pathway. All ECG electrodes must be equipped with protective impedances or HF choke coils.
- 21- In operations where the high-frequency current is passing from the very minor parts of the body, it is preferred to use the bipolar technique instead of monopolar use to avoid unwanted coagulation.
- 22- In Monopolar electrosurgery, patient laying position, the use and the proper placement of a dispersive (neutral) electrode is a key element in the safe and effective use, particularly in the prevention of burns.
- 23- Apparent low output or failure of the generator to function correctly at the normal operating settings may indicate faulty application of the dispersive (neutral) electrode or poor contact in its connections. In this case, prior to increasing the output power, check the adherence of the dispersive (neutral) electrode and its connections.
- 24- For patients with electrically conductive implants, a possible hazard exists due to the concentration or re-direction of HF currents. In case of doubt, qualified advice should be obtained.
- 25- Do not use with hybrid trocar systems, i.e., a combination of metal and plastic, when using monopolar active components. This may result in alternate site burns due to capacitive coupling. Use only all-metal or all-plastic trocar systems.
- 26- Check this instruction manual, the label, or the current product catalogue for the maximum rated voltage of the HF cable. In any combination with another electrosurgical accessory, the maximum rated voltage of the combination corresponds to the lowest rated voltage of the accessories used. Do not use HF cables with a generator set to output voltage above the specified maximum accessory-rated voltage.
- 27- As a rule, HF electrodes combined with standard resectoscopes should only be used at a maximum output voltage of 2.0 kVp in the usual cut and coagulation modes.
- 28- The Contact Quality Monitoring System (CQMS, REM, ARM, NESSY, etc.) will not function if the incorrect dispersive cable is used. The type of generator will determine which cable should be used.
- 29- Neuromuscular stimulation can cause unintentional movement of the patient.
- 30- Reprocessing and mechanical stress damage the HF cable depending on the intensity of use. All reusable HF cables should not exceed 50 uses. Damage caused by misuse, overuse, mishandling, or normal wear may cause electrical shocks or burns to the patient or user.
- 31- If this product is used on a patient with or suspected of having CJD (Creutzfeldt-Jacob disease), vCJD (variant Creutzfeldt-Jacob disease), BSE (Bovine Spongiform Encephalopathy), and TSE (Transmissible Spongiform Encephalopathy) the product cannot be reused and must be destroyed to eliminate the risk of cross-contamination. Such pathogens cannot be killed by standard disinfection and sterilization methods.
- 32- Never attempt to perform repairs. Do not modify the HF cables in any way. Modifications may affect patient safety, and performance, and void the warranty.
- 33- Do not use damaged bipolar/monopolar cables. Dispose of damaged cables according to national and local laws and regulations for medical waste.
- 34- Refer to the manual of the high-frequency generator for additional cautions and warnings.

Manual Cleaning and Disinfection:

Manual Cleaning: The HF cable must be cleaned and disinfected immediately after each use. Remove heavy debris from the HF cable by wiping it with an appropriate, moist, lint-free cloth, sponge, or gauze pad. Prepare the detergent solution with deionized or distilled water according to the manufacturer's instructions. Saturate a synthetic felt or sponge pad with the cleaning solution. Wipe all surfaces of the HF cable until all visible contamination is removed. Wash the pad in the cleaning solution and repeat wiping the cable. Use a plastic cleaning brush for adhered particles or dried secretions. At least allow a contact time of 2 minutes. If visible contamination is still present on the cable at the end of the cleaning repeat the cleaning process. Moist a clean cloth, sponge, or gauze pad with deionized or distilled water and thoroughly wipe the HF cable to remove any debris and all traces of the cleaning solution. Dry manually by using a lint-free cloth or absorbent paper and, in particular, for drying cavities, sterile compressed air.

Manual Disinfection: Prepare the disinfectant solution according to the manufacturer's instructions. Only disinfect clean products that have dried slightly. Saturate a synthetic felt or sponge pad with the disinfectant solution. Wipe all surfaces of the HF cable. Wash the pad in the disinfectant solution and repeat wiping the cable. At least allow a contact time of 6 minutes. Moist a clean cloth, sponge, or gauze pad with deionized or distilled water and thoroughly wipe the HF cable to remove any residual debris and all traces of disinfectant solution. Dry manually by using a lint-free cloth or absorbent paper and, in particular, for drying cavities, sterile compressed air.

HF cable should never be immersed in cleaning, disinfecting, or rinsing solutions, or placed in an ultrasonic cleaning device. Soaking will reduce the life of HF cable by oxidation of metal parts and hardening of the cable jacket. Avoid using excessive force when cleaning or disinfecting by hand. It is not recommended to use saline solutions as a final rinse as saline solutions may interfere with subsequent disinfection or sterilization steps. Please do not clean with hard or frictional materials, abrasive cleaners, or solvents. Cleaning agents containing phenol or chloride are not suitable. These instructions were validated using the below detergents and disinfectants.

Cleaning agent (manual cleaning): Enzo® / Cidezyme®, enzymatic detergent, Johnson&Johnson (add 8ml enzymatic detergent per 1 liter distilled water, for cable with dried-on organic matter use 16 ml enzymatic detergent per 1 liter distilled water/warm water)

Disinfectant (manual disinfection): Cidex OPA®, 0.55% ortho-phthalaldehyde solution, Johnson&Johnson (at a minimum temperature of 20°C/68°F).

Reverse osmosis/distilled water.

Machine Cleaning and Disinfection:

Place the pre-cleaned HF cables in a suitable container and put them into the washer and disinfectant machine (WDM). Make sure that HF cables do not touch each other and are not kinked or crushed. Do not clean together with sharp-edged or pointed objects.

Only a washing and disinfecting machine with a CE mark (or FDA approval) and validated according to ISO15883, with a tested program for thermal disinfection (at least 90°C for 5 minutes or A₀ value > 3000) should be used. It is recommended to use of neutral to mild alkaline cleaning agents. Discoloration of metal surfaces is possible when alkaline cleaners are used (pH 9.5 - 11.5). However, this does not impair the functional capability of the cable. Strong alkaline cleaners (pH>11.5) must be avoided. Do not use neutralization agents.

Select the program according to the WDM manufacturer's instructions. The following cycle is recommended with these parameters.

- 1- Pre-rinse: < 45°C for 2 minutes with softened water
- 2- Cleaning: at 55°C for 5 minutes with demineralized water and 0.5 % (5 ml/liter) mild alkaline detergent (MetalClean, Ecolab Ltd.)
- 3- Intermediate rinse: < 45°C for 2 minutes with demineralized water
- 4- Neutralization / Last rinse: > 40°C for minimum 5 minutes with demineralized water
- 5- Thermal disinfection / Final rinse: at 90°C for 5 minutes with demineralized water
- 6- Drying: 80°C for minimum 18 minutes

At the end of the program, remove the HF cables from the WDM. Dry the HF cables manually with a clean, soft cloth or sterile compressed air (< 2 bar) if necessary. Check all the surfaces of HF cables for residues. If there is visible residue, repeat the cleaning process. Immediately pack the HF cable in a single-use sterilization pouch or put it in a sterilization container.

Caution: During the cleaning and disinfection of HF cable disposable gloves, protective eyewear, and filtration mask should be used to reduce the risk of transmission of infectious agents by splashing or the creation of aerosols.

Warning: FMT high-frequency bipolar/monopolar cables are sold NON-STERILE. Clean and disinfect the cables before the first and every use.

Warning: Before cleaning or disinfecting the cable, disconnect it from the high-frequency generator.

Warning: Avoid contact with strong, aromatic, chlorinated, ketone, ether, or ester solvents. Do not use bleach on electrical contacts or connectors.

Warning: Avoid using flammable cleaning and disinfecting agents. If flammable agents are used, allow them to evaporate before HF surgery.

Sterilization: The HF cables must be cleaned and disinfected before sterilization. When preparing the HF cable for sterilization, do not coil the cable into a circle of less than 10cm (4") in diameter. Pack the cable in disposable sterilization packaging (single or double packaging) made from paper/polyethylene and/or put it in a sterilization container. Packaging must comply with the ISO 11607-1 and EN 868-5 standards for packaging sterilized instruments. Sterilization of the HF cable should be done by steam sterilization using the fractioned pre-vacuumed method. The cables must be sterilized at a minimum of 134°C [273°F] and a maximum of 137°C [278°F] in saturated steam during a holding (exposure) time of at least 4 minutes to at most 20 minutes, then dried in a vacuum for at least 20 minutes. When sterilizing several products at the same time in a steam autoclave, ensure that the maximum load capacity of the autoclave specified by the manufacturer is not exceeded. Sterilization must be performed in accordance with EN ISO 17665 (Sterilization of Medical Devices in Moist Heat).

The above sterilization process has been validated to assure a Sterility Assurance Level of 10⁻⁶ (SAL ≤ 1:10⁶) and the proof of the effectiveness of steam sterilization has been provided by an independent, accredited test laboratory (test report number AYA-SVR100620-01 dated 10.06.2020). The use of other sterilization methods (e.g. ethylene oxide gas, formaldehyde, gamma radiation, or low-temperature plasma sterilization) or other autoclave sterilization cycles is outside the manufacturer's responsibility.

Note: Before use, the HF cable must be cooled to room temperature. It is the sole responsibility of the user to safeguard and maintain the sterile condition of the cable after the sterilization process.

Note: The user can track the number of reprocessing cycles by punching a hole every time the product is autoclaved with the permanently attached autoclave record tag on the cable.

Warning: Do not sterilize with hot air. Do not use "Flash" autoclave procedures.

Compatibility: To ensure compatibility, the HF cable should only be used with the equipment and instrument for which they have been designed and is specified on the label on the product's primary packaging. For product models and compatible devices, please refer to the catalog at metkomedical.com/PDF/06ElectrosurgeryCablesandAccessories.pdf

Environmental Requirements:

Operating Temperature : 5°C - 45°C (41°F - 113°F)

Storage Temperature : 0°C - 50°C (32°F - 122°F)

Relative Humidity : 20 - 80% (not condensing)

Store the cables in a clean, cool, and dry place. During storage, the products should be protected from sunlight. The storage of the product in its original packaging until its first use is recommended. Used cables should be loosely coiled in 10 to 15 cm diameter (4 to 6 inches) and must be stored individually in a protective container or sterile packaging.

Safety:

Degree of protection from electric shocks: type CF

Degree of protection against the ingress of water: Keep Dry

FMT reusable high-frequency bipolar/monopolar cables are not made with natural rubber latex and PVC.

Refer to the separate leaflet "Description of Symbols" located within the product package for the explanation of symbols.

Reporting of serious incidents: Any serious incident related to the use of this product should be reported to both the manufacturer and the health authority/competent authority where the product is installed. Either contact your local representative or report to: metko@metkomedical.com

Caution: Federal Law (U.S.A.) restricts this device to sale by or on the order of a physician.

Warranty: Each reusable high-frequency cable is designed to withstand a minimum of 50 sterilization cycles when adequately cared for and sterilized according to instructions. All reusable high-frequency bipolar/monopolar cables are warranted for 50 uses against defects in workmanship and materials. Damage caused by misuse, overuse, wrong sterilization, or normal wear is not covered under warranty. After the first use, the product has a service life of 12 months. Dispose of the HF cable after 12 months.

FMT® is a registered trademark of Metko A.Ş.

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