# Reusable Electrosurgical Adaptor

English

**Indications:** Reusable electrosurgical adaptors are intended for providing the electrical connection between the electrosurgical generator (high-frequency unit) and foot control or hand switch electrosurgical pencils, or other monopolar or bipolar cables, by adapting the physical form of one connector to another. Electrosurgical adaptors transfer high-frequency electrical power/energy without changing. Reusable electrosurgical adaptors are intended to be used with appropriate electrosurgical equipment by trained and qualified medical personnel in hospitals and hospital-type facilities.

**Contraindications:** While not an absolute contraindication to electrosurgery, implantable electrical devices (IEDs) such as cardiac and gastric pacemakers, implantable cardiac defibrillators (ICDs), cochlear implants, and deep brain, nerve, spinal cord, or bone stimulators create unique risks. HF current during electrosurgery can interfere with the functioning of these devices and damage them. Patients with these devices require a thorough pre-procedure evaluation and may require intraoperative monitoring and post-procedure device assessment.

**Adverse Reactions / Events:** Adverse events have been reported when using electrosurgical devices, such as tissue damage at the wrong site caused by unintentional activation, fires involving surgical drapes and other flammable materials, burns caused by alternating current paths where the patient or user is in contact with exposed metal, explosions caused by an electrosurgical spark in the flammable gas environment, organ perforation and sudden massive bleeding.

**Clinical Benefits:** Enables conducting of high-frequency current during electrosurgery.

**Product Description:** FMT reusable electrosurgical adaptors are manufactured from Thermoplastic Elastomer (TPE) and silicone. Electrosurgical adaptors have different instrument-side connectors, device-side connectors, and connector wiring configurations for different electrosurgical generators (high-frequency units). Therefore, the same appearance does not guarantee full compatibility with the electrosurgical generator. Always consider the information on the label of the product and/or its primary packaging.

- ESU-ADP/01**, Monopolar adaptor for hand switch pencils with 3-pin Valleylab type connector, compatible with Erbe T Series generators
- ESU-ADP/02**, Monopolar adaptor for hand switch pencils with 3-pin Valleylab type connector, compatible with Erbe ACC/ICC Series generators
- ESU-ADP/03**, Monopolar adaptor for hand switch pencils with 3-pin Valleylab type connector, compatible with Berchtold - Martin generators
- ESU-ADP/04**, Bipolar adaptor for bipolar cable with Ø4.0mm banana plug, compatible with Alsa Excel generators
- ESU-ADP/05**, Bipolar adaptor for bipolar cable with Ø4.0mm banana plug, compatible with KLS Martin generators

Maximum accessory rated voltage	(U <sub>max</sub> ):
Electrosurgical bipolar adaptors	0.8 kVp
Electrosurgical monopolar adaptors	4.0 kVp

**Directions for Use:** Electrosurgical adaptors and their connectors should be inspected for visual defects. Loose or damaged connector mating parts (broken or missing components), corrosion of metal contacts, damage in the adaptor's insulation (i.e. cuts, cracks, brittle or solid, permanent bends, burnt or melted), wires or metal parts that have come out of the cable/connector insulation, loose connectors, and not completely dried and ventilated adaptors should not be used.

- 1- Electrosurgical adaptors transmit the HF current required to cut or coagulate tissues to high-frequency electrodes, pencils, forceps, or instruments. Therefore, a suitable HF electrode or instrument or HF cable should be connected to the adaptor's socket side. The adaptor's plug connects to the high-frequency unit/generator.
- 2- Check whether the adaptor plug is suitable for the correct socket on the generator. Careless pushing or pushing with force may damage the plug and/or the socket on the generator.
- 3- The instrument connector(s) should be plugged into the socket(s) of the electrosurgical adaptor, covering all the metal contacts of the instrument connector.
- 4- When connecting the adaptor, it is important to check whether good and continuous contact is obtained. This is obtained by pushing the adaptor plug into the mechanical borders of the socket on the generator, covering all the metal contacts.
- 5- Before disconnecting the adaptors, always switch the HF unit off. When connecting and disconnecting adaptors, hold only the plastic part of the plug. Otherwise, it may result in a burn or electric shock in patients or users. If it is available on the adaptor, do not pull on the cable. In preparation, transport, sterilization, and storage, adaptor cables should not be bent, twisted tightly, or folded. Such movements reduce the life of the product. Do not load heavy or sharp objects on the adaptor. This may damage the adaptor insulation. Damage to the insulation may result in burns with the unwanted electric current passes.
- 6- Activation of the HF current is done by the foot or hand switch or, if available, by the unit's autostart mode. Directives in the HF unit's user manual should be followed.

**Special Information:** The socket on the unit for bipolar handles with 2-line (bipolar) cables can be mixed with other bipolar coagulation instruments. For example, the bipolar coagulation connection cable in ERBE fits the monopolar cut and coagulation socket of high-frequency units of AESCULAP, BERCHTOLD, INTEGRA, and MARTIN. For this reason before use make sure that the product is correct and suitable for the HF unit.

#### Warnings:

- 1- The user of electrosurgical adaptors should have the necessary knowledge, expertise, and education to use the products. Read all instruction manuals of the electrosurgical adaptors, generators, and connected instruments thoroughly. An insufficient understanding of the dangers, warnings, cautions, and information in these manuals can result in death, serious injury, or equipment damage.
- 2- Before the first and each use, an electrical continuity test must be performed and the insulation of the electrosurgical adaptor must be inspected for damage. If there is a break in the cable wire or the adaptor becomes otherwise electrically discontinuous, arcing may occur in the patient-return or active circuit and may burn the patient or user or create a fire. Visually inspect the entire adaptor. Visual inspection alone may not be sufficient to ensure that the insulation is intact. Inspection may be done visually under magnification or with a high-voltage insulation testing device. Do not use electrosurgical adaptors with brittle or defective insulation or exhibiting pressure points.
- 3- Reusable electrosurgical adaptors should be inspected, cleaned, sterilized, and dried/aerated prior to use per this instruction manual prior to first and each subsequent use. The user is recommended to implement a program to track the number of uses of the adaptors. The serial number on the product allows the user to keep track of the adaptor and record reprocessing cycles or the service life of each adaptor.
- 4- Properly reprocess the electrosurgical adaptor before the first and each subsequent use following the instructions in this manual. Improper and/or incomplete reprocessing can cause infection in the patient and/or medical personnel.
- 5- Do not perform electrosurgery in the presence of flammable anesthetics or oxidizing gases, other flammable gases, or in close proximity to volatile solvents, fluids, or objects, or in the oxygen-enriched environments, as explosions and fire could result. Do not place electrosurgical instruments near or in contact with flammable materials (such as gauze or surgical drapes). Instruments that are activated or hot from use may cause a fire.
- 6- High-frequency accessories connected to the electrosurgery unit but not temporarily used should not be touched to the patient and held in a different position. The active electrode should be placed in a clean, dry, well-insulated safety holster when not in use to prevent inadvertent activation or contamination. Inadvertent contact with the patient may result in burns.
- 7- During the use of high-frequency equipment in patients with cardiac pacemakers, implanted heart defibrillators, or pace electrodes there is a danger to ruin the function or the device itself. The pacemaker or electrodes should never be touched with high-frequency electrodes. Before using advice from cardiology departments must always be taken.
- 8- Electrosurgical adaptors should only be connected to the high-frequency generator when the unit is switched off. Make sure that the adaptor is not connected when preparing the electrosurgical instrument (working element, etc.). Failure to do so may result in injury or electrical shock to the patient or user.
- 9- Avoid unnecessary or prolonged activation. This may cause injuries due to capacitive coupling with other surgical equipment. Activate only when touching or immediately next to the intended point of contact. Keep the voltage/power settings at the lowest possible level to achieve the desired effect. Keep the active electrodes clean. When the effectiveness of the instrument is reduced, do not increase the power of the generator. Instead, clean the active electrode contacts.
- 10- Aspirate fluid from the area before activating the HF instrument. Conductive fluids (e.g., blood or saline) in direct contact with or in close proximity to an active electrode may carry electrical current or heat away from target tissues, which may cause unintended burns to the patient. Flammable solutions pooling under the patient or in body depressions may cause burns to the patient.
- 11- Potentially hazardous conditions may exist when accessories are not appropriate for the type of generator output used. If unusual operating characteristics develop or are observed during procedures, the use of the product should be discontinued.
- 12- Connect bipolar adaptors to the bipolar receptacle and monopolar adaptors to the monopolar receptacle only. Improper connection of adaptors may result in inadvertent accessory activation or other potentially hazardous conditions.
- 13- Make sure that the patient is electrically isolated against electrical conductive elements and does not come into direct contact with other metal parts or electric conductive instruments which are grounded or which have an appreciable capacitance to ground. Metal devices can provide an alternate pathway for the electrical current and be responsible for patient burns. Avoid skin-to-skin contact of the patient.
- 14- Always provide a safe working environment. Avoid sharp edges and trip hazards. Avoid kinking or sharply bending the cables or damaging the insulation. Failure to do so may result in damage to the insulation and result in burns with the unwanted electric current passes.
- 15- Always plug and unplug the adaptors by holding on to the connector. Do not touch the metal parts of the connector. Never pull at the cable.
- 16- Adaptors whose connectors have not been thoroughly rinsed and dried/aerated may cause electrical shocks or burns to the patient or user.
- 17- The patient monitoring electrodes, such as ECG electrodes, should be placed as far apart from the surgical site as feasible. Do not use ECG needle electrodes for monitoring. Needle monitoring electrodes provide an alternate pathway for the electrical current placing the patient at risk for burns. All ECG electrodes must be equipped with protective impedances or HF choke coils.
- 18- In operations where the high-frequency current is passing from the very minor parts of the body, it is preferred to use the bipolar technique instead of monopolar use to avoid unwanted coagulation.
- 19- Apparent low output or failure of the generator to function correctly at the normal operating settings may indicate faulty application of the dispersive (neutral) electrode or poor contact in its connections. In this case, prior to increasing the output power, check the adherence of the dispersive (neutral) electrode and its connections.
- 20- For patients with electrically conductive implants, a possible hazard exists due to the concentration or re-direction of HF currents. In case of doubt, qualified advice should be obtained.
- 21- Check this instruction manual, the label, or the current product catalog for the maximum rated voltage of the electrosurgical adaptor. If anything is unclear, contact the manufacturer. In any combination with another electrosurgical accessory, the maximum rated voltage of the combination corresponds to the lowest rated voltage of the accessories used. Do not use electrosurgical adaptors with a generator set to output voltage above the specified maximum accessory-rated voltage.
- 22- Reprocessing and mechanical stress damage the electrosurgical adaptor depending on the intensity of use. All reusable electrosurgical adaptors should not exceed 50 uses. Damage caused by misuse, overuse, mishandling, or normal wear may cause electrical shocks or burns to the patient and/or user.
- 23- If this product is used on a patient with or suspected of having CJD (Creutzfeldt-Jakob disease), vCJD (variant Creutzfeldt-Jakob disease), BSE (Bovine Spongiform Encephalopathy), or TSE (Transmissible spongiform Encephalopathy), it cannot be reused and must be destroyed to eliminate the risk of cross-contamination. Standard disinfection and sterilization methods cannot kill such pathogens.
- 24- In Monopolar electrosurgery, patient laying position, the use and the proper placement of a dispersive (neutral) electrode is a key element in the safe and effective use, particularly in the prevention of burns.
- 25- The Contact Quality Monitoring System (CQMS, REM, ARM, NESSY, etc.) will not function if the incorrect dispersive cable is used. The type of generator will determine which cable should be used.
- 26- Never attempt to perform repairs. Do not modify the adaptors in any way. Modifications may affect patient safety, and performance, and void the warranty.
- 27- Do not use damaged electrosurgical adaptors. Dispose of damaged adaptors according to national and local laws and regulations for medical waste.
- 28- Refer to the manual of the high-frequency generator for additional cautions and warnings.

#### Manual Cleaning and Disinfection:

**Manual Cleaning:** The electrosurgical adaptors must be cleaned and disinfected immediately after each use. Remove heavy debris from the adaptor by wiping with an appropriate, moist lint-free cloth, sponge, or gauze pad. Prepare the detergent solution with deionized or distilled water according to the manufacturer's instructions. Saturate a synthetic felt or sponge pad with the cleaning solution. Wipe all surfaces of the adaptor until all visible contamination is removed. Wash the pad in the cleaning solution and repeat wiping the adaptor. Use a plastic cleaning brush for adhered particles or dried secretions. Wipe or brush all surfaces at least six strokes back and forth. At least allow a contact time of 2 minutes. If visible contamination is still present on the adaptor at the end of the cleaning

repeat the cleaning process or safely dispose of the adaptor. Moist a clean cloth, sponge, or gauze pad with deionized or distilled water and thoroughly wipe the adaptor to remove any debris and all traces of the cleaning solution. Dry manually by using a lint free cloth or absorbent paper and, in particular, for drying cavities, sterile compressed air.

**Manual Disinfection:** Prepare the disinfectant solution according to the manufacturer's instructions. Only disinfect clean products that have dried slightly. Saturate a synthetic felt or sponge pad with the disinfectant solution. Wipe all surfaces of the electrosurgical adaptor. Wash the pad in the disinfectant solution and repeat wiping the adaptor. Wipe or brush all surfaces at least six strokes back and forth. At least allow a contact time of 6 minutes. Moist a clean cloth, sponge, or gauze pad with deionized or distilled water and thoroughly wipe the adaptor to remove any residual debris and all traces of disinfectant solution. Dry manually by using a lint-free cloth or absorbent paper and, in particular, for drying cavities, sterile compressed air.

Electrosurgical adaptors should never be immersed in cleaning, disinfecting, or rinsing solutions, or placed in an ultrasonic cleaning device. Soaking will reduce the life of the adaptor by oxidation of metal parts. Avoid using excessive force when cleaning or disinfecting by hand. This may reduce the adaptor life. It is not recommended to use saline solutions as a final rinse as saline solutions may interfere with subsequent disinfection or sterilization steps. Please do not clean with hard or frictional materials, abrasive cleaners, or solvents. Cleaning agents containing phenol or chloride are not suitable. These instructions were validated using the below detergents and disinfectants.

Cleaning agent (manual cleaning): Enzol® / Cidezyme®, enzymatic detergent, Johnson&Johnson (add 8ml enzymatic detergent per 1 liter distilled water, for cable with dried-on organic matter use 16 ml enzymatic detergent per 1 liter distilled water / warm water)

Disinfectant (manual disinfection): Cidex OPA®, 0.55% ortho-phthalaldehyde solution, Johnson&Johnson (at a minimum temperature of 20°C/68°F).

Reverse osmosis / distilled water.

**Caution:** During the cleaning and disinfection of electrosurgical adaptors disposable gloves, protective eyewear, and filtration mask should be used to reduce the risk of transmission of infectious agents by splashing or the creation of aerosols.

**Warning:** FMT electrosurgical adaptors are sold NON-STERILE. Clean and disinfect the adaptors before the first and every use.

**Warning:** Before cleaning or disinfecting the adaptor, disconnect it from the high-frequency generator.

**Warning:** Electrosurgical adaptor should be cleaned of overburden prior to disinfection to improve the effectiveness.

**Warning:** Avoid contact with strong, aromatic, chlorinated, ketone, ether, or ester solvents. Do not use bleach on electrical contacts or connectors.

**Warning:** Avoid using flammable cleaning and disinfecting agents. If flammable agents are used, allow them to evaporate before HF surgery.

**Sterilization:** The electrosurgical adaptors must be cleaned and disinfected prior to sterilization. Pack the adaptor in disposable sterilization packaging (single or double packaging) made from paper/polyethylene and/or put it in a sterilization container. Packaging must comply with the ISO 11607-1 and EN 868-5 standards for packaging sterilized instruments. Sterilization of the electrosurgery adaptors should be done by steam sterilization using the fractionated pre-vacuumed method. The adaptor must be sterilized at a minimum of 134°C [273°F] and a maximum of 137°C [278°F] in saturated steam during a holding (exposure) time of at least 4 minutes to at most 20 minutes, then dried in a vacuum for at least 20 minutes. When sterilizing several products at the same time in a steam autoclave, ensure that the maximum load capacity of the autoclave specified by the manufacturer is not exceeded. Sterilization must be performed in accordance with EN ISO 17665 (Sterilization of Medical Devices in Moist Heat).

The above sterilization process has been validated to assure a Sterility Assurance Level of  $10^{-6}$  ( $SAL \leq 1:10^6$ ) and the proof of the effectiveness of steam sterilization has been provided by an independent, accredited test laboratory (test report number AYA-SVR100620-01 dated 10.06.2020). The use of other sterilization methods (e.g. ethylene oxide gas, formaldehyde, gamma radiation, or low-temperature plasma sterilization) or other autoclave sterilization cycles is outside the manufacturer's responsibility.

**Note:** The electrosurgical adaptors must be cooled to room temperature before use. After the sterilization process, it is the user's sole responsibility to safeguard and maintain the adaptors sterile condition.

**Warning:** Do not sterilize with hot air. Do not use "Flash" autoclave procedures.

**Compatibility:** In order to ensure compatibility, the reusable electrosurgical adaptor should only be used with the equipment and instrument they have been designed and is specified in a label on the product's primary packaging. For product models and compatible devices, please refer to the catalog at [metkomedical.com/PDF/06ElectrosurgeryCablesandAccessories.pdf](http://metkomedical.com/PDF/06ElectrosurgeryCablesandAccessories.pdf)

#### Environmental Requirements:

Operating Temperature : 5°C - 45°C (41°F - 113°F)

Storage Temperature : 0°C - 50°C (32°F - 122°F)

Relative Humidity : 20 - 80% (not condensing)

Store the adaptors in a clean, cool, and dry place. During storage, the product should be protected from sunlight. The storage of the product in its original packaging until its first use is recommended.

Used adaptors should be stored individually in a protective container or sterile packaging.

#### Safety:

Degree of protection from electric shocks: type CF

Degree of protection against the ingress of water: Keep Dry

FMT reusable electrosurgical adaptors are sold non-sterile. They must be cleaned, disinfected, and sterilized before initial use and each subsequent use.

FMT reusable electrosurgical adaptors are not made with natural rubber latex and PVC.

Refer to the separate leaflet "Description of Symbols" located within the product package for the explanation of symbols.

**Reporting of Serious Incidents:** Any serious incident related to the use of this product should be reported to both the manufacturer and the health authority/competent authority where the product is installed. Either contact your local representative or report to: [metko@metkomedical.com](mailto:metko@metkomedical.com)

**Caution:** Federal Law (U.S.A.) restricts this device to sale by or on the order of a physician.

**Warranty:** Each reusable electrosurgical adaptor is designed to withstand a minimum of 50 sterilization cycles when properly cared for and sterilized according to given instructions. All reusable electrosurgical adaptors are warranted for 50 uses against defects in workmanship and materials. Damage caused by misuse, overuse, mishandling, or normal wear is not covered under warranty. After the first use, the product has a service life of 12 months.

FMT® is a registered trademark of Metko A.Ş.

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