**Indications:** Disposable pulse oximetry (SpO<sub>2</sub>) sensors are indicated for the continuous noninvasive monitoring of functional oxygen saturation of arterial hemoglobin (SpO<sub>2</sub>), pulse rate, and plethysmographic pulse wave of the patient. Disposable (SpO<sub>2</sub>) sensors are intended to be used with appropriate pulse oximeters or patient monitors by trained and qualified medical personnel in hospitals and hospital-type facilities.

**Contraindications:** Disposable SpO<sub>2</sub> sensors are contraindicated for use on patients with allergic reactions to adhesive tape or for prolonged periods of use. The sensor site must be inspected at least every four (4) hours or sooner, and if the circulatory condition or skin integrity is compromised, the sensor should be applied to a different site. Disposable pulse oximeter sensors are also contraindicated if they are used on other than the specified patient size (patient weight) and application site.

**Clinical Benefits:** A wired LED and PD sensor enables continuous measurement of patient functional oxygen saturation.

**Product Description:** The compatible disposable SpO<sub>2</sub> sensors are intended for use with major pulse oximeter brands and are designed to match the specifications of the original equipment manufacturer. There are three types of disposable SpO<sub>2</sub> sensors, according to the structures covering optoelectronic sensors; M type microfoam elastic foam tape with an adhesive layer, F type woven fabric tape with an adhesive layer, and DVF type foam sponge with Velcro. Disposable SpO<sub>2</sub> sensors are four sizes according to patient weight. DAF sensors are used with Adult patients, DPF sensors with Pediatric patients, and DIF sensors with Infant patients. DNF and DVF sensors can be used in Neonatal or Adult patients depending on the application site. The disposable SpO<sub>2</sub> sensor is for multiple uses on a single patient and is packaged sterile. Disposable SpO<sub>2</sub> sensors can be used on the same patient as long as the adhesive tape attaches without slippage.


Compatible Brand or Pulse Oximeter Technology	Adult		Pediatric		Infant		Neonatal /Adult	
	> 40 kg	10-50 kg	10-50 kg	3-20 kg	Neonatal < 3 kg	Adult > 40 kg	Neonatal Hand/Foot,	Adult Finger Application
	Finger Application	Finger Application	Finger Application	Toe/Thumb Application	Neonatal Hand/Foot,	Adult Finger Application	Neonatal Hand/Foot,	Adult Finger Application
BCI	FMT-DAF / BCI-X	FMT-DPF / BCI-X	FMT-DIF / BCI-X	FMT-DNF / BCI-X	FMT-DNF / BCI-X	FMT-DVF / BCI		
CSI (CRITICARE)	FMT-DAF / CSI-X	FMT-DPF / CSI-X	FMT-DIF / CSI-X	FMT-DNF / CSI-X	FMT-DNF / CSI-X	FMT-DVF / CSI		
DATEx	FMT-DAF / DTX-X	FMT-DPF / DTX-X	FMT-DIF / DTX-X	FMT-DNF / DTX-X	FMT-DNF / DTX-X	FMT-DVF / DTX		
GE OXY TIP (UN Connector)	FMT-DAF / OXY-X	FMT-DPF / OXY-X	FMT-DIF / OXY-X	FMT-DNF / OXY-X	FMT-DNF / OXY-X	FMT-DVF / OXY		
GE TRUSIGNAL	FMT-DAF / OXN-X	FMT-DPF / OXN-X	FMT-DIF / OXN-X	FMT-DNF / OXN-X	FMT-DNF / OXN-X	FMT-DVF / OXN		
MASIMO (LNCS)	FMT-DAF / MSM-X	FMT-DPF / MSM-X	FMT-DIF / MSM-X	FMT-DNF / MSM-X	FMT-DNF / MSM-X	FMT-DVF / MSM		
MASIMO (LNOP)	FMT-DAF / MSD-X	FMT-DPF / MSD-X	FMT-DIF / MSD-X	FMT-DNF / MSD-X	FMT-DNF / MSD-X	FMT-DVF / MSD		
MEK	FMT-DAF / MEK-X	FMT-DPF / MEK-X	FMT-DIF / MEK-X	FMT-DNF / MEK-X	FMT-DNF / MEK-X	FMT-DVF / MEK		
NELLCOR	FMT-DAF / NLC-X	FMT-DPF / NLC-X	FMT-DIF / NLC-X	FMT-DNF / NLC-X	FMT-DNF / NLC-X	FMT-DVF / NLC		
NELLCOR OXIMAX	FMT-DAF / NLO-X	FMT-DPF / NLO-X	FMT-DIF / NLO-X	FMT-DNF / NLO-X	FMT-DNF / NLO-X	FMT-DVF / NLO		
NIHON KOHDEN	FMT-DAF / NKN-X	FMT-DPF / NKN-X	FMT-DIF / NKN-X	FMT-DNF / NKN-X	FMT-DNF / NKN-X	FMT-DVF / NKN		
NOVAMETRIX	FMT-DAF / NVM-X	FMT-DPF / NVM-X	FMT-DIF / NVM-X	FMT-DNF / NVM-X	FMT-DNF / NVM-X	FMT-DVF / NVM		
NONIN	FMT-DAF / NON-X	FMT-DPF / NON-X	FMT-DIF / NON-X	FMT-DNF / NON-X	FMT-DNF / NON-X	FMT-DVF / NON		
OHMEDA	FMT-DAF / OHM-X	FMT-DPF / OHM-X	FMT-DIF / OHM-X	FMT-DNF / OHM-X	FMT-DNF / OHM-X	FMT-DVF / OHM		
Adhesive sensors (X: material code of the sensors, M for microfoam medical tape, and F for fabric medical tape).							Non-adhesive sensor. Foam sponge sensors (DVF)	

**Comparative value measurement in % saturation:** Between 70% and 100% SpO<sub>2</sub> levels, SpO<sub>2</sub> accuracy (Functional Oxygen Saturation Arms) is ±2%, and SpO<sub>2</sub> accuracy is not specified under the 70% SpO<sub>2</sub> levels. Pulse Rate accuracy (Arms) from 35 to 240 bpm is ± 2 bpm.

**Directions for Use:**


**Applying the DAF (Adult) and DPF (Pediatric) sensors to the patient:**

- 1- Open the pouch and remove the sensor. Remove the plastic backing from the DAF/DPF sensor and locate transparent windows (a) on the adhesive side. Windows cover optical components.
- 2- Select a site that least restricts a conscious patient's movements. The patient's index finger of the non-dominant hand is the preferred location. Alternatively, the middle finger, ring finger, or little finger of the hand may be used. The thumb or the great toe may be used as secondary alternatives.
- 3- Orient the DAF/DPF sensor so the dashed line in the middle of the sensor is centered on the tip of the digit. Wrap adhesive flaps on the non-cable end around the digit. Note that the cable must be positioned on the top of the hand.
- 4- Fold the cable end over the top of the digit so that the windows are directly opposite each other. Wrap adhesive flaps securely around the sides of the digit.
- 5- Plug the DAF/DPF sensor connector into the adapter & extension cable socket or the pulse oximeter and verify proper operation as described in the oximeter operator's manual.


**Applying the DIF (Infant) sensor to the patient:**

- 1- Open the pouch and remove the sensor. Remove the plastic backing from the DIF sensor and locate transparent windows (b) on the adhesive side. Windows cover optical components. Note corresponding alignment marks (a) on the non-adhesive side and the dashed line (c) midway between the marks.
- 2- The preferred sensor site is the great toe for infants. Alternatively, apply the sensor to another digit of similar size, for example, a thumb.
- 3- Orient the DIF sensor so the window next to the cable is aligned on the bottom of the great toe as shown. The cable should extend towards the heel.
- 4- Wrap the DIF sensor firmly, but not too tightly around the toe. Windows must directly oppose each other.
- 5- Wrap any excess tape loosely around the toe. Secure the sensor cable with medical tape, preferably across the bottom of the foot. Make sure that the tape securing the cable does not restrict blood circulation.
- 6- Plug the DIF sensor connector into the adapter & extension cable socket or the pulse oximeter and verify proper operation as described in the oximeter operator's manual.


**Applying the DNF (Neonatal / Adult) sensor to the patient:**

- 1- Open the pouch and remove the sensor. Remove the plastic backing from the DNF sensor and locate transparent windows (b) on the adhesive side. Windows cover optical components. Note corresponding alignment marks (a) on the non-adhesive side and the dashed line (c) midway between the marks.
- 2- Orient the DNF sensor so the dashed line is on the lateral edge of the site (c).
- Neonate:** The preferred sensor site is the fleshy part of a foot for neonates. Alternatively, apply the sensor to a hand palm. The window next to the cable goes sole of the foot as shown.
- Adult:** The preferred sensor site is an index finger for adults. Alternatively, apply the sensor to another finger. The window next to the cable goes on the nail side, distal to the first joint. Do not place the sensor on a joint. Note that the cable must be positioned on the top of the hand.
- 3- Wrap the DNF sensor firmly but not too tightly around the foot or finger. Windows must directly oppose each other.
- 4- Plug the DNF sensor connector into the adapter & extension cable socket or the pulse oximeter and verify proper operation as described in the oximeter operator's manual.

**Applying the DVF Non-Adhesive Sponge (Neonatal / Adult) sensor to the patient:**

- 1- Open the pouch and remove the DVF sensor. Locate transparent windows (b) on the blue-colored sponge side. Windows cover optical components.
- 2- Position the optical window that aligns with the sensor cable over the selected extremity of the patient.
- Neonate:** The preferred sensor site is the fleshy part of a foot for neonates. Alternatively, apply the sensor to a hand palm. The window next to the cable goes sole of the foot as shown.
- Adult:** The preferred sensor site is an index finger for adults. Alternatively, apply the sensor to another finger. The window next to the cable goes on the nail side, distal to the first joint. Do not place the sensor on a joint. Note that the cable must be positioned on the top of the hand.
- 3- Wrap the DVF sensor around the extremity firmly but not too tightly, so that both optical windows oppose each other. Use the Velcro® at the end of the sensor (d) to secure it in place. If necessary, secure the cable to the hand or foot with medical tape.
- 4- Plug the DVF sensor connector into the adapter & extension cable socket or the pulse oximeter and verify proper operation as described in the oximeter operator's manual.

DVF-Neonatal / Adult

**Warnings:**

- 1- All SpO2 sensors and adapter cables are designed for use with specific patient monitors or pulse oximeters. Verify the compatibility of the monitor, adapter cable, and sensor and all functions including alarming perform properly before use, otherwise degraded performance and/or patient injury can result.
- 2- The application of disposable SpO2 sensors should be done by trained medical personnel following established procedures.
- 3- All wire-lead patient-connected transducer assemblies are subject to reading error, local heating, and possible damage from high-intensity sources of RF energy. Electrosurgical equipment represents one such source since capacitively-coupled currents may seek alternative paths to ground through sensor cables and isolated instruments. Patient burns may result. If possible, remove the sensor from patient contact before activating the surgical unit or other RF source. If sensors must be used simultaneously with electrosurgical apparatus, the instruments to which the sensors are connected should be checked for adequate isolation from electrical grounds at radio frequencies. Hazards can be reduced by selecting a monitoring point located away from the expected RF current path between the active probe and the return pad.
- 4- The sensor should be free of visible defects, discoloration and damage. If the sensor is wet, discolored, or damaged, discontinue use. Never use a damaged sensor or one with exposed electrical circuitry.
- 5- The SpO2 sensor is in direct contact with the patient's skin and therefore allergic reactions may occur in susceptible individuals. In case of allergic or skin reactions discontinue use and consult a specialist.
- 6- All SpO2 sensors may be used on the same site for a maximum of four (4) hours, provided the site is inspected routinely per clinical protocol to ensure skin integrity, circulatory status, and correct optical alignment. Skin erosion, tissue ischemia, and/or pressure necrosis can be caused when the sensor is not frequently moved, applied too tightly, or becomes too tight due to edema. Assess the site as frequently as every (1) hour and move the sensor or use a different type of sensor if there are signs of pyrexia or loss of skin integrity and/or loss of circulation or perfusion.
- 7- As with all medical equipment, carefully route cables to reduce the possibility of patient entanglement or strangulation.
- 8- Venous congestion may cause under the reading of actual arterial oxygen saturation. Therefore, assure proper venous outflow from the monitored site. The sensor should not be below heart level (e.g. sensor on the hand of a patient in a bed with an arm dangling to the floor).
- 9- Venous pulsations may cause inaccurate readings (e.g. tricuspid valve regurgitation, Trendelenburg position).
- 10- Do not use the sensor during Magnetic Resonance Imaging (MRI) scanning or in an MRI environment as it may result in physical harm.
- 11- Misapplied sensors or sensors that become partially dislodged may cause incorrect readings.
- 12- Inaccurate readings may be caused by EMI radiation interference or strong electromagnetic sources such as electrosurgery equipment.
- 13- Abnormal fingers, intravascular dyes such as indocyanine green or methylene blue, or externally applied coloring and texture such as fingernail polish, or artificial fingernails may lead to inaccurate or no readings.
- 14- The application site should be cleaned of debris prior to sensor placement. Do not apply the sensor to the injured skin.
- 15- Elevated levels of Carboxyhemoglobin (COHb), Methemoglobin (MetHb), and total bilirubin may lead to inaccurate SpO2 readings. High levels of COHb or MetHb may occur with a seemingly normal SpO2. When elevated levels of COHb or MetHb are suspected, laboratory analysis (CO-Oximetry) of a blood sample should be performed.
- 16- Inaccurate SpO2 readings may be caused by severe anemia, very low arterial perfusion, vasospastic disease such as Raynaud's, and peripheral vascular disease at the monitored site or extreme motion artifact.
- 17- Inaccurate SpO2 readings may be caused by elevated levels of dyshemoglobin, hypocapnic or hypercapnic conditions, and severe vasoconstriction or hypothermia.
- 18- Hemoglobinopathies and synthesis disorders such as thalassemias, Hb s, Hb c, sickle cell, etc. may cause inaccurate SpO2 readings.
- 19- If using pulse oximetry during full body irradiation, keep the sensor out of the radiation field. If the sensor is exposed to the radiation, the reading might be inaccurate or not provided or the unit might read zero for the duration of the active radiation period.
- 20- The pulsations from intra-aortic balloon support can be additive to the pulse rate on the oximeter pulse rate display. Verify pulse rate against the ECG heart rate.
- 21- Sensors should not compress the limb. A tight sensor will compromise venous return and may give rise to inaccurate readings. When selecting a sensor type it is important to check if the patient's weight is in the recommended patient weight range.
- 22- High ambient light sources such as surgical lights (especially those with a xenon light source), bilirubin lamps, fluorescent lights, infrared heating lamps, and direct sunlight can interfere with the performance of the sensor. To prevent interference from ambient light, ensure that the sensor is properly applied, and cover the sensor site with opaque material if required. Failure to take this precaution in high ambient light conditions may result in inaccurate readings.
- 23- At elevated ambient temperatures, the patient's skin could be severely burned after prolonged sensor application at sites that are not well perfused. To prevent this condition, be sure to check patient application sites frequently. All listed sensors operate without risk of exceeding 41°C on the skin if the initial skin temperature does not exceed 35°C.
- 24- When selecting the sensor site, priority should be given to an extremity free of an arterial catheter, blood pressure cuff, or intravascular infusion line.
- 25- Do not use tape to secure the sensor to the site; this can restrict blood flow and cause inaccurate readings. The use of tape can cause skin damage, and/or pressure necrosis or damage to the sensor.
- 26- The performance of the sensor is compromised by motion. In such cases, try to keep the patient still, or change the sensor site to one with less motion.
- 27- Do not twist unnecessarily or use excessive force when using, connecting, disconnecting, or storing the sensor. Do not pull on the cable.
- 28- High oxygen concentration may predispose a premature infant to retinopathy. Therefore, the upper alarm limit for oxygen saturation must be carefully selected in accordance with accepted clinical standards. The accuracy range of the oximeter being used must be considered.
- 29- Portable and mobile RF communications equipment can affect medical electrical equipment and should be used no closer to any part of the equipment.
- 30- The sensor is not protected against the effect of the discharge of a cardiac defibrillator. Remove the sensor before cardiac defibrillation.
- 31- Keep the sensor dry. Do not use a wet sensor. This may cause burns during the application of high-frequency devices.
- 32- Never attempt to perform repairs. Do not modify the SpO2 sensors in any way. Modifications adversely affect patient safety, performance, and/or accuracy, and void the warranty.
- 33- Do not attempt to re-sterilize disposable sensors. In the event of damage to the sterile packaging, do not use or re-sterilize.
- 34- Open the single-use product sterile package just before application. Do not use dropped sterile products on the patient.
- 35- Sterile disposable sensors should be used before the expiry date has been exceeded.
- 36- Do not re-use disposable SpO2 sensors on multiple patients. Dispose of used or damaged sensors according to local laws and regulations for medical waste.
- 37- Refer to the monitor's operation manual for additional cautions and warnings. For operation errors and screen error messages, refer to the troubleshooting table provided in the operation manual of the monitor.

**Cleaning, Disinfection and Sterilization:** Disposable SpO2 sensors are intended for multiple uses on a single patient. Single patient use SpO2 sensors probes should be disposed of after use on a single patient, so do not attempt to clean, disinfect or sterilize them.

**Measurement Wavelengths and Output Power (\*)**

LED Type	Red Peak Wavelength	Red Max. Output Power	IR Peak Wavelength	IR Max. Output Power
2-Leads	663nm	@ 1.2mW nominal	890nm	@ 1.0mW nominal
3-Leads	661nm	@ 1.2mW nominal	940nm	@ 1.2mW nominal
4-Leads	660nm	@ 1.2mW nominal	905/940nm	@ 1.0mW nominal

\*This information is especially useful for clinicians performing photodynamic therapy.

**Environmental Requirements:**

Operating Temperature : 5°C - 45°C (41°F - 113°F)  
 Storage Temperature : 0°C - 50°C (32°F - 122°F)  
 Relative Humidity : 20 - 80% (not condensing)

During the storage the products should be protected from sunlight. It is recommended to store the products in original sterile packages until the first use.

**Compatibility:** In order to ensure compatibility and claimed accuracy, the disposable SpO2 sensor should only be used with the equipment for which they have been designed and is specified in the label on the product's primary packaging. The disposable SpO2 sensor is designed to be compatible with the patient monitors and pulse oximeters using the appropriate adapter and extension cable. For product models and compatible devices, please refer to the catalog at [metkomedical.com/PDF/03PulseOximetry\(SpO2\)Sensors.pdf](http://metkomedical.com/PDF/03PulseOximetry(SpO2)Sensors.pdf)

**Safety:**

Degree of protection from electric shocks: type BF  
 Degree of protection against the ingress of water: IPX2  
 The FMT disposable SpO2 sensors are sold STERILE. Sterilized by Ethylene Oxide (ETO) sterilization method.  
 The FMT disposable SpO2 sensors are not made with natural rubber latex.  
 The FMT disposable SpO2 sensors are not made with Diethylhexyl phthalate (DEHP) plasticizer.  
 The FMT disposable SpO2 sensors have not been evaluated for safety in the MR environment. Scanning a patient who has this device may result in patient injury.  
 Refer to the separate leaflet "Description of Symbols" located within the product package for the explanation of symbols.  
**Note:** Functional Tester cannot be used to assess the accuracy of a pulse oximeter sensor or a pulse oximeter monitor.

**Reporting of serious incidents:** Any serious incident related to the use of this product should be reported to both the manufacturer and the health authority/competent authority where the product is installed. Either contact your local representative or report to: [metko@metkomedical.com](mailto:metko@metkomedical.com)

**Caution:** Federal Law (U.S.A.) restricts this device to sale by or on the order of a physician.

**Warranty:** FMT disposable SpO2 sensors are under six (6) months warranty against material and workmanship defects from the date of original purchase. In warranty period, METKO will be responsible for changing the sensor free of charge if the defect is proven. Single use products are warranted for single patient use only. This warranty does not extend to any product that has been subject to misuse, neglect or accident; or that has been damaged by causes external to the product; or that has been used in violation of the operating instructions supplied with the product. Product life is 3 years.

FMT® is a registered trademark of Metko A.Ş.

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