**Indications:** FMT disposable neonatal Non-Invasive Blood Pressure (NIBP) cuffs are intended for use with automated oscillometric NIBP devices to measure the systolic and diastolic arterial pressure of the patient and should be disposed of after use on a single patient. NIBP cuffs are intended to be used by qualified medical personnel in hospitals and hospital-type facilities. Neonatal NIBP cuffs can be used on neonates and infants depending on selection of the appropriate cuff size.

**Contraindications:** Measurement of blood pressure at the brachial artery is a generally benign procedure. However, there are some circumstances in which obtaining readings from a particular arm may not be appropriate; such circumstances include the presence of an arterial-venous shunt, recent axillary node dissection, or any deformity or surgical history that interferes with proper access or blood flow to the upper arm. If these relative contraindications are present, blood pressure should be assessed in the opposite arm. There may also be pre-existing conditions that can interfere with the accuracy or interpretation of readings, such as aortic coarctation, arterial-venous malformation, occlusive arterial disease, or the presence of an antecubital bruit. If neither arm can be used, then measurement of blood pressure in a leg may be indicated.

**Adverse Reactions:** No allergic reactions are observed at places where NIBP cuffs are in contact with the intact skin of a patient. However, the such reaction cannot be completely excluded in people particularly susceptible to allergies; it should be noted that no such events have been reported to the manufacturer. Incorrect cleaning or disinfection, or a failure to perform them, as well as failure to observe hygiene procedures may cause the transfer of pathogens and infections.

**Clinical Benefits:** Enables NIBP monitoring with patient monitoring equipment.


**Product Description:** FMT disposable neonatal NIBP cuffs are manufactured from Thermoplastic Polyurethane (TPU) and in bladderless form. Cuff sizes are marked with different animal figures and color-coded to allow easy size selection. There are five sizes depending on patient size. The choice of single-tube or dual-tube cuff depends on the NIBP monitoring system. Connectors attached to the end of the cuff tube (or tubes) allow the cuff to be connected to a compatible NIBP monitor via the NIBP hose.

**FMT- DSX/BPYY** Disposable neonatal cuff, Single tube


**FMT- DDX/BPYYYY** Disposable neonatal cuff, Double tube

(X is cuff size 1: Neonatal Small, 2: Neonatal Medium, 3: Neonatal Large, 4: Neonatal Extra Large, and 5: Infant size, YY is connector number)

Cuff Text Color and Symbol	Navy Blue - Sea Horse 	Light Blue - Penguin 	Green - Seal 	Red - Duck 	Orange - Horse 
	#1 Neonatal Small	#2 Neonatal Medium	#3 Neonatal Large	#4 Neonatal Extra Large	#5 Infant
Cuff Range (Circumference)	3 - 6 cm	4 - 8 cm	6 - 11 cm	7 - 13 cm	8 - 15 cm
Cuff Dimensions (L x W)	10.0 x 2.5 cm	11.9 x 3.1 cm	14.6 x 4.1 cm	18.1 x 4.8 cm	19.1 x 5.6 cm

**Directions for Use:**

- 1- Before use, check that the neonatal NIBP cuff, cuff tubing, and hose are intact and clean. Replace cuff when aging, tearing, or weak closure is apparent. Do not inflate the cuff when unwrapped. Do not use cuff if structural integrity is suspect.
- 2- Select the appropriate blood pressure measurement site. Inspect the patient's limb prior to application. Do not apply the cuff to areas where the skin is not intact or tissue is injured. Because normative values are generally based on this site and as a matter of convenience, the upper arm is the preferred site for the placement of the cuff. The forearm, thigh, and ankle can also be used. If the largest cuff does not fit the upper arm adequately, use an alternate site. When using the upper arm choose the correct cuff accordingly, taking into account the patient's cardiovascular status and the effect of the alternative site on blood pressure values, proper cuff size, and comfort.
- 3- Ensure that cuffed limb is supported to keep the cuff level with the patient's heart. If the cuff is not at heart level, the difference in systolic and diastolic values due to the hydrostatic effect must be considered.
- 4- Select the appropriate cuff size. Measure the length of the limb in centimeters (cm) from joint to joint. The width of the cuff that you select should be in the range of 55% to 65% of the length of the limb. Avoid using a cuff with a width that extends over a joint. The cuff should be long enough to encircle the circumference of the limb completely. Applicable RANGE printed in centimeters and marked on the cuff. When cuff sizes overlap for a specified circumference, choose the larger size cuff. Accuracy depends on the use of the proper size cuff.
- 5- Remove clothing from the limb before positioning the cuff. Place the cuff on a bare limb. Place the opened cuff to the limb such that the side with the animal figure is away from the patient. Ensure the side with the warning "THIS SIDE TO PATIENT" is toward the patient.
- 6- Wrap the cuff snugly around the limb ensuring the INDEX LINE on the edge of the cuff falls within the range indicated by the RANGE lines on the cuff, then fasten. The INDEX LINE at the edge of the cuff should reach into, but not extend beyond, the RANGE marking on the cuff. If they do not line up, use a different size cuff. Make sure to align the artery mark (  ) over the patient's brachial (or femoral for ankle) artery.
- 7- The cuff should not be so tight as to prevent venous return between determinations. Ensure that hook and loop closures are properly engaged so that pressure is evenly distributed throughout the cuff.
- 8- Connect the cuff to the NIBP hose ensuring the proper engagement of the connectors. Ensure that the connection is secure and that tubing is not kinked and the hose is routed to avoid kinking or compression.
- 9- Start NIBP monitoring. Confirm the artery marker continues pointing to the brachial artery after each inflation. Reposition the cuff to the correct position if it slips.

**Warnings:**

- 1- Neonatal NIBP cuffs are designed for use with compatible automated oscillometric NIBP devices to measure arterial pressure on the neonatal patient's upper arm. Verify the compatibility of the NIBP cuff, NIBP hose, and NIBP monitoring device all functions perform properly before use, otherwise degraded performance and/or patient injury can result.
- 2- Limb circumference ranges are printed on each cuff. Be sure the cuff is correctly matched to patient size. Cuffs that are too small may produce false high blood pressure measurements; cuffs that are too large may produce false low-pressure measurements.
- 3- Before each use, visually inspect the NIBP cuff and connected NIBP hose. Use only monitor hoses identified for use with neonatal cuffs. The use of other hoses can prevent monitors from obtaining readings. The NIBP cuffs should be free of visible defects and damage. Replace cuffs if hook and loop fail to hold during inflation, artery, index, or range markings are not visible, or the monitor displays a message that a leak has been detected. Never use damaged, clogged, or visibly contaminated NIBP cuffs. Dispose of damaged NIBP cuffs according to hospital disposal practices, local environmental and waste disposal laws, and regulations.
- 4- Do not start NIBP measurements unless you are sure that the cuff inflation and alarm settings are appropriate for the patient. Incorrect settings can result in patient injury or inaccurate measurements. Do not apply the higher adult inflation, overpressure limits, and measurement duration to neonatal patients.
- 5- Do not connect the cuff to intravascular fluid systems which could allow air to be pumped into a blood vessel, which could lead to serious patient injury.
- 6- Do not apply cuffs to areas where the skin is not intact or tissue is injured, or where dermal disruption is at greater risk. Ensure that the rough side of the closure does not contact with skin; contact may cause irritation.
- 7- Do not obtain determinations more frequently than clinically indicated, weighing the benefits of frequent measurement against risk. Rapidly repeating measurements can impair circulation in the monitored limb.

- 8- Check NIBP cuff/adaptor/NIBP hose, cuff site, and cuff limb frequently. Signs of skin irritation or impeded blood flow, especially when monitoring at frequent intervals and/or extended periods of time should be checked frequently. Rotate the site if appropriate.
- 9- Do not allow an NIBP cuff to remain on the patient for more than 10 minutes when inflated above 10 mm Hg. This may cause patient distress, disturb blood circulation, and contribute to the injury of peripheral nerves. Remove the cuff from the patient when monitoring is not in process.
- 10- Devices that exert pressure on tissue have been associated with purpura, skin avulsion, compartmental syndrome, ischemia, and/or neuropathy.
- 11- Do not apply the cuff to a limb used for intravenous infusion, arterial monitoring, or any other intravascular access, therapy, or an arterio-venous (A-V) shunt, where AV fistulas are present or in areas where circulation is compromised. The cuff inflation can temporarily block blood flow, potentially causing harm to the patient. Assess limb for the risk of lymphedema (due to mastectomy, etc.).
- 12- Use care when placing the cuff on limbs used to monitor other patient parameters. Do not apply the cuff to a limb being used for SpO2 monitoring. Cuff inflation may cause inaccurate SpO2 measurements.
- 13- Do not apply external pressure against the cuff and avoid contact with the cuff while monitoring. This may cause inaccurate blood pressure values.
- 14- Use cuff only under direct supervision by a trained healthcare professional when attached to automated monitors without fail-to-read alarms.
- 15- To optimize performance and accuracy, minimize limb movement/cuff motion.
- 16- Avoid placing the cuff on a patient in a manner that can lead to the cuff tube becoming kinked. Tube kinking may cause inaccurate readings.
- 17- Do not modify the Neonatal NIBP cuffs in any way. Modifications may affect patient safety, performance, and accuracy and void the warranty.
- 18- As with all medical equipment, carefully route hoses to reduce the possibility of patient entanglement or strangulation.
- 19- Do not use disposable cuffs on multiple patients.
- 20- For further information and warnings read the instruction manual of the NIBP monitoring device or patient monitor before using the NIBP cuffs.

**Cleaning and Disinfection:** Disposable neonatal NIBP Cuffs are intended for multiple uses on a single patient. Disposable cuffs should be disposed of after use on a single patient, or if they become soiled, so do not attempt to clean, disinfect or sterilize them.

**Warning:** Do not clean, or disinfect disposable neonatal NIBP cuffs. Do not sterilize by autoclave, radiation or steam.

**Warning:** Care should be taken to prevent liquid from entering the cuff tubing. The liquid in the tubing may affect blood pressure determination accuracy and damage automatic or manual monitors.

**Environmental Requirements:**

Operating Temperature : 5°C - 45°C (41°F - 113°F)

Storage Temperature : 0°C - 50°C (32°F - 122°F)

Relative Humidity : 20 - 80% (not condensing)

During storage, the product should be protected from sunlight. It is recommended to store the product in its original package until the first use.

**Compatibility:** In order to ensure compatibility, the Neonatal NIBP cuff should only be used with the equipment for which they have been designed and is specified in the label on the product's primary packaging. For product models and compatible devices, please refer to the catalog at [metkomedical.com/PDF/05NIBPCuffsandHoses.pdf](http://metkomedical.com/PDF/05NIBPCuffsandHoses.pdf)

**Safety:**

Degree of protection from electric shocks: type BF

Degree of protection against the ingress of water: Keep dry

The FMT disposable neonatal NIBP cuffs are sold NON-STERILE.

The FMT disposable neonatal NIBP cuffs are not made with natural rubber latex and PVC.

The FMT disposable neonatal NIBP cuffs have not been evaluated for safety in the MR environment. Scanning a patient who has this device may result in patient injury.

Refer to the separate leaflet "Description of Symbols" located within the product package for the explanation of symbols.

**Reporting of serious incidents:** Any serious incident related to the use of this product should be reported to both the manufacturer and the health authority/competent authority where the product is installed. Either contact your local representative or report to: [metko@metkomedical.com](mailto:metko@metkomedical.com)

**Caution:** Federal Law (U.S.A.) restricts this device to sale by or on the order of a physician.

**Warranty:** FMT disposable neonatal NIBP cuffs are under six (6) months warranty against material and workmanship defects from the date of original purchase. In the warranty period, METKO will be responsible for repairing or replacing the cuff free of charge if the defect is proven. This warranty does not cover any product that has been subjected to misuse, neglect, or accident and that has been damaged by causes external to the product, or that has been used in violation of the operating instructions supplied with the product. The product life is 3 years of the production date.

FMT® is a registered trademark of Metko A.Ş.

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