Indications:

The disposable skin temperature probe is indicated for routine continuous monitoring of the skin temperature of infant patients. The disposable temperature probes are intended to be used with compatible incubators, infant warmers, or phototherapy units by trained and qualified medical personnel in hospitals and hospital-type facilities.

Contraindications:

The use of a disposable skin temperature probe is contraindicated for usage over traumatized and hypo-perfused areas.

Adverse Reactions:

Adverse reactions reported during applications of skin temperature probes include skin sensitizations or irritations due to adhesive cover.

Clinical Benefits:

A wired temperature sensor enables continuous measurement of patient skin temperature in a warmer.

Product Description:

FMT disposable probes are intended for use with compatible electronic thermometers. The use of the temperature probe with an incompatible electronic thermometer can result in degraded performance and accuracy. The operator is responsible for checking the compatibility of the monitor and the temperature probe before its use. Probes are accurate and interchangeable with any other probe in the same series at specified measurement temperatures. Within the measurement range of 25 - 45 °C, the accuracy for the probes is ± 0.1°C. The minimum measurement time for accurate temperature readings is 150 seconds. The temperature reaches a stabilized temperature in 5 minutes, depending on the body site. The system accuracy is the sum of the probe and instrument accuracies.

The disposable skin temperature probe is designed for placement on the skin surface. In order to fix the sensor and minimize the effect of the ambient temperatures on temperature measurement, the sensor is covered with a heat-reflective adhesive cover. The heat-reflective adhesive cover is made of metallic-looking Mylar film, closed-cell foam, and backing paper to protect the adhesive surface. The temperature sensor is placed in the center of the foam cut out to prevent pressure on the skin below the sensor. The adhesive on the foam surface holds the skin sensor in situ. The disposable temperature probe is for single-patient use only and is individually packaged.


FMT-INC/SD Skin temperature probe for Air Shields / Hillrom - Drager Incubators (OEM P/N 68 209 20 – MU06933)

FMT-INC/SD2 Skin temperature probe for Drager Incubators (OEM P/N MX11000, MX11001)

FMT-INC/SD3 Skin temperature probe for GE Giraffe Incubators (OEM P/N 6600 0873 700)

Directions for Use:

- 1- Inspect the packaging and product for any damage, contamination, or missing parts prior to use.
- 2- Remove the probe from its package.
- 3- Make sure that the probe is compatible with the device to be connected and that the plug of the probe matches the temperature socket of the device. To connect the probe to the device, align the plug of the probe with the temperature socket and push firmly to ensure full contact. Forced mating of connectors without proper alignment can result in damage to connectors and loss of electrical continuity.
- 4- Clean and dry skin completely in the area intended for probe placement.
- 5- When the infant is on its back, the probe should be placed on the abdomen, halfway between the xiphoid and the umbilicus (1). When the infant is in the prone position, the probe should be on the infant's back. If the baby is placed in the lateral position, ensure the skin probe is located on the side of the infant and clear of the infant's arm (2). The probe must never be placed under the infant or used rectally.
- 6- Peel off the backing paper under the probe cover and apply the adhesive side to the skin. Fix the skin probe to the selected skin site by gently pressing the edges of the probe cover.
- 7- A piece of medical tape can be placed over the probe wire approximately 4 or 5 centimeters from the probe tip to stabilize the attached probe.
- 8- Wait to see the complete measurement reading after the response time is reached. Follow the monitoring device manufacturer's instructions of use and your hospital's protocols on temperature monitoring of the patients.
- 9- When temperature monitoring is completed, first remove the heat-reflective adhesive cover then carefully remove the skin probe tip from the patient's skin then disconnect the probe from the device. Hold the plastic plug firmly and pull. Do not pull on the cable.
- 10- This probe is for single-patient use. Discard after use according to facility protocol.

Warnings:

- 1- All temperature probes are designed for use with specific patient monitors or electronic thermometers. Verify the compatibility of the monitor, extension cable, and probe and all functions including alarming perform properly before use, otherwise degraded performance, accuracy, and/or patient injury can result.
- 2- The application of temperature probes should be done by trained medical personnel following established procedures.
- 3- All wire-lead patient-connected transducer assemblies are subject to reading error, local heating, and possible damage from high-intensity sources of RF energy. Electrosurgical equipment represents one such source since capacitively-coupled currents may seek alternative paths to ground through probe cables and isolated instruments. Patient burns may result. If possible, remove the probe from patient contact before activating the surgical unit or other RF source. If probes must be used simultaneously with electrosurgical apparatus, the instruments to which the probes are connected should be checked for adequate isolation from electrical grounds at radio frequencies. Hazards can be reduced by selecting a temperature monitoring point located away from the expected RF current path between the active probe and the return pad.
- 4- It is not recommended to use the temperature probes within the range of Magnetic Resonance Imaging (MRI), Computed Tomography (CT), etc. applications. Conducted current may cause burns.
- 5- The temperature probe is in direct contact with the patient's skin and therefore allergic reactions may occur in susceptible individuals. In case of allergic or skin reactions discontinue use and consult a specialist.
- 6- As with all medical equipment, carefully route temperature probe cables to reduce the possibility of patient entanglement or strangulation.
- 7- Assure that probe cable is not positioned where they could become entangled around IV tubes and restrict IV flow.
- 8- The medical thermometer cannot measure the patient's temperature accurately if the probe is not correctly positioned or the probe tip is not firmly contacting the body site. Regularly check the probe remains correctly positioned on the patient and probe is always in direct contact with the patient's skin.
- 9- Make sure the skin temperature probe is free from pressure. Never place skin temperature probes underneath the patient.

- 10- Do not use a skin probe without a reflective probe cover. Do not remove the reflective foil from the reflective probe cover. Do not use any type of lotion or dressing between patient skin and adhesive probe cover. Replace the probe cover when repositioning the skin probe or when adhesive strength degrades.
- 11- Place skin temperature probe on soft muscle area. Avoid any bony area.
- 12- Do not apply the temperature probe to a skin surface that has been damaged by abrasion, laceration, incision, or burning.
- 13- When monitoring the infants, if the infant is on its back, the skin probe should be placed on the abdomen, halfway between the xyphoid and the umbilicus. When the infant is in the prone position, the skin probe should be on the infant's back. If the baby is placed in the lateral position, ensure the skin probe is located on the side of the infant and clear of the infant's arm.
- 14- To ensure the continued comfort of the infant, avoid pressing directly on the skin probe when attaching it to the infant. Always use caution when applying or removing a temperature probe from a patient.
- 15- Do not leave the patient unattended when using the incubator. Check the patient's temperature periodically to ensure the comfort and safety of the patient.
- 16- Do not place the incubator in direct sunlight or in proximity to any other source of radiant heat. External heat sources can cause elevated interior temperatures that may overheat the infant. Phototherapy units located too close to the incubator may affect the hood wall temperature, incubator temperature, and infant skin temperature. Patient injury could occur. Do not place phototherapy units too close to the incubator.
- 17- When using phototherapy lamps, the probe must be directly in the path of the radiant heat of the lamp, do not place the probe in an area shielded from the lamp's light. The phototherapy lamp may raise the infant's skin temperature.
- 18- Do not use an auxiliary or rectal temperature measurement to control the infant's skin temperature. It can result in overheating or underheating of the infant.
- 19- Mishandling of the probes could result in damage to internal wires and loss of electrical isolation or improper temperature readings. Don't tie the probes in knots or bunch tightly. Don't wrap probes around equipment to avoid damaging internal wires. Don't leave the probes dangling from the instrument.
- 20- Care must be taken to ensure that the probe cable is not subjected to mechanical stress, such as pulling, bending, and squeezing.
- 21- Portable and mobile RF communications equipment can affect medical electrical equipment and should be used no closer to any part of the equipment.
- 22- The probe is not protected against the effect of the discharge of a cardiac defibrillator. Remove the probe before cardiac defibrillation.
- 23- Never attempt to perform repairs. Do not modify the temperature probes in any way. Modifications may affect patient safety, and performance, and void the warranty.
- 24- Do not attempt to sterilize disposable probes
- 25- Do not re-use disposable probes on multiple patients. Dispose of used or damaged probes according to local laws and regulations for medical waste.
- 26- For further information and warnings read the instructions of use accompanying the temperature monitoring device.

Cleaning, Disinfection, and Sterilization:

Disposable temperature probes are intended for multiple uses on a single patient. Disposable temperature probes should be disposed of after use on a single patient, so do not attempt to clean, disinfect or sterilize them.

Environmental Requirements:

Operating Temperature : 5°C - 45°C (41°F - 113°F)

Storage Temperature : 0°C - 50°C (32°F - 122°F)

Relative Humidity : 20 - 80% (not condensing)

During the storage the products should be protected from sunlight. It is recommended to store the products in original packages until the first use.

Compatibility:

In order to ensure compatibility and claimed accuracy, the temperature probe should only be used with the equipment for which they have been designed and is specified in the label on the product's primary packaging. For product models and compatible devices, please refer to the catalog at metkomedical.com/PDF/04MedicalTemperatureProbes.pdf

Safety:

Degree of protection from electric shocks: type BF

Degree of protection against the ingress of water: IPX2

The FMT disposable temperature probes are sold NON-STERILE.

The FMT disposable temperature probes are not made with natural rubber latex.

The FMT disposable temperature probes are not made with a Diethylhexyl phthalate (DEHP) plasticizer.

The FMT disposable temperature probes have not been evaluated for safety in the MR environment. Scanning a patient who has this device may result in patient injury.

Refer to the separate leaflet "Description of Symbols" located within the product package for the explanation of symbols.

Reporting of serious incidents:

Any serious incident related to the use of this product should be reported to both the manufacturer and the health authority/competent authority where the product is installed. Either contact your local representative or report to: metko@metkomedical.com

Caution: Federal Law (U.S.A.) restricts this device to sale by or on the order of a physician.

Warranty:

FMT disposable temperature probes are under six (6) months warranty against material and workmanship defects from the date of original purchase. In warranty period, METKO will be responsible for changing the probe free of charge if the defect is proven. Single-use products are warranted for single-patient use only. This warranty does not extend to any product that has been subject to misuse, neglect, or accident or that has been damaged by causes external to the product, or that has been used in violation of the operating instructions supplied with the product. Product life is 3 years.

FMT® is a registered trademark of Metko A.Ş.

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