Reusable Diagnostic EKG Cable

English

Indications: The standard 12-channel electrocardiogram represents the heart's electrical activity with electrical signals measured from ECG electrodes placed on the patient's skin. Diagnostic EKG cables are used to transmit a patient's electrocardiographic signals from electrodes to an electrocardiograph, a device that processes these signals. EKG cables are intended to be used by qualified medical personnel on infants, children, and adults in hospitals and hospital-type facilities.

Contraindications: No known contraindications to the use of the product.

Adverse Reactions: No allergic reactions are observed at places where EKG cables are in contact with a patient's intact skin. However, such reactions cannot be completely excluded in people particularly susceptible to allergies. It should be noted that no such events have been reported to the manufacturer. Incorrect cleaning or disinfection, or a failure to perform them, as well as failure to observe hygiene procedures, may cause the transfer of pathogens and infections.

Clinical Benefits: Enables ECG recording with diagnostic EKG equipment.

Product Description: FMT EKG cables have different device-side connectors and connector wiring configurations for different electrocardiograph manufacturers. Therefore, the same appearance does not guarantee full compatibility with the electrocardiograph. Always consider the information on the product's label and/or primary packaging. There are Ø4.0mm banana, Ø3.2 mm straight pin, snap, or grabber-type patient terminations available.

E100/xxxY/Z, Diagnostic EKG cable with fixed patient leadwires.

E101/xxxY/Z, Diagnostic EKG cable with detachable patient leadwires.

E103/xxxY/Z, Diagnostic EKG cable with detachable patient leadwire set (GE Multi-Link type).

E104/xxxY/Z, Diagnostic EKG cable with detachable patient leadwire set (GE Multi-Link 12SL type).

E105/xxxY/Z, Diagnostic EKG cable with detachable patient leadwire set (PHILIPS M1949A type).

E106/xxxY/Z, Diagnostic EKG cable with detachable patient leadwire set (Medtronic/PhysioControl type).

E107/xxxY/Z, Diagnostic EKG cable with detachable patient leadwire set (Schiller type).

(xxx shows the compatible brand. Y shows the patient termination type B=4.0mm Banana, P=3.2mm Straight Pin, S=Snap, G=Grabber, O=trunk cable w/o leadwires. Z shows the color code I=IEC, A=AHA/AAMI)

Directions for Use:

- 1- Before use, check that the product is intact and clean.
- 2- Connect the EKG cable connector to the socket of the compatible electrocardiograph and check for correct operation.
- 3- Before applying the electrodes, ensure the patient's skin is clean, dry, and properly prepared. Always apply the electrodes to intact skin.
- 4- Attach leadwire terminations to the electrodes. Ensure the leadwires do not pass twist stress onto electrodes.
- 5- Apply the electrodes firmly on prepared sites using your hospital's protocol or the chart given below.
- 6- Ensure the leads conform to the body contours and no strain is placed on the electrodes.
- 7- Verify that the ECG signal appears on the electrocardiograph screen. Please refer to the electrocardiograph instructions for use for detailed information on ECG recording.

IEC Color Code		AHA/AAMI Color Code		Location (Resting EKG)
R	Red	RA	White	Right Arm, avoiding bony prominences.
L	Yellow	LA	Black	Left Arm, avoiding bony prominences.
F	Green	LL	Red	Left Leg, avoiding bony prominences.
N	Black	RL	Green	Right Leg, avoiding bony prominences.
C1	Red	V1	Red	Fourth intercostal space at the right margin of sternum.
C2	Yellow	V2	Yellow	Fourth intercostal space at the left margin of sternum.
C3	Green	V3	Green	Midway between position C2 (V2) and position C4 (V4).
C4	Brown	V4	Blue	Fifth intercostal space at the junction of the left midclavicular line.
C5	Black	V5	Orange	At the horizontal level of position C4 (V4) at the left anterior axillary line.
C6	Purple	V6	Purple	At the horizontal level of position C4 (V4) at the left mid-axillary line.

Warnings:

- 1- All EKG cables are designed for use with a specific electrocardiograph. Before use, verify the compatibility of the electrocardiograph and EKG cable and that all functions function properly; otherwise, degraded performance may result.
- 2- Before each use, visually inspect the EKG cable. The EKG cable should be free of visible defects and damage. If there are breaks, scratches, or cracks on the cable jacket or plastic insulation of the connector or patient terminations or visibly bent, damaged pins of the connector or patient terminal, discontinue using the cable. Never use a damaged cable or one with exposed electrical circuitry. Even if there is no exposed electrical circuit, the scratches or cracks in the cable jacket or plastic insulation may cause the accumulation of pathogens and reduce the effectiveness of cleaning and disinfection of the product. Dispose of damaged cables according to hospital disposal practices, and local environmental and waste disposal laws and regulations.
- 3- Do not use the EKG cable during Magnetic Resonance Imaging (MRI) scanning or in an MRI environment, as it may cause physical harm.
- 4- Keep the EKG cable connectors as far as possible from electromagnetic sources, such as an electromedical device's power transformer. Also, avoid placing the EKG cable alongside the device's power cord. Strong electromagnetic fields may adversely interfere with the function of an EKG cable.
- 5- As with all medical equipment, carefully route cables to reduce the possibility of patient entanglement or strangulation.
- 6- Do not connect the EKG cable to any power outlet to prevent shock hazards.
- 7- To avoid a shock hazard and interface with nearby electrical equipment, keep electrodes and patient cables away from grounded metal and other electrical equipment.
- 8- When connecting and disconnecting the EKG cable and the patient terminations, pull on the connector or plug itself, not on the cable.
- 9- Do not modify the EKG cable in any way. Modifications may affect patient safety, performance, and accuracy and void the warranty.
- 10- Mishandling of the EKG cables could result in damage to internal wires or loss of electrical isolation. Do not tie the EKG cable and leadwires in knots or bunch tightly. Do not wrap the EKG cable around the equipment to avoid damaging internal wires. Do not leave the EKG cable dangling from the instrument. Keep EKG cable leadwires away from moving parts that might pinch or crush them. When not in use, loosely coil the EKG cable to a bend diameter of at least 10 cm and store it at room temperature.
- 11- If all the leadwires are not connected correctly with the patient, that may cause reading errors.
- 12- Dirty, compromised or unprepared skin may cause incorrect ECG readings. Using previously opened or out-of-date electrodes may degrade the ECG signal quality. ECG electrodes should be removed from their sealed package immediately prior to use.
- 13- Always ensure that the patient is kept motionless during 12 lead ECG signal acquisition and analysis.
- 14- Mating connectors with force without proper alignment may damage the connectors and disrupt the electrical continuity.
- 15- Before using the EKG cable, read the instruction manual of the EKG device for further information and warnings.

Manual Cleaning and Disinfection:

Manual Cleaning: The EKG cables must be cleaned and disinfected immediately after each use. Remove heavy debris from the EKG cable by wiping it with an appropriate, moist, lint-free cloth, sponge, or gauze pad. Prepare the detergent solution with reverse osmosis or distilled water according to the manufacturer's instructions. Saturate a synthetic felt or sponge pad with the cleaning solution. Wipe all surfaces of the EKG cable until all visible contamination is removed. Wash the pad in the cleaning solution and repeat wiping the cable. Use a plastic cleaning brush for adhered particles or dried

secretions. Wipe or brush all surfaces at least six strokes back and forth. At least allow a contact time of 2 minutes. If visible contamination is still present on the cable at the end of the cleaning repeat the cleaning process. Moist a clean cloth, sponge, or gauze pad with tap water or reverse osmosis and thoroughly wipe the EKG cable to remove any debris and all traces of the cleaning solution. Dry manually by using a lint-free cloth or absorbent paper and, in particular, for drying cavities, sterile compressed air.

Manual Disinfection: Prepare the disinfectant solution according to the manufacturer's instructions. Only disinfect clean products that have dried slightly. Saturate a synthetic felt or sponge pad with the disinfectant solution. Wipe all surfaces of the EKG cable. Wash the pad in the disinfectant solution and repeat wiping the cable. Wipe or brush all surfaces at least six strokes back and forth. At least allow a contact time of 6 minutes. Moist a clean cloth, sponge, or gauze pad with reverse osmosis or distilled water and thoroughly wipe the EKG cable to remove any residual debris and all traces of disinfectant solution. Dry manually by using a lint-free cloth or absorbent paper and, in particular, for drying cavities, sterile compressed air.

EKG cable, especially metallic parts of the plugs and connector, should never be immersed in cleaning, disinfecting, or rinsing solutions or placed in an ultrasonic cleaning device. Soaking will reduce the life of the EKG cable by oxidation of metal parts and hardening of the cable jacket. Avoid using excessive force when cleaning or disinfecting by hand, which may damage the cable's internal wires and reduce the product's life. Please do not clean with hard or frictional materials, abrasive cleaners, or solvents. Cleaning agents containing phenol or chloride are not suitable. The EKG cable should be completely dry before use, do not use wet EKG cables. These instructions were validated using the below detergents and disinfectants.

Cleaning agent (manual cleaning): Enzol® / Cidezyme®, enzymatic detergent, Johnson&Johnson (add 8ml enzymatic detergent per 1 liter distilled water, for cable with dried-on organic matter use 16 ml enzymatic detergent per 1 liter distilled water/ warm water)
Disinfectant (manual disinfection): Cidex OPA®, 0.55% ortho-phthalaldehyde solution, Johnson&Johnson (at a minimum temperature of 20°C/68°F).
Reverse osmosis/ distilled water.

Caution: During the cleaning and disinfection of EKG cables, disposable gloves, protective eyewear, and a filtration mask should be used to reduce the risk of transmission of infectious agents by splashing or creating aerosols.

Warning: FMT EKG cables are sold NON-STERILE. Clean and disinfect the cables before the first and every use.

Warning: Before cleaning or disinfecting the cable, disconnect it from the EKG device.

Warning: EKG cable should be cleaned of overburden prior to disinfection to improve the effectiveness.

Warning: Do not sterilize EKG cables and leadwires by autoclave, radiation, or steam.

Warning: Avoid contact with strong, aromatic, chlorinated, ketone, ether, or ester solvents. Do not use bleach on electrical contacts or connectors.

Environmental Requirements:

Operating Temperature : 5°C - 45°C (41°F - 113°F)

Storage Temperature : 0°C - 50°C (32°F - 122°F)

Relative Humidity : 20 - 80% (not condensing)

The product should be protected from sunlight during storage. It is recommended that it be stored in its original package until its first use. Used cables should be loosely coiled in 10 to 15 cm diameter (4 to 6 inches) and must be stored individually in a protective container or packaging.

Compatibility: In order to ensure compatibility, the EKG / ECG cable should only be used with the equipment they have been designed for and is specified in a label on the product's primary packaging. For product models and compatible devices, please refer to the catalog at metkomedical.com/PDF/01DiagnosticEKGCablesandAccessories.pdf

Safety:

Degree of protection from electric shocks: type CF

Degree of protection against the ingress of water: Keep dry

FMT EKG / ECG cables are sold NON-STERILE.

FMT EKG / ECG cables are not made with natural rubber latex and PVC.

FMT reusable EKG / ECG cables have not been evaluated for safety in the MR environment. Scanning a patient who has this device may result in patient injury.

For an explanation of symbols, refer to the separate leaflet "Description of Symbols" located within the product package.

Reporting of Serious Incidents: Any serious incident related to the use of this product should be reported to both the manufacturer and the health authority/competent authority where the product is installed. Either contact your local representative or report to: metko@metkomedical.com

Caution: Federal Law (U.S.A.) restricts this device to sale by or on the order of a physician.

Warranty: FMT EKG cables are under six (6) months warranty against material and workmanship defects from the original purchase date. In the warranty period, METKO will be responsible for repairing or changing the cable free of charge if the defect is proven. This warranty does not extend to any product that has been subject to misuse, wrong disinfection, neglect, or accident, or that has been damaged by causes external to the product or that has been used in violation of the operating instructions supplied with the product. The product life is 3 years from the production date.

FMT® is a registered trademark of Metko A.Ş.

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