Indications: Reusable pulse oximetry (SpO₂) sensors are indicated for the continuous noninvasive monitoring of functional oxygen saturation of arterial hemoglobin (SpO₂), pulse rate, and plethysmographic pulse wave of the patient during no motion conditions. Reusable SpO₂ sensors are intended to be used with appropriate pulse oximeters or patient monitors by trained and qualified medical personnel in hospitals and hospital-type facilities.

Contraindications: Reusable SpO₂ sensors are contraindicated for use on active patients or for prolonged periods of use. The sensor site must be inspected at least every four (4) hours or sooner, and if the circulatory condition or skin integrity is compromised, the sensor should be applied to a different site. Reusable SpO₂ sensors are also contraindicated if they are used on other than the specified patient size (patient weight) and application site.

Clinical Benefits: A wired LED and PD sensor enables continuous measurement of patient functional oxygen saturation.

Product Description:


Sensor Type	Sensor Housing Type	Patient Size (Body Weight)	Application Site
FMT-RAF/XXX (Reusable clip SpO ₂ sensor)	Solid type finger clip	> 40 kg Adult	Index or other finger
FMT-RAF(B)/XXX (Reusable clip SpO ₂ sensor)	Solid type finger clip	> 40 kg Adult	Index or other finger
FMT-RPF/XXX (Reusable clip SpO ₂ sensor)	Solid type finger clip	10 - 40 kg Pediatric	Index or other finger
FMT-REC/XXX (Reusable ear clip SpO ₂ sensor)	Solid type ear clip	> 40 kg Adult	Ear lobe
FMT-RAS/XXX (Reusable soft SpO ₂ sensor)	Soft type silicone	> 40 kg Adult	Index or other finger
FMT-RAS(B)/XXX (Reusable soft SpO ₂ sensor)	Soft type silicone	> 40 kg Adult	Index or other finger
FMT-RPS/XXX (Reusable soft SpO ₂ sensor)	Soft type silicone	10 - 40 kg Pediatric	Index or other finger
FMT-RIS/XXX (Reusable soft SpO ₂ sensor)	Soft type silicone	3 - 15 kg Infant	Index or other finger
FMT-RWS/XXX (Reusable Wrap Y SpO ₂ Sensor)	Soft Type Silicone	> 40 kg Adult	Index or other finger
		10 - 40 kg Pediatric	Index or other finger
		3 - 15 kg Infant	Great toe
		1 - 4 kg Neonatal	Foot
		1 - 4 kg Neonatal	Foot or hand palm
FMT-RYS/XXX (Reusable Y SpO ₂ Sensor)	Soft Type Silicone	> 40 kg Adult	Index or other finger
		> 40 kg Adult	Earlobe
		10 - 40 kg Pediatric	Index or other finger
		3 - 15 kg Infant	Great toe
		1 - 4 kg Neonatal	Foot or hand palm
FMT-RNS/XXX (Reusable Wrap Soft SpO ₂ Sensor)	Soft Type Silicone	1 - 4 kg Neonatal	Foot or hand palm

XXX is the compatible device code. See the compatibility section.

Comparative value measurement in % saturation: Between 70% and 100% SpO₂ levels, SpO₂ accuracy (Functional Oxygen Saturation Arms) is $\pm 2\%$, and SpO₂ accuracy is not specified under the 70% SpO₂ levels. Pulse Rate accuracy (Arms) from 35 to 240 bpm is ± 2 bpm.


Directions for Use:

FMT-RAF/XXX, FMT-RPF/XXX, FMT-RAS/XXX, FMT-RPS/XXX and FMT-RIS/XXX sensor application: Select a site that least restricts a conscious patient's movements. The patient's index finger of the non-dominant hand is the preferred location. Alternatively, the middle finger, ring finger, or little finger of the hand may be used. Always choose a site that will completely cover the sensor's detector window. Thumb, the great toe, or the long toe (next to the great toe) may be used as secondary alternatives. Fit the sensor as illustrated. The patient's finger must be inserted until the tip of the finger reaches the finger stop. If the fingernail is long, it may extend over and pass the finger stop. Direct the cable along the patient's finger and parallel to the arm. Connect the sensor cable to the monitor or to the adaptor & extension cable and verify proper operation as described in the oximeter's operator manual. To disconnect, hold the connector and pull straight out.

FMT-RYS/XXX sensor application (as ear clip): Open the ear clip by pressing at the ends with the thumb and index finger. Slide each sensor head into ear clip slots with both optical component sides facing inwards. The appropriate sensor sites are either the ear lobe (a) or the ear pinna (b). For fast response, rub the application site in order to increase the circulation before applying the sensor. Clip the sensor onto the ear so that the sensor cable runs down the side of the patient's face and body. To remove the ear clip and sensor from the patient's ear, press open the ear clip ends. Remove the sensor from the ear clip by gently sliding the sensor head out of the ear clip slot. **Note:** If the sensor is not positioned properly, light may bypass the ear tissue and result in SpO₂ inaccuracies. If measurement cannot be obtained try another site on the same ear or on the other ear. If measurement still cannot be obtained, consider using the RYS sensor on another body site with adhesive tape instead of an ear clip.

FMT-REC/XXX sensor application: Rub the application site on the ear in order to increase the circulation before applying the sensor. Apply the sensor cable to the soft silicone fixing mechanism as shown in the picture. Open the ear clip by pressing at the ends with the thumb and index finger. The appropriate sensor site is the ear lobe (a). Clip sensor onto fleshy part of the lobe. Do not position the sensor on cartilage or press against the head. Hang the soft silicone fixing mechanism (ear hanger) on top of the ear. The fixing mechanism minimizes artifacts caused by patient motion. Remove the ear clip sensor from the patient's ear by pressing from the ends.

FMT-RYS/XXX and FMT-RWS/XXX sensor application with medical tape: Cut the adhesive medical tape with the appropriate length for the selected sensor site. Position the back of the sensor heads on the adhesive tape. Both optical component sides should face outwards. The preferred sensor sites are the index finger for adult and pediatric patients, the great toe or thumb for infants, and the fleshy part of the foot or the hand palm for Neonatals. Ensure that sensor heads are directly opposite each other through the tissue. Once positioned on the patient, wrap the tape smoothly. Do not wrap the tape around the limb so tightly that normal blood flow may be restricted. Secure the "Y" junction to the patient with additional tape. Be sure to leave slack in the wires between the sensor heads and the "Y" junction. When removing the sensor, first gently remove the tape from the patient's skin. Then separate the adhesive tape from the sensor heads by peeling it away.

FMT-RNS/XXX sensor application: The preferred sensor sites are fleshy part of the foot or hand palm for Neonatal patients. Apply the sensor on the patient. Position the inner face of the "C" shaped sensor on the limb. Insert the tip of the belt into first slot on the back of the sensor and draw out the belt. Then put the belt tip through the second slot to fasten the sensor. Do not pull the belt too much since this may restrict the normal blood flow. Ensure that optical components of the sensor are directly opposite each other through the tissue. When removing the sensor, first remove the belt from second slot then loosen the sensor.

Warnings:

- All SpO₂ sensors and adapter cables are designed for use with specific patient monitors or pulse oximeters. Verify the compatibility of the monitor, adapter cable, and sensor and all functions including alarming perform properly before use, otherwise degraded performance and/or patient injury can result.
- The application of reusable SpO₂ sensors should be done by trained medical personnel following established procedures.
- Clean or disinfect the sensors prior to reuse on multiple patients.
- All wire-lead patient-connected transducer assemblies are subject to reading error, local heating, and possible damage from high-intensity sources of RF energy. Electrosurgical equipment represents one such source since capacitively-coupled currents may seek alternative paths to ground through sensor cables and isolated instruments. Patient burns may result. If possible, remove the sensor from patient contact before activating the surgical unit or other RF source. If sensors must be used simultaneously with electrosurgical apparatus, the instruments to which the sensors are connected should be checked for adequate isolation from electrical grounds at radio frequencies. Hazards can be reduced by selecting a monitoring point located away from the expected RF current path between the active probe and the return pad.
- The sensor should be free of visible defects, discoloration and damage. If the sensor is wet, discolored or damaged, discontinue use. Never use a damaged sensor or one with exposed electrical circuitry.
- The SpO₂ sensor is in direct contact with the patient's skin and therefore allergic reactions may occur in susceptible individuals. In case of allergic or skin reactions discontinue use and consult a specialist.
- All SpO₂ sensors may be used on the same site for a maximum of four (4) hours, provided the site is inspected routinely per clinical protocol to ensure skin integrity, circulatory status, and correct optical alignment. Skin erosion, tissue ischemia, and/or pressure necrosis can be caused when the sensor is not frequently moved, applied too tightly, or becomes too tight due to edema. Assess the site as frequently as every (1) hour and move the sensor or use a different type of sensor if there are signs of pyrexia or loss of skin integrity and/ or loss of circulation or perfusion.
- As with all medical equipment, carefully route cables to reduce the possibility of patient entanglement or strangulation.
- Venous congestion may cause under reading of actual arterial oxygen saturation. Therefore, assure proper venous outflow from the monitored site. The sensor should not be below heart level (e.g. sensor on the hand of a patient in a bed with an arm dangling to the floor).
- Non-venous pulsations may cause inaccurate readings (e.g. tricuspid valve regurgitation, Trendelenburg position).
- Do not use the sensor during Magnetic Resonance Imaging (MRI) scanning or in an MRI environment as it may result in physical harm.
- Misapplied sensors or sensors that become partially dislodged may cause incorrect readings.
- Inaccurate readings may be caused by EMI radiation interference or strong electromagnetic sources such as electrosurgery equipment.
- Abnormal fingers, intravascular dyes such as indocyanine green or methylene blue, or externally applied coloring and texture such as fingernail polish, or artificial fingernails may lead to inaccurate or no readings.
- The application site should be cleaned of debris prior to sensor placement. Do not apply the sensor to the injured skin.
- Elevated levels of Carboxyhemoglobin (COHb), Methemoglobin (MetHb), and total bilirubin may lead to inaccurate SpO₂ readings. High levels of COHb or MetHb may occur with a seemingly normal SpO₂. When elevated levels of COHb or MetHb are suspected, laboratory analysis (CO-Oximetry) of a blood sample should be performed.
- Inaccurate SpO₂ readings may be caused by severe anemia, very low arterial perfusion, vasospastic disease such as Raynaud's, and peripheral vascular disease at the monitored site or extreme motion artifact.
- Inaccurate SpO₂ readings may be caused by elevated levels of dysmethemoglobin, hypoxic or hypercapnic conditions, and severe vasoconstriction or hypothermia.
- Hemoglobinopathies and synthesis disorders such as thalassemias, Hb s, Hb c, sickle cell, etc. may cause inaccurate SpO₂ readings.

- 20- If using pulse oximetry during full body irradiation, keep the sensor out of the radiation field. If the sensor is exposed to the radiation, the reading might be inaccurate or not provided or the unit might read zero for the duration of the active radiation period.
- 21- The pulsations from intra-aortic balloon support can be additive to the pulse rate on the oximeter pulse rate display. Verify pulse rate against the ECG heart rate.
- 22- Sensors should not compress the limb. A tight sensor will compromise venous return and may give rise to inaccurate readings. When selecting a sensor type it is important to check if the patient's weight is in the recommended patient weight range.
- 23- Sensors that are unapplied or fallen due to the loose application may cause readings. To avoid misdiagnosis, verify the sensor match the patient weight and the sensor is applied to the patient correctly.
- 24- High ambient light sources such as surgical lights (especially those with a xenon light source), bilirubin lamps, fluorescent lights, infrared heating lamps, and direct sunlight can interfere with the performance of the sensor. To prevent interference from ambient light, ensure that the sensor is properly applied, and cover the sensor site with opaque material if required. Failure to take this precaution in high ambient light conditions may result in inaccurate readings.
- 25- At elevated ambient temperatures, the patient's skin could be severely burned after prolonged sensor application at sites that are not well perfused. To prevent this condition, be sure to check patient application sites frequently. All listed sensors operate without risk of exceeding 41°C on the skin if the initial skin temperature does not exceed 35°C.
- 26- When selecting the sensor site, priority should be given to an extremity free of an arterial catheter, blood pressure cuff, or intravascular infusion line.
- 27- Do not use tape to secure the sensor to the site; this can restrict blood flow and cause inaccurate readings. The use of tape can cause skin damage, and/or pressure necrosis or damage to the sensor.
- 28- The performance of the sensor is compromised by motion. Finger clip type sensor is contraindicated for use on active (mobile) patients. Soft type silicone sensors tolerate a moderate amount of patient activity. Use soft type silicone or disposable sensors for active patients.
- 29- Do not twist unnecessarily or use excessive force when using, connecting, disconnecting, or storing the sensor. Do not pull on the cable.
- 30- High oxygen concentration may predispose a premature infant to retinopathy. Therefore, the upper alarm limit for the oxygen saturation must be carefully selected in accordance with accepted clinical standards. The accuracy range of the oximeter being used must be considered.
- 31- Portable and mobile RF communications equipment can affect medical electrical equipment and should be used no closer to any part of the equipment.
- 32- The sensor is not protected against the effect of the discharge of a cardiac defibrillator. Remove the sensor before cardiac defibrillation.
- 33- To prevent damage, do not expose connector pins to cleaning solutions, and do not soak or immerse the sensor in any liquid solution.
- 34- Finger clip type sensors are not intended for use on the thumb or across an infant or Neonatal patient's hand or foot.
- 35- Keep the ear clip of the Y sensor out of the reach of children under the age of 3, it is small enough to be swallowed and may block the trachea.
- 36- Discard the ear clip sensor if loss of spring tension allows for slippage or movement of the ear clip from its proper position on the ear lobe or pinna. That may result in inaccurate readings.
- 37- Do not use the Y sensor with ear clip or ear clip sensor on any site other than those on the ear. This may result in inaccurate readings.
- 38- Remove earrings from the patient's ear before applying the ear clip sensor. Ear clip sensors are not recommended for pediatric patients with less than 40kg or for neonatal use.
- 39- When using wrap Neonatal sensors, do not pull the belt too tight as normal blood circulation may be affected.
- 40- Keep the sensor dry. Do not use a wet sensor. This may cause burns during the application of high-frequency devices.
- 41- Never attempt to perform repairs. Do not modify the SpO2 sensors in any way. Modifications adversely affect patient safety, performance, and/or accuracy, and void the warranty.
- 42- Do not attempt to sterilize by irradiation, steam, autoclave or ethylene oxide as it will damage the sensor.
- 43- Dispose damaged sensors according to local laws and regulations for medical waste.
- 44- Refer to the monitor's operation manual for additional cautions and warnings. For operation errors and screen error messages, refer to the trouble shooting table provided in the operation manual of the monitor.

Manual Cleaning and Disinfection:

Manual Cleaning: The reusable SpO2 sensors must be cleaned and disinfected immediately after each use. Remove heavy debris from the SpO2 sensor by wiping with an appropriate, moist lint-free cloth, sponge, or gauze pad. Prepare the detergent solution with reverse osmosis or distilled water according to the manufacturer's instructions. Saturate a synthetic felt or sponge pad with the cleaning solution. Wipe all surfaces of the SpO2 sensor until all visible contamination is removed. Wash the pad in the cleaning solution and repeat wiping the sensor. Use a plastic cleaning brush for adhered particles or dried secretions. Wipe or brush all surfaces at least six strokes back and forth. At least allow a contact time of 2 minutes. If visible contamination is still present on the sensor at the end of the cleaning repeat the cleaning process. Moist a clean cloth, sponge, or gauze pad with tap water or reverse osmosis and thoroughly wipe the SpO2 sensor to remove any debris and all traces of cleaning solution. Dry manually by using a lint free cloth or absorbent paper and, in particular, for drying cavities, sterile compressed air.

Manual Disinfection: Prepare the disinfectant solution according to the manufacturer's instructions. Only disinfect clean products that have dried slightly. Saturate a synthetic felt or sponge pad with the disinfectant solution. Wipe all surfaces of the SpO2 sensor. Wash the pad in the disinfectant solution and repeat wiping the cable. Wipe or brush all surfaces at least six strokes back and forth. At least allow a contact time of 6 minutes. Moist a clean cloth, sponge, or gauze pad with reverse osmosis or distilled water and thoroughly wipe the SpO2 sensor to remove any residual debris and all traces of disinfectant solution. Dry manually by using a lint free cloth or absorbent paper and, in particular, for drying cavities, sterile compressed air.

SpO2 sensor, especially metallic parts of connector should never be immersed in cleaning, disinfecting, or rinsing solutions, or placed in an ultrasonic cleaning device. Soaking will reduce the life of SpO2 sensor by oxidation of metal parts and hardening of the cable jacket. Avoid using excessive force when cleaning or disinfecting by hand that may damage the internal wires of the cable and reduces the product life. Please do not clean with hard or frictional materials, abrasive cleaners, or solvents. Cleaning agents containing phenol or chloride are not suitable. The SpO2 sensor should be completely dry before use, do not use wet SpO2 sensors. These instructions were validated using the below detergents and disinfectants.

Cleaning agent (manual cleaning): Enzo® / Cidezyme®, enzymatic detergent, Johnson&Johnson (add 8ml enzymatic detergent per 1 liter distilled water, for cable with dried-on organic matter use 16 ml enzymatic detergent per 1 liter distilled water/ warm water)

Disinfectant (manual disinfection): Cidex OPA®, 0.55% ortho-phthalaldehyde solution, Johnson&Johnson (at a minimum temperature of 20°C/68°F).

Reverse osmosis/ distilled water.

Caution: During the cleaning and disinfection of SpO2 sensors disposable gloves, protective eyewear, and filtration mask should be used to reduce the risk of transmission of infectious agents by splashing or the creation of aerosols.

Warning: FMT reusable SpO2 sensors are sold NON-STERILE. Clean and disinfect the sensors before first and every use.

Warning: Before cleaning or disinfecting the sensor, disconnect it from the pulse oximeter.

Warning: SpO2 sensors should be cleaned of overburden prior to disinfection to improve the effectiveness.

Warning: Do not sterilize SpO2 sensors by autoclave, radiation or steam.

Warning: Avoid contact with strong, aromatic, chlorinated, ketone, ether, or ester solvents. Do not use bleach on electrical contacts or connectors.

Measurement Wavelengths and Output Power (*)

LED Type	Red Peak Wavelength	Red Max. Output Power	IR Peak Wavelength	IR Max. Output Power
2-Leads	663nm	@1.2mW nominal	890nm	@1.0mW nominal
3-Leads	661nm	@1.2mW nominal	940nm	@1.2mW nominal
4-Leads	660nm	@1.2mW nominal	905/940nm	@1.0mW nominal

*This information is especially useful for clinicians performing photodynamic therapy.

Environmental Requirements:

Operating Temperature : 5°C - 45°C (41°F - 113°F)

Storage Temperature : 0°C - 50°C (32°F - 122°F)

Relative Humidity : 20 - 80% (not condensing)

During the storage the product should be protected from sunlight. It is recommended to store the product in its original package until the first use. Used sensors should be loosely coiled in 10 to 15 cm diameter (4 to 6 inches) and must be stored individually in a protective container or packaging.

Compatibility: In order to ensure compatibility and claimed accuracy, the reusable SpO2 sensors should only be used with the equipment for which they have been designed and is specified in the label on the product and/or the product's primary packaging. The reusable SpO2 sensor is designed to be compatible with the patient monitors using the appropriate adapter and extension cable. For product models and compatible devices, please refer to the catalog at [metkomedical.com/PDF/03PulseOximetry\(SpO2\)Sensors.pdf](http://metkomedical.com/PDF/03PulseOximetry(SpO2)Sensors.pdf)

Safety:

Degree of protection from electric shocks: type BF

Degree of protection against the ingress of water: IPX2

The FMT reusable SpO2 sensors are sold NON-STERILE.

The FMT reusable SpO2 sensors are not made with natural rubber latex and PVC.

The FMT reusable SpO2 sensors have not been evaluated for safety in the MR environment. Scanning a patient who has this device may result in patient injury.

Refer to the separate leaflet "Description of Symbols" located within the product package for the explanation of symbols.

Note: Functional Tester cannot be used to assess the accuracy of a pulse oximeter sensor or a pulse oximeter monitor.

Reporting of serious incidents: Any serious incident related to the use of this product should be reported to both the manufacturer and the health authority/competent authority where the product is installed. Either contact your local representative or report to: metko@metkomedical.com

Caution: Federal Law (U.S.A.) restricts this device to sale by or on the order of a physician.

Warranty: FMT reusable SpO2 sensors are under twelve (12) months warranty against material and workmanship defects from the date of original purchase. In warranty period, METKO will be responsible for repairing or changing the sensor free of charge if the defect is proven. This warranty does not extend to any product that has been subject to misuse, neglect or accident; or that has been damaged by causes external to the product, or that has been used in violation of the operating instructions supplied with the product. Product life is 3 years.

FMT® is registered trademark of Metko A.Ş.

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