Reusable Esophageal / Rectal Temperature Probe

English

Indications:

The reusable esophageal/rectal temperature probe is intended for use in routine continuous monitoring of esophageal or rectal temperature as an indicator of the patient's core body temperature. Reusable esophageal/rectal temperature probes are intended to be used with appropriate electronic thermometers by trained and qualified medical personnel in hospitals and hospital-type facilities.

Contraindications:

The use of an esophageal probe may be contraindicated in patients with esophageal diverticulum or stenosis, in neonates or small infants undergoing neck surgery, and in patients undergoing a tracheostomy or insertion of an internal jugular catheter. Rectal placement is contraindicated in neonates and small infants due to the possibility of colon perforation.

Adverse Reactions:

Adverse reactions reported during applications of esophageal / rectal temperature probes include: Accidental tracheal or bronchial placement accompanied with an airway obstruction, esophageal abrasion and/or perforation, pharyngeal abrasion, rectal abrasion and Tissue burns due to aberrant electro-cautery current radio-frequency current pathways.

Clinical Benefits:

A wired temperature sensor enables continuous measurement of patient core (central) temperature.

Product Description:

The use of the temperature probe with incompatible electronic thermometers can result in degraded performance and accuracy. The operator is responsible for checking the compatibility of the monitor, the temperature probe, and the extension cable before use. Probes are accurate and interchangeable with any other probe with matching temperature/resistance characteristics at specified measurement temperatures. Within the measurement range of 25 - 45 °C, the accuracy for the probes is $\pm 0.1^{\circ}\text{C}$. The minimum measurement time for accurate temperature readings is 150 seconds. The temperature reaches a stabilized temperature in 5 minutes, depending on the body site. The system accuracy is the sum of the probe and instrument accuracies.

FMT400/AOR, FMT400/AOR/Z (3m), **FMT400/AOR-15, FMT400/AOR/Z-15** (1.5m), Esophageal / Rectal probe, Adult
FMT400/POR, FMT400/POR/Z (3m), **FMT400/POR-15, FMT400/POR/Z-15** (1.5m), Esophageal / Rectal probe, Pediatric
FMT400/IOR, FMT400/IOR/Z (3m), **FMT400/IOR-15, FMT400/IOR/Z-15** (1.5m), Esophageal / Rectal probe, Infant
FMT400/AOR-A, FMT400/AOR/Z-A (3m), **FMT400/AOR-A15, FMT400/AOR/Z-A15** (1.5m), Esophageal / Rectal probe, Adult, Autoclavable
FMT400/POR-A, FMT400/POR/Z-A (3m), **FMT400/POR-A15, FMT400/POR/Z-A15** (1.5m), Esophageal / Rectal probe, Pediatric, Autoclavable
(Z for Straight 6.35mm Phone Plug)

FMT400/AOR/XX (3m), **FMT400/AOR/XX-15** (1.5m), Esophageal / Rectal probe, Adult

FMT400/POR/XX (3m), **FMT400/POR/XX-15** (1.5m), Esophageal / Rectal probe, Pediatric

FMT400/IOR/XX (3m), **FMT400/IOR/XX-15** (1.5m), Esophageal / Rectal probe, Infant

FMT400/AOR/XX-A (3m), **FMT400/AOR/XX-A15** (1.5m), Esophageal / Rectal probe, Adult, Autoclavable

FMT400/POR/XX-A (3m), **FMT400/POR/XX-A15** (1.5m), Esophageal / Rectal probe, Pediatric, Autoclavable

(XX variable is for compatible monitor: EDN (Edan), HP (HP & Philips), MND (Mindray), S (Drager/Siemens), BLT (Biolight), GE (GE), SW (S&W), THT (Tecotherm), MMM (Belmont), GH (GE))

Directions for Use:

The reusable temperature probe and its plug should be inspected for visual defects. Probes with a loose or damaged plug or contact pins or corrosion of metal contacts, damage in the insulation of cable (i.e. cuts, cracks, brittle or solid, permanent bends, burnt or melted), wires or metal parts that have come out of the insulation of the cable or plug, or not completely dried and ventilated probes should not be used.

- 1- Before use, check that the product is intact and clean. Inspect the temperature probe cable for wear, breakage or fraying. Replace if necessary.
- 2- Make sure that the probe is compatible with the monitor to be connected and that the plug of the probe matches the socket on the monitor. To connect the probe to the temperature monitor, align the plug of the probe with the socket of the monitor or extension cable and push firmly to ensure full contact. Forced mating of connectors without proper alignment can result in damage to connectors and loss of electrical continuity.
- 3- Apply the probe to the rectum or the esophagus according to accepted clinical procedures. If used in the esophagus, remove the probe prior to intubation or extubation of the tracheal tube following accepted medical technique. Verify position of the probe by medical procedures.
- 4- Wait to see the complete measurement reading after the response time is reached. Follow the monitoring device manufacturer's instructions of use and your hospital's protocols on temperature monitoring of the patients.
- 5- When temperature monitoring is completed, first carefully remove the esophageal / rectal probe from the patient according to accepted medical techniques then disconnect the probe from the monitor. Hold the plastic plug firmly and pull. Do not pull on the cable.

Warnings:

- 1- All temperature probes are designed for use with specific patient monitors or electronic thermometers. Verify the compatibility of the monitor, extension cable, and probe and all functions including alarming perform properly before use, otherwise degraded performance, accuracy, and/or patient injury can result.
- 2- Application of temperature probes should be done by trained medical personnel following established procedures.
- 3- All wire-lead patient-connected transducer assemblies are subject to reading error, local heating, and possible damage from high-intensity sources of RF energy. Electrosurgical equipment represents one such source since capacitively-coupled currents may seek alternative paths to ground through probe cables and isolated instruments. Patient burns may result. If possible, remove the probe from patient contact before activating the surgical unit or other RF source. If probes must be used simultaneously with electrosurgical apparatus, the instruments to which the probes are connected should be checked for adequate isolation from electrical grounds at radio frequencies. Hazards can be reduced by selecting a temperature monitoring point located away from the expected RF current path between the active probe and the return pad.
- 4- It is not recommended to use the temperature probes within the range of Magnetic Resonance Tomography (MRI), Computed Tomography (CT) etc. applications. Conducted current may cause burns.
- 5- The temperature probe is in direct contact with the patient's skin and therefore allergic reactions may occur in susceptible individuals. In case of allergic or skin reactions discontinue use and consult a specialist.
- 6- As with all medical equipment, carefully route temperature probe cables to reduce the possibility of patient entanglement or strangulation.
- 7- Assure that probe and extension cable are not positioned where they could become entangled around IV tubes and restrict IV flow.
- 8- The medical thermometer cannot measure the patient's temperature accurately if the Esophageal / Rectal probe is not correctly positioned or probe tip is not in direct contact with the patient's skin. Regularly check the skin probe remain correctly positioned on the patient and probe is always in direct contact with the patient's skin
- 9- Use only water or medical-grade water based lubricants to lighten insertion of the esophageal / rectal probe. Do not use petroleum based lubricants
- 10- When used with external heating and cooling systems, do not position the probe that can be affected by the heating or cooling source.
- 11- If the esophageal/rectal probe is to be used esophageal, intubate the patient with the endotracheal tube before inserting the probe. Inform the surgeon that a probe has been inserted. Remove the probe prior to extubation of the tracheal tube following the accepted medical technique.
- 12- Do not administer the esophageal/rectal probe in the tracheal or bronchial way as it may cause airway obstruction.
- 13- After long term use the jacket of the Esophageal / Rectal probe may lose flexibility due to effect of disinfection agents. In medical applications, the user must determine that a probe is suitable and sufficient flexible for esophageal or rectal use.
- 14- Always use caution when applying, inserting, or removing a temperature probe from a patient.
- 15- Inspect and clean or disinfect the probe before attaching to a new patient.
- 16- Do not use wet probes. This may cause burns during application of high frequency devices.
- 17- Portable and mobile RF communications equipment can affect medical electrical equipment and should be used no closer to any part of the equipment.
- 18- The probe is not protected against the effect of the discharge of a cardiac defibrillator. Remove the probe before cardiac defibrillation.
- 19- Mishandling of the probes could result in damage to internal wires and loss of electrical isolation or improper temperature readings. Don't tie the probes in knots or bunch tightly. Don't wrap probes around equipment to avoid damaging internal wires. Don't leave the probes dangling from the instrument. When not in use, probes and leads should be loosely coiled and stored in room temperature.
- 20- Care must be taken to ensure that the probe cable is not subjected to mechanical stress, such as pulling, bending and squeezing.
- 21- Never boil temperature probes or clean with ultrasonic cleaner. Do not immerse the probe in liquid.
- 22- Never attempt to perform repairs. Do not modify the temperature probes in any way. Modifications may affect patient safety, performance, and void the warranty.
- 23- Do not use damaged probes. Dispose damaged probes according to local laws and regulations for medical waste.
- 24- For further information and warnings read the instructions of use accompanying the temperature monitoring device.

Manual Cleaning and Disinfection:

Manual Cleaning:

The temperature probes must be cleaned and disinfected immediately after each use. Remove heavy debris from the temperature probe by wiping with an appropriate, moist lint-free cloth, sponge, or gauze pad. Prepare the detergent solution with reverse osmosis or distilled water according to the manufacturer's instructions. Saturate a synthetic felt or sponge pad with the cleaning solution. Wipe all surfaces of the temperature probe cable until all visible contamination is removed. Wash the pad in the cleaning solution and repeat wiping the cable. Use a plastic cleaning brush for adhered particles or dried secretions. Wipe or brush all surfaces at least six strokes back and forth. At least allow a contact time of 2 minutes. If visible contamination is still present on the probe at the end of the cleaning repeat the cleaning process. Moisten a clean cloth, sponge, or gauze pad with tap water or reverse osmosis and thoroughly wipe the temperature probe to remove any debris and all traces of cleaning solution. Dry manually by using a lint free cloth or absorbent paper and, in particular, for drying cavities, sterile compressed air.

Manual Disinfection:

Prepare the disinfectant solution according to the manufacturer's instructions. Only disinfect clean products that have dried slightly. Saturate a synthetic felt or sponge pad with the disinfectant solution. Wipe all surfaces of the temperature probe. Wash the pad in the disinfectant solution and repeat wiping the probe cable. Wipe or brush all surfaces at least six strokes back and forth. At least allow a contact time of 6 minutes. Moist a clean cloth, sponge, or gauze pad with reverse osmosis or distilled water and thoroughly wipe the temperature probe to remove any residual debris and all traces of disinfectant solution. Dry manually by using a lint free cloth or absorbent paper and, in particular, for drying cavities, sterile compressed air.

Temperature probe, especially metallic parts of connector should never be immersed in cleaning, disinfecting, or rinsing solutions, or placed in an ultrasonic cleaning device. Soaking will reduce the life of temperature probe by oxidation of metal parts and hardening of the cable jacket. Avoid using excessive force when cleaning or disinfecting by hand that may damage the internal wires of the probe and reduces the product life. Please do not clean with hard or frictional materials, abrasive cleaners, or solvents. Cleaning agents containing phenol or chloride are not suitable. The temperature probe should be completely dry before use, do not use wet temperature probes. These instructions were validated using the below detergents and disinfectants.

Cleaning agent (manual cleaning): Enzol® / Cidezyme®, enzymatic detergent, Johnson&Johnson (add 8ml enzymatic detergent per 1 liter distilled water, for cable with dried-on organic matter use 16 ml enzymatic detergent per 1 liter distilled water/ warm water).

Disinfectant (manual disinfection): Cidex OPA®, 0.55% ortho-phthalaldehyde solution, Johnson&Johnson (at a minimum temperature of 20°C/68°F).
Reverse osmosis/ distilled water.

Caution: During the cleaning and disinfection of temperature probes disposable gloves, protective eyewear, and filtration mask should be used to reduce the risk of transmission of infectious agents by splashing or the creation of aerosols.

Warning: FMT temperature probes are sold NON STERILE. Clean and disinfect the probes before first and every use.

Warning: Before cleaning or disinfecting the probes, disconnect it from the medical thermometer.

Warning: Temperature probes should be cleaned of overburden prior to disinfection to improve the effectiveness.

Warning: Do not sterilize temperature probes by autoclave (except blue probes), radiation or steam.

Warning: Avoid contact with strong, aromatic, chlorinated, ketone, ether, or ester solvents. Do not use bleach on electrical contacts or connectors.

Sterilization:

The temperature probes must be cleaned and disinfected prior to sterilization. When preparing the probe for sterilization, do not coil the cable into a circle of less than 10cm (4") in diameter. Pack the probe in disposable sterilization packaging (single or double packaging) made from paper/polyethylene and/or put it in a sterilization container. Packaging must comply with the ISO 11607-1 and EN 868-5 standards for packaging for sterilized instruments. Sterilization of the probes should be done by steam sterilization using the fractionated pre vacuumed method. The probes must be sterilized at a minimum of 134°C [273°F] and a maximum of 137°C [278°F] in saturated steam during a holding (exposure) time of at least 4 minutes to at most 20 minutes, then dried in vacuum for at least 20 minutes. When sterilizing several products at the same time in a steam autoclave, ensure that the maximum load capacity of the autoclave specified by the manufacturer is not exceeded. Sterilization must be performed in accordance with EN ISO 17665 (Sterilization of Medical Devices in Moist Heat).

The above sterilization process has been validated to assure a Sterility Assurance Level of 10^{-6} ($SAL \leq 1:10^6$) and the proof of the effectiveness of steam sterilization has been provided by an independent, accredited test laboratory (test report number AYA-SVR100620-01 dated 10.06.2020). The use of other sterilization methods (e.g. ethylene oxide gas, formaldehyde, gamma radiation, or low-temperature plasma sterilization) or other autoclave sterilization cycles is outside the manufacturer's responsibility.

Note: Before use, the temperature probes must be cooled to room temperature. It is the sole responsibility of the user to safeguard and maintain the sterile condition of the probes after the sterilization process.

Warning: Only blue colored temperature probes can be autoclaved. Do not autoclave grey colored non-autoclavable standard temperature probes. Such an attempt will destroy the probe.

Warning: Do not sterilize with hot air. Do not use "Flash" autoclave procedures.

Environmental Requirements:

Operating Temperature : 5°C - 45°C (41°F - 113°F)

Storage Temperature : 0°C - 50°C (32°F - 122°F)

Relative Humidity : 20 - 80% (not condensing)

During the storage the product should be protected from sunlight. It is recommended to store the product in its original package until the first use. Used probes should be loosely coiled in 10 to 15 cm diameter (4 to 6 inches) and must be stored individually in a protective container or sterile packaging.

Compatibility:

In order to ensure compatibility and claimed accuracy, the temperature probe should only be used with the equipment for which they have been designed and is specified in the label on the product's primary packaging. All the temperature probes except EDAN compatible EDN coded probes are compatible with patient monitors that use 400-Series technology. For product models and compatible devices, please refer to the catalog at metkomedical.com/PDF/04MedicalTemperatureProbes.pdf

Safety:

Degree of protection from electric shocks: type BF

Degree of protection against the ingress of water: IPX2

The FMT reusable temperature probes are sold NON-STERILE.

The FMT reusable temperature probes are not made with natural rubber latex and PVC.

The FMT reusable temperature probes have not been evaluated for safety in the MR environment. Scanning a patient who has this device may result in patient injury.

Refer to the separate leaflet "Description of Symbols" located within the product package for the explanation of symbols.

Reporting of serious incidents:

Any serious incident related to the use of this product should be reported to both the manufacturer and the health authority/competent authority where the product is installed.

Either contact your local representative or report to: metko@metkomedical.com

Caution: Federal Law (U.S.A.) restricts this device to sale by or on the order of a physician.

Warranty:

FMT reusable standard non-autoclavable temperature probes are under twelve (12) months and autoclavable temperature probes are under six (6) months warranty against material and workmanship defects from the date of original purchase. Autoclavable reusable temperature probes designed to withstand a minimum of 50 sterilization cycles when properly cared for and sterilized according to given instructions. In warranty period, METKO will be responsible for repairing or changing the probe free of charge if the defect is proven. This warranty does not extend to any product that has been subject to misuse, wrong sterilization, neglect or accident; or that has been damaged by causes external to the product; or that has been used in violation of the operating instructions supplied with the product. The product life is 3 years of the production date.

FMT® is a registered trademark of Metko A.Ş.

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FMT400/AOR, FMT400/AOR/Z (3m), **FMT400/AOR-15, FMT400/AOR/Z-15** (1.5m), Özofagual / Rektal prob, Yetişkin
FMT400/POR, FMT400/POR/Z (3m), **FMT400/POR-15, FMT400/POR/Z-15** (1.5m), Özofagual / Rektal prob, Pediatrik
FMT400/OR, FMT400/OR/Z (3m), **FMT400/OR-15, FMT400/OR/Z-15** (1.5m), Özofagual / Rektal prob, İnfant
FMT400/AOR-A, FMT400/AOR/Z-A (3m), **FMT400/AOR-A15, FMT400/AOR/Z-A15** (1.5m), Özofagual / Rektal prob, Yetişkin, Otoklavlanabilir
FMT400/POR-A, FMT400/POR/Z-A (3m), **FMT400/POR-A15, FMT400/POR/Z-A15** (1.5m), Özofagual / Rektal prob, Pediatrik, Otoklavlanabilir
(Düz 6.35mm Telefon Fişi için /Z)
FMT400/AOR/XX (3m), **FMT400/AOR/XX-15** (1.5m), Özofagual / Rektal prob, Yetişkin
FMT400/POR/XX (3m), **FMT400/POR/XX-15** (1.5m), Özofagual / Rektal prob, Pediatrik
FMT400/OR/XX (3m), **FMT400/OR/XX-15** (1.5m), Özofagual / Rektal prob, İnfant
FMT400/AOR/XX-A (3m), **FMT400/AOR/XX-A15** (1.5m), Özofagual / Rektal prob, Yetişkin, Otoklavlanabilir
FMT400/POR/XX-A (3m), **FMT400/POR/XX-A15** (1.5m), Özofagual / Rektal prob, Pediatrik, Otoklavlanabilir
(XX değişkeni uyumlu monitör içindir: EDN (Edan), HP (HP & Philips), MND (Mindray), S (Dräger/Siemens), BLT (Biolight), GE (GE), SW (S&W), THT (Tecotherm), MMM (Belmont), GH (GE))

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FMT400/AOR, FMT400/AOR/Z (3m), **FMT400/AOR-15, FMT400/AOR/Z-15** (1.5m), Sonde œsophagienne / rectale, Adulte
FMT400/POR, FMT400/POR/Z (3m), **FMT400/POR-15, FMT400/POR/Z-15** (1.5m), Sonde œsophagienne / rectale, Pédiatrique
FMT400/IOR, FMT400/IOR/Z (3m), **FMT400/IOR-15, FMT400/IOR/Z-15** (1.5m), Sonde œsophagienne/rectale, Nourrisson
FMT400/AOR-A, FMT400/AOR/Z-A (3m), **FMT400/AOR-A15, FMT400/AOR/Z-A15** (1.5m), Sonde œsophagienne / rectale, Adulte, Autoclavable
FMT400/POR-A, FMT400/POR/Z-A (3m), **FMT400/POR-A15, FMT400/POR/Z-A15** (1.5m), Sonde œsophagienne / rectale, Pédiatrique, Autoclavable
(Z pour prise téléphonique droite de 6.35 mm)

FMT400/AOR/XX (3m), **FMT400/AOR/XX-15** (1.5m), Sonde œsophagienne / rectale, Adulte
FMT400/POR/XX (3m), **FMT400/POR/XX-15** (1.5m), Sonde œsophagienne / rectale, Pédiatrique

FMT400/IOR/XX (3m), **FMT400/IOR/XX-15** (1.5m), Sonde œsophagienne/rectale, Nourrisson

FMT400/AOR/XX-A (3m), **FMT400/AOR/XX-A15** (1.5m), Sonde œsophagienne / rectale, Adulte, Autoclavable

FMT400/POR/XX-A (3m), **FMT400/POR/XX-A15** (1.5m), Sonde œsophagienne / rectale, Pédiatrique, Autoclavable

(La variable **XX** est pour un moniteur compatible: EDN (Edan), HP (HP & Philips), MND (Mindray), S (Drager /Siemens), BLT (Biolight), GE (GE), SW (S&W), THT (TecoTherm), MMM (Belmont), GH (GE))

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FMT400/AOR, FMT400/AOR/Z (3m), FMT400/AOR-15, FMT400/AOR/Z-15 (1.5m), Sonda Esofagea/Rettale, Adulto
FMT400/POR, FMT400/POR/Z (3m), FMT400/POR-15, FMT400/POR/Z-15 (1.5m), Sonda Esofagea/Rettale, Pediatrica
FMT400/IOR, FMT400/IOR/Z (3m), FMT400/IOR-15, FMT400/IOR/Z-15 (1.5m), Sonda Esofagea/Rettale, Neonato
FMT400/AOR-A, FMT400/AOR/Z-A (3m), FMT400/AOR-A15, FMT400/AOR/Z-A15 (1.5m), Sonda Esofagea/Rettale, Adulto, Autoclavabile
FMT400/POR-A, FMT400/POR/Z-A (3m), FMT400/POR-A15, FMT400/POR/Z-A15 (1.5m), Sonda Esofagea/Rettale, Pediatrica, Autoclavabile
(Z per spina telefonica dritta da 6.35 mm)
FMT400/AOR/XX (3m), FMT400/AOR/XX-15 (1.5m), Sonda Esofagea/Rettale, Adulto
FMT400/POR/XX (3m), FMT400/POR/XX-15 (1.5m), Sonda Esofagea/Rettale, Pediatrica
FMT400/IOR/XX (3m), FMT400/IOR/XX-15 (1.5m), Sonda Esofagea/Rettale, Neonato
FMT400/AOR/XX-A (3m), FMT400/AOR/XX-A15 (1.5m), Sonda Esofagea/Rettale, Adulto, Autoclavabile
FMT400/POR/XX-A (3m), FMT400/POR/XX-A15 (1.5m), Sonda Esofagea/Rettale, Pediatrica, Autoclavabile
(XX variabile è per monitor compatibile: EDN (Edan), HP (HP & Philips), MND (Mindray), S (Drager/Siemens), BLT (Biolight), GE (GE), SW (S&W), THT (Tecotherm), MMM (Belmont), GH (GE))

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FMT400/AOR, FMT400/AOR/Z (3m), FMT400/AOR-15, FMT400/AOR/Z-15 (1.5m), sonda esofágica/rectal, adulto
FMT400/POR, FMT400/POR/Z (3m), FMT400/POR-15, FMT400/POR/Z-15 (1.5m), sonda esofágica/rectal, pediátrica
FMT400/IOR, FMT400/IOR/Z (3m), FMT400/IOR-15, FMT400/IOR/Z-15 (1.5m), sonda esofágica/rectal, lactante
FMT400/AOR-A, FMT400/AOR-Z-A (3m), FMT400/AOR-A15, FMT400/AOR-Z-A15 (1.5m), sonda esofágica/rectal, adulto, esterilizable en autoclave
FMT400/POR-A, FMT400/POR-Z-A (3m), FMT400/POR-A15, FMT400/POR-Z-A15 (1.5m), sonda esofágica/rectal, pediátrica, esterilizable en autoclave
(Z para conector de teléfono recto de 6.35 mm)
FMT400/AOR/XX (3m), FMT400/AOR/XX-15 (1.5m), sonda esofágica/rectal, adulto
FMT400/POR/XX (3m), FMT400/POR/XX-15 (1.5m), sonda esofágica/rectal, pediátrica
FMT400/IOR/XX (3m), FMT400/IOR/XX-15 (1.5m), sonda esofágica/rectal, lactante
FMT400/AOR/XX-A (3m), FMT400/AOR/XX-A15 (1.5m), sonda esofágica/rectal, adulto, esterilizable en autoclave
FMT400/POR/XX-A (3m), FMT400/POR/XX-A15 (1.5m), sonda esofágica/rectal, pediátrica, esterilizable en autoclave
(XX variable es para monitor compatible: EDN (Edan), HP (HP y Philips), MND (Mindray), S (Drager/Siemens), BLT (Biolight), GE (GE), SW (S&W), THT (TecoTherm), MMM (Belmont), GH (GE))

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FMT400/AOR, FMT400/AOR/Z (3m), **FMT400/AOR-15, FMT400/AOR/Z-15** (1.5m), sonda esofágica/retal, adulto
FMT400/POR, FMT400/POR/Z (3m), **FMT400/POR-15, FMT400/POR/Z-15** (1.5m), sonda esofágica/retal, pediátrica
FMT400/IOR, FMT400/IOR/Z (3m), **FMT400/IOR-15, FMT400/IOR/Z-15** (1.5m), sonda esofágica/retal, infantil
FMT400/AOR-A, FMT400/AOR-Z-A (3 m), **FMT400/AOR-A15, FMT400/AOR-Z-A15** (1.5m), sonda esofágica/retal, adulto, autoclavável
FMT400/POR-A, FMT400/POR-Z-A (3 m), **FMT400/POR-A15, FMT400/POR-Z-A15** (1.5m), sonda esofágica/retal, pediátrica, autoclavável
(Z para plugue de telefone reto de 6.35 mm)

FMT400/AOR/XX (3m), **FMT400/AOR/XX-15** (1.5m), sonda esofágica/retal, adulto
FMT400/POR/XX (3m), **FMT400/POR/XX-15** (1.5m), sonda esofágica/retal, pediátrica
FMT400/IOR/XX (3m), **FMT400/IOR/XX-15** (1.5m), sonda esofágica/retal, infantil

FMT400/AOR/XX-A (3m), **FMT400/AOR/XX-A15** (1.5m), sonda esofágica/retal, adulto, autoclavável
FMT400/POR/XX-A (3m), **FMT400/POR/XX-A15** (1.5m), sonda esofágica/retal, pediátrica, autoclavável

(A variável **XX** é para monitor compatível: EDN (Edan), HP (HP & Philips), MND (Mindray), S (Drager/Siemens), BLT (Biolight), GE (GE), SW (S&W), THT (Tecotherm), MMM (Belmont), GH (GE))

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