**Indications:**

The reusable tubular temperature probe is intended for use in instantaneous spot-check measurement of the rectal temperature of adult or pediatric patients as an indicator of the patient's core body temperature and laboratory temperature measurements by liquid immersion. Reusable tubular temperature probes are intended to be used with YSI 400 series compatible electronic thermometers by trained and qualified medical personnel in hospitals and hospital-type facilities.

Contraindications:

The reusable tubular probe is contraindicated for esophageal or nasal use, long-term temperature monitoring and use simultaneously with electrosurgical devices. Rectal placement is contraindicated in neonates and small infants due to possibility of colon perforation.

Adverse Reactions:

No known adverse reactions to the use of the product.

Clinical Benefits:

A wired temperature sensor enables continuous measurement of patient core (central) temperature.

Product Description:

The use of the temperature probe with incompatible electronic thermometers can result in degraded performance and accuracy. The operator is responsible for checking the compatibility of the monitor and the temperature probe before its use. Probes are accurate and interchangeable with any other probe in the same series at specified measurement temperatures. Within the measurement range of 25 - 45 °C, the accuracy for the probes is $\pm 0.1^{\circ}\text{C}$ ($\pm 0.1^{\circ}\text{C}$ accuracy in the measurement range of 0°C to 70°C for the laboratory). The minimum measurement time for accurate temperature readings is 150 seconds. The temperature reaches a stabilized temperature in 5 minutes. The system accuracy is the sum of the probe and instrument accuracies.

FMT400/ATB (3m), Ø3.8 mm, Tubular temperature probe, Rectal, Adult
FMT400/PTB (3m), Ø3.2 mm, Tubular temperature probe, Rectal, Pediatric

Directions for Use:

- 1- Before use, check that the product is intact and clean. Inspect the temperature probe cable for wear, breakage or fraying. Replace if necessary.
- 2- Make sure that the probe is compatible with the electronic thermometer to be connected and that the plug of the probe matches the socket on the device. To connect the probe to the thermometer, align the plug of the probe with the socket of the thermometer and push firmly to ensure full contact. Forced mating of connectors without proper alignment can result in damage to connectors and loss of electrical continuity.
- 3- Apply the probe to the rectum according to accepted clinical procedures. Verify position of the probe by medical procedures.
- 4- Wait to see the complete measurement reading after the response time is reached. Follow the electronic thermometer device manufacturer's instructions of use and your hospital's protocols on temperature measurement of the patients.
- 5- When temperature measurement is completed, first carefully remove the tubular probe from the patient according to accepted medical techniques then disconnect the probe from the thermometer. Hold the plastic plug firmly and pull. Do not pull on the cable.

Warnings:

- 1- All temperature probes are designed for use with specific patient monitors or electronic thermometers. Verify the compatibility of the monitor, extension cable, and probe and all functions including alarming perform properly before use, otherwise degraded performance, accuracy, and/or patient injury can result.
- 2- Application of temperature probes should be done by trained medical personnel following established procedures.
- 3- Tubular temperature probes are contraindicated to long-term continuous temperature monitoring. Use only for instantaneous spot-check measurements.
- 4- Never use tubular temperature probes simultaneously with electrosurgical devices. Patient burns may result. Remove the probe from patient contact before activating the surgical unit or other RF source.
- 5- It is not recommended to use the temperature probes within the range of Magnetic Resonance Tomography (MRI), Computed Tomography (CT), etc. applications. Conducted current may cause burns.
- 6- The temperature probe is in direct contact with the patient's skin and therefore allergic reactions may occur in susceptible individuals. In case of allergic or skin reactions discontinue use and consult a specialist.
- 7- As with all medical equipment, carefully route temperature probe cables to reduce the possibility of patient entanglement or strangulation.
- 8- Assure that probe cable is not positioned where they could become entangled around IV tubes and restrict IV flow.
- 9- The medical thermometer cannot measure the patient's temperature accurately if the tubular probe is not correctly positioned or the probe tip is not in direct contact with the patient's skin. Check the probe tip is correctly positioned on the patient and that the probe is always in direct contact with the patient's skin during measurement.
- 10- Use only water or medical-grade water-based lubricants to lighten the insertion of the tubular probe. Do not use petroleum-based lubricants
- 11- Never use a tubular probe on the esophageal or nasal way.
- 12- Never use a tubular probe for neonate and infant patients.
- 13- When using the probe always verify the position of the tubular probe with medical procedures.
- 14- Always use caution when applying, inserting, or removing a temperature probe from a patient.
- 15- Inspect and clean or disinfect the probe before attaching it to a new patient.
- 16- Do not use wet probes. This may cause burns during the application of high-frequency devices.
- 17- Mishandling of the probes could result in damage to internal wires and loss of electrical isolation or improper temperature readings. Don't tie the probes in knots or bunch tightly. Don't wrap probes around equipment to avoid damaging internal wires. Don't leave the probes dangling from the instrument. When not in use, probes and leads should be loosely coiled and stored at room temperature.
- 18- In laboratory use junction between the probe and probe cable should not be immersed in the liquid.
- 19- Care must be taken to ensure that the probe cable is not subjected to mechanical stress, such as pulling, bending, and squeezing.
- 20- Never boil temperature probes or clean with an ultrasonic cleaner. Do not immerse the probe in a liquid.
- 21- Never attempt to perform repairs. Do not modify the temperature probes in any way. Modifications may affect patient safety, and performance, and void the warranty.
- 22- Do not use damaged probes. Dispose of damaged probes according to local laws and regulations for medical waste.
- 23- For further information and warnings read the instructions of use accompanying the temperature monitoring device.

Manual Cleaning and Disinfection:**Manual Cleaning:**

The tubular temperature probes must be cleaned and disinfected immediately after each use. Remove heavy debris from the temperature probe by wiping with an appropriate, moist lint-free cloth, sponge, or gauze pad. Prepare the detergent solution with reverse osmosis or distilled water according to the manufacturer's instructions. Saturate a synthetic felt or sponge pad with the cleaning solution. Wipe all surfaces of the temperature probe cable until all visible contamination is removed. Wash the pad in the cleaning solution and repeat wiping the cable. Use a plastic cleaning brush for adhered particles or dried secretions. Wipe or brush all surfaces at least six strokes back and forth. At least allow a contact time of 2 minutes. If visible contamination is still present on the probe at the end of the cleaning repeat the cleaning process. Moist a clean cloth, sponge, or gauze pad with tap water or reverse osmosis and thoroughly wipe the temperature probe to remove any debris and all traces of cleaning solution. Dry manually by using a lint free cloth or absorbent paper and, in particular, for drying cavities, sterile compressed air.

Manual Disinfection:

Prepare the disinfectant solution according to the manufacturer's instructions. Only disinfect clean products that have dried slightly. Saturate a synthetic felt or sponge pad with the disinfectant solution. Wipe all surfaces of the tubular temperature probe. Wash the pad in the disinfectant solution and repeat wiping the probe cable. Wipe or brush all surfaces at least six strokes back and forth. At least allow a contact time of 6 minutes. Moist a clean cloth, sponge, or gauze pad with reverse osmosis or distilled water and thoroughly wipe the temperature probe to remove any residual debris and all traces of disinfectant solution. Dry manually by using a lint free cloth or absorbent paper and, in particular, for drying cavities, sterile compressed air.

Temperature probe, especially metallic parts of connector should never be immersed in cleaning, disinfecting, or rinsing solutions, or placed in an ultrasonic cleaning device. Soaking will reduce the life of temperature probe by oxidation of metal parts and hardening of the cable jacket. Avoid using excessive force when cleaning or disinfecting by hand that may damage the internal wires of the probe and reduces the product life. Please do not clean with hard or frictional materials, abrasive cleaners, or solvents. Cleaning agents containing phenol or chloride are not suitable. The temperature probe should be completely dry before use, do not use wet temperature probes. These instructions were validated using the below detergents and disinfectants.

Cleaning agent (manual cleaning): Enzol® / Cidezime®, enzymatic detergent, Johnson&Johnson (add 8ml enzymatic detergent per 1 liter distilled water, for cable with dried-on organic matter use 16 ml enzymatic detergent per 1 liter distilled water/ warm water).

Disinfectant (manual disinfection): Cidex OPA®, 0.55% ortho-phthalaldehyde solution, Johnson&Johnson (at a minimum temperature of 20°C/68°F). Reverse osmosis/ distilled water.

Caution: During the cleaning and disinfection of temperature probes disposable gloves, protective eyewear, and filtration mask should be used to reduce the risk of transmission of infectious agents by splashing or the creation of aerosols.

Warning: FMT temperature probes are sold NON-STERILE. Clean and disinfect the probes before first and every use.

Warning: Before cleaning or disinfecting the probes, disconnect it from the medical thermometer.

Warning: Temperature probes should be cleaned of overburden prior to disinfection to improve the effectiveness.

Warning: Do not sterilize temperature probes by autoclave, radiation or steam.

Warning: Avoid contact with strong, aromatic, chlorinated, ketone, ether, or ester solvents. Do not use bleach on electrical contacts or connectors.

Environmental Requirements:

Operating Temperature : 5°C - 45°C (41°F - 113°F)

Storage Temperature : 0°C - 50°C (32°F - 122°F)

Relative Humidity : 20 - 80% (not condensing)

During the storage the products should be protected from sunlight. It is recommended to store the products in original packages until the first use. Used probes should be loosely coiled and must be stored individually in a protective container or sterile packaging.

Compatibility:

In order to ensure compatibility and claimed accuracy, the temperature probe should only be used with the equipment for which they have been designed and is specified in the label on the product's primary packaging. All tubular temperature probes are compatible with patient monitors that use 400-Series technology. For product models and compatible devices, please refer to the catalog at metkomedical.com/PDF/04MedicalTemperatureProbes.pdf

Safety:

Degree of protection from electric shocks: type BF

Degree of protection against the ingress of water: Keep dry

FMT reusable tubular temperature probes are sold NON-STERILE.

FMT reusable tubular temperature probes are not made with natural rubber latex.

Refer to the separate leaflet "Description of Symbols" located within the product package for the explanation of symbols.

Reporting of Serious Incidents:

Any serious incident related to the use of this product should be reported to both the manufacturer and the health authority/competent authority where the product is installed. Either contact your local representative or report to: metko@metkomedical.com

Caution: Federal Law (U.S.A.) restricts this device to sale by or on the order of a physician.

Warranty:

FMT reusable tubular temperature probes are under twelve (12) months warranty against material and workmanship defects from the date of original purchase. In warranty period, METKO will be responsible for repairing or changing the probe free of charge if the defect is proven. This warranty does not extend to any product that has been subject to misuse, wrong disinfection, neglect or accident, or that has been damaged by causes external to the product, or that has been used in violation of the operating instructions supplied with the product. The product life is 3 years of the production date.

FMT® is a registered trademark of Metko A.Ş.

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