## Indications:

The reusable oxygenator temperature probe is intended for use in continuous, non-contact blood temperature monitoring by the YSI400 electronic thermometer on an extracorporeal circulation device. The reusable oxygenator temperature probes are intended to be used with compatible electronic thermometers and oxygenators by trained and qualified medical personnel in hospitals and hospital-type facilities.

## Contraindications:

No known contraindications to the use of the product.

## Adverse Reactions:

No known adverse reactions to the use of the product.

## Clinical Benefits:

A wired temperature sensor enables continuous measurement of blood temperature in an oxygenator.

## Product Description:


FMT reusable probes are intended for use with YSI400 compatible electronic thermometers. The use of the temperature probe with an incompatible electronic thermometer can result in degraded performance and accuracy. The operator is responsible for checking the compatibility of the electronic thermometer and the temperature probe before its use. Probes are accurate and interchangeable with any other probe in the same series at specified measurement temperatures. Within the measurement range of 25 - 45 °C, the accuracy for the probes is  $\pm 0.2^{\circ}\text{C}$ . The minimum measurement time for accurate temperature readings is 150 seconds. The temperature reaches a stabilized temperature in 5 minutes.

### FMT400/OXY

(Dideco, Euroset), (3m)  
Tube diameter:  $\varnothing$  2.4 mm  
Tube length: 35 mm

### FMT400/OXY-2

(Maquet), (3m)  
Tube diameter:  $\varnothing$  3.8 mm  
Tube length: 30 mm

## Directions for Use:

- 1- Before use, check that the product is intact and clean. Inspect the temperature probe and cable for rust, wear, breakage or fraying.
- 2- Connect the temperature probe plug to the temperature socket of the electronic thermometer.
- 3- Align the metal BNC type connector to the arterial (or venous) probe port of the oxygenator, insert the stainless steel tube into the hole on the port then twist to lock the BNC connector into the port.
- 4- Follow the oxygenator and heart-lung machine manufacturer's instructions for use and your hospital's protocols on temperature monitoring.

## Warnings:

- 1- All temperature probes are designed for use with specific patient monitors or electronic thermometers. Verify the compatibility of the electronic thermometer, extension cable, and probe and all functions including alarming perform properly before use, otherwise degraded performance, accuracy, and/or patient injury can result.
- 2- Application of temperature probes should be done by trained medical personnel following established procedures.
- 3- All wire-lead patient-connected transducer assemblies are subject to reading error, local heating, and possible damage from high-intensity sources of RF energy. Electrosurgical equipment represents one such source since capacitively-coupled currents may seek alternative paths to ground through probe cables and isolated instruments.
- 4- It is not recommended to use the temperature probes within the range of Magnetic Resonance Tomography (MRI), Computed Tomography (CT), etc. applications.
- 5- Always use caution when applying, inserting, or removing a temperature probe.
- 6- As with all medical equipment, carefully route temperature probe cables to reduce the possibility of patient entanglement or strangulation.
- 7- Temperatures are most accurate when the probe tip firmly contacts the probe port. Make sure the probe tube is fully inserted into the port and locked.
- 8- Mishandling of the probes could result in damage to internal wires and loss of electrical isolation or improper temperature readings. Don't tie the probes in knots or bunch tightly. Don't wrap probes around equipment to avoid damaging internal wires. Don't leave the probes dangling from the instrument. When not in use, probes and leads should be loosely coiled and stored at room temperature.
- 9- Care must be taken to ensure that the probe cable is not subjected to mechanical stress, such as pulling, bending, and squeezing.
- 10- Never boil temperature probes or clean with an ultrasonic cleaner. Do not immerse the probe in a liquid.
- 11- Never attempt to perform repairs. Do not modify the temperature probes in any way. Modifications may affect patient safety, and performance, and void the warranty.
- 12- Do not use damaged probes. Dispose of damaged probes according to local laws and regulations for medical waste.
- 13- For further information and warnings read the instructions of use accompanying the temperature monitoring device and the oxygenator.

## Manual Cleaning and Disinfection:

### Manual Cleaning:

The oxygenator temperature probes must be cleaned and disinfected immediately after each use. Remove heavy debris from the temperature probe by wiping with an appropriate, moist lint-free cloth, sponge, or gauze pad. Prepare the detergent solution with reverse osmosis or distilled water according to the manufacturer's instructions. Saturate a synthetic felt or sponge pad with the cleaning solution. Wipe all surfaces of the temperature probe cable until all visible contamination is removed. Wash the pad in the cleaning solution and repeat wiping the cable. Use a plastic cleaning brush for adhered particles or dried secretions. Wipe or brush all surfaces at least six strokes back and forth. At least allow a contact time of 2 minutes. If visible contamination is still present on the probe at the end of the cleaning repeat the cleaning process. Moist a clean cloth, sponge, or gauze pad with tap water or reverse osmosis and thoroughly wipe the temperature probe to remove any debris and all traces of the cleaning solution. Dry manually by using a lint-free cloth or absorbent paper and, in particular, for drying cavities, sterile compressed air.

### Manual Disinfection:

Prepare the disinfectant solution according to the manufacturer's instructions. Only disinfect clean products that have dried slightly. Saturate a synthetic felt or sponge pad with the disinfectant solution. Wipe all surfaces of the oxygenator temperature probe. Wash the pad in the disinfectant solution and repeat wiping the probe cable. Wipe or brush all surfaces at least six strokes back and forth. At least allow a contact time of 6 minutes. Moist a clean cloth, sponge, or gauze pad with reverse osmosis or distilled water and thoroughly wipe the temperature probe to remove any residual debris and all traces of disinfectant solution. Dry manually by using a lint-free cloth or absorbent paper and, in particular, for drying cavities, sterile compressed air.

Temperature probe, especially metallic parts of connectors should never be immersed in cleaning, disinfecting, or rinsing solutions, or placed in an ultrasonic cleaning device. Soaking will reduce the life of the temperature probe by oxidation of metal parts and hardening of the cable jacket. Avoid using excessive force when cleaning or disinfecting by hand that may damage the internal wires of the probe and reduces the product life. Please do not clean

with hard or frictional materials, abrasive cleaners, or solvents. Cleaning agents containing phenol or chloride are not suitable. The temperature probe should be completely dry before use, do not use wet temperature probes. These instructions were validated using the below detergents and disinfectants.

Cleaning agent (manual cleaning): Enzol® / Cidezyme®, enzymatic detergent, Johnson&Johnson (add 8ml enzymatic detergent per 1 liter distilled water, for cable with dried-on organic matter use 16 ml enzymatic detergent per 1 liter distilled water/ warm water).

Disinfectant (manual disinfection): Cidex OPA®, 0.55% ortho-phthalaldehyde solution, Johnson&Johnson (at a minimum temperature of 20°C/68°F).

Reverse osmosis/ distilled water.

**Caution:** During the cleaning and disinfection of temperature probes disposable gloves, protective eyewear, and filtration mask should be used to reduce the risk of transmission of infectious agents by splashing or the creation of aerosols.

**Warning:** FMT temperature probes are sold NON-STERILE. Clean and disinfect the probes before first and every use.

**Warning:** Before cleaning or disinfecting the probes, disconnect it from the medical thermometer.

**Warning:** Temperature probes should be cleaned of overburden prior to disinfection to improve the effectiveness.

**Warning:** Do not sterilize temperature probes by autoclave, radiation, or steam.

**Warning:** No automated method is recommended. Do not use dishwasher detergents.

**Warning:** Avoid contact with strong, aromatic, chlorinated, ketone, ether, or ester solvents. Do not use bleach on electrical contacts or connectors.

#### Environmental Requirements:

Operating Temperature : 5°C - 45°C (41°F - 113°F)

Storage Temperature : 0°C - 50°C (32°F - 122°F)

Relative Humidity : 20 - 80% (not condensing)

During the storage the product should be protected from sunlight. It is recommended to store the product in its original package until the first use. Used probes should be loosely coiled in 10 to 15 cm diameter (4 to 6 inches) and must be stored individually in a protective container or packaging.

#### Compatibility:

In order to ensure compatibility and claimed accuracy, the temperature probe should only be used with the equipment for which they have been designed and is specified in the label on the product's primary packaging. All FMT oxygenator temperature probes are compatible with electronic thermometers that use 400-Series technology. For product models and compatible devices, please refer to the catalog at [metkomedical.com/PDF/04MedicalTemperatureProbes.pdf](http://metkomedical.com/PDF/04MedicalTemperatureProbes.pdf)

#### Safety:

Degree of protection from electric shocks: type B

Degree of protection against the ingress of water: Keep dry

FMT reusable oxygenator temperature probes are sold NON-STERILE.

FMT reusable oxygenator temperature probes are not made with natural rubber latex.

FMT reusable oxygenator temperature probes are not made with PVC.

Refer to the separate leaflet "Description of Symbols" located within the product package for the explanation of symbols.

#### Reporting of serious incidents:

Any serious incident related to the use of this product should be reported to both the manufacturer and the health authority/competent authority where the product is installed. Either contact your local representative or report to: [metko@metkomedical.com](mailto:metko@metkomedical.com)

**Caution:** Federal Law (U.S.A.) restricts this device to sale by or on the order of a physician.

#### Warranty:

FMT reusable oxygenator temperature probes are under twelve (12) months warranty against material and workmanship defects from the date of original purchase. In warranty period, METKO will be responsible for repairing or changing the probe free of charge if the defect is proven. This warranty does not extend to any product that has been subject to misuse, wrong disinfection, neglect or accident, or that has been damaged by causes external to the product, or that has been used in violation of the operating instructions supplied with the product. The product life is 3 years of the production date.

FMT® is a registered trademark of Metko A.Ş.

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