Indications: FMT Reusable Non-Invasive Blood Pressure (NIBP) cuffs are intended for use with automated oscillometric NIBP devices to measure the systolic and diastolic arterial pressure of the patient. NIBP cuffs are intended to be used by qualified medical personnel in hospitals and hospital-type facilities. NIBP cuffs can be used in infant, pediatric and adult patients, depending on the selection of the appropriate cuff size.

Contraindications: Measurement of blood pressure at the brachial artery is a generally benign procedure. However, there are some circumstances in which obtaining readings from a particular arm may not be appropriate; such circumstances include the presence of an arterial-venous shunt, recent axillary node dissection, or any deformity or surgical history that interferes with proper access or blood flow to the upper arm. If these relative contraindications are present, blood pressure should be assessed in the opposite arm. There may also be pre-existing conditions that can interfere with the accuracy or interpretation of readings, such as aortic coarctation, arterial-venous malformation, occlusive arterial disease, or the presence of an antecubital bruit. If neither arm can be used, then measurement of blood pressure in a leg may be indicated.

Adverse Reactions: No allergic reactions are observed at places where NIBP cuffs are in contact with a patient's intact skin. However, such reactions cannot be completely excluded in people particularly susceptible to allergies; it should be noted that no such events have been reported to the manufacturer. Incorrect cleaning or disinfection, or a failure to perform them, as well as failure to observe hygiene procedures, may cause the transfer of pathogens and infections.

Clinical Benefits: Enables NIBP monitoring with patient monitoring equipment.

Product Description: FMT reusable NIBP cuffs are manufactured from Nylon and are available in bladdered (with internal inflation bag) and bladderless types. Both types have six sizes, depending on the patient's limb size. The choice of single-tube or double-tube cuff depends on the NIBP monitoring system. Connectors attached to the end of the cuff tube (or tubes) allow the cuff to be connected to a compatible NIBP monitor via the NIBP hose.

FMT-SX/BPYY Cuff with bladder, Single Tube

FMT-DX/BPYYY Cuff with bladder, Double Tube


FMT-SBX/BPYY Bladderless Cuff, Single Tube

FMT-DBX/BPYYY Bladderless Cuff, Double Tube

(X is cuff size; 6: Thigh, 5: Large Adult, 4L: Adult Extra Long, 4: Adult, 3: Small Adult, 2: Pediatric, and 1: Infant size, YY is connector number)

Directions for Use:

- 1- Before use, check that the NIBP cuff, cuff tubing, and hose are intact and clean. Replace the cuff when ageing, tearing, or weak closure is apparent. Do not inflate the cuff when unwrapped. Do not use a cuff if its structural integrity is suspect.
- 2- Select the appropriate blood pressure measurement site. Inspect the patient's limb prior to application. Do not apply the cuff to areas where the skin is not intact, or tissue is injured. Because normative values are generally based on this site and as a matter of convenience, the upper arm is the preferred site for the placement of the cuff. The forearm, thigh, and ankle can also be used. If the largest cuff does not fit the upper arm adequately, use an alternate site. When using the upper arm, choose the correct cuff accordingly, taking into account the patient's cardiovascular status and the effect of the alternative site on blood pressure values, proper cuff size, and comfort.
- 3- If the patient is standing, sitting, or inclined, ensure that the cuffed limb is supported to keep the cuff level with the patient's heart. If the cuff is not at heart level, the difference in systolic and diastolic values due to the hydrostatic effect must be considered.
- 4- Select the appropriate cuff size. Measure the patient's limb circumference and select an appropriately sized cuff according to the applicable RANGE printed in centimeters and marked on the cuff. When cuff sizes overlap for a specified circumference, choose the larger size cuff. Accuracy depends on the use of the proper size cuff.
- 5- Remove clothing from the limb before positioning the cuff. Place the cuff on a bare limb. Apply the cuff to the limb such that the "THIS SIDE AGAINST SKIN" printed side contacts the patient. Wrap the cuff around the limb, ensuring the INDEX LINE (INDEX LINE) falls between the RANGE (RANGE) marks on the cuff. Make sure to align the arrow marked ARTERY (ARTERY) over the patient's brachial (or femoral for ankle) artery and the bottom edge of the cuff 1-2 cm above the antecubital fossa.
- 6- Wrap the cuff snugly around the limb but still allow space for a finger between the patient and the cuff. The cuff should not be so tight that it prevents venous return between determinations. Ensure that hook-and-loop closures are properly engaged so that pressure is evenly distributed throughout the cuff.
- 7- Connect the cuff to the NIBP hose, ensuring the proper engagement of the connectors. Ensure that the connection is secure, the tubing is not kinked, and the hose is routed to avoid kinking or compression.
- 8- Start NIBP monitoring. Confirm that the artery marker continues pointing to the brachial artery after each inflation. If it slips, reposition the cuff to the correct position.

Warnings:

- 1- NIBP cuffs are designed for use with compatible automated oscillometric NIBP devices to measure radial arterial pressure on the patient's upper arm. Verify the compatibility of the NIBP cuff, NIBP hose, and NIBP monitoring device all functions perform properly before use, otherwise degraded performance and/or patient injury can result.
- 2- Limb circumference ranges are printed on each cuff. Be sure the cuff is correctly matched to the patient's size. Cuffs that are too small may produce false high-blood-pressure measurements, and cuffs that are too large may produce false low-pressure measurements.
- 3- Before each use, visually inspect the NIBP cuff and connected NIBP hose. The NIBP cuffs should be free of visible defects and damage. Replace the cuff if the hook and loop fail to hold during inflation, artery, index, or range markings are not visible, or the monitor displays a message that a leak has been detected. Never use damaged, clogged, or visibly contaminated NIBP cuffs. Dispose of damaged NIBP cuffs according to hospital disposal practices, local environmental and waste disposal laws, and regulations.
- 4- Do not connect the cuff to intravascular fluid systems. Doing so could allow air to be pumped into a blood vessel, which could lead to serious patient injury.
- 5- Do not apply cuffs to areas where the skin is not intact, tissue is injured, or dermal disruption is at greater risk. Ensure that the rough side of the closure does not contact the skin; contact may cause irritation.
- 6- Do not obtain determinations more frequently than clinically indicated, weighing the benefits of frequent measurement against the risk. Rapidly repeating measurements can impair circulation in the monitored limb.
- 7- Frequently check the NIBP cuff/adaptor/NIBP hose, cuff site, and cuff limb for signs of skin irritation or impeded blood flow, especially when monitoring at frequent intervals and/or extended periods of time. Rotate the site if appropriate.
- 8- Do not allow an NIBP cuff inflated above 10 mm Hg to remain on the patient for more than 10 minutes. This may cause patient distress, disturb blood circulation, and contribute to the injury of peripheral nerves. Remove the cuff from the patient when monitoring is not in process.
- 9- Devices that exert pressure on tissue have been associated with purpura, skin avulsion, compartmental syndrome, ischemia, and/or neuropathy.
- 10- Do not apply the cuff to a limb used for intravenous infusion, arterial monitoring, or any other intravascular access, therapy, or an arterio-venous (A-V) shunt where AV fistulas are present or in areas where circulation is compromised. The cuff inflation can temporarily block blood flow, potentially harming the patient. Assess the limb for the risk of lymphedema (due to mastectomy, etc.).
- 11- Use care when placing the cuff on limbs used to monitor other patient parameters. Do not apply the cuff to a limb being used for SpO2 monitoring. Cuff inflation may cause inaccurate SpO2 measurements.
- 12- Do not apply external pressure against the cuff and avoid contact with it while monitoring. This may cause inaccurate blood pressure values.
- 13- Use the cuff only under direct supervision by a trained healthcare professional when attached to automated monitors without fail-to-read alarms.
- 14- To optimize performance and accuracy, minimize limb movement/cuff motion.
- 15- Avoid placing the cuff on a patient in a manner that can cause the cuff tube to become kinked. Tube kinking may cause inaccurate readings.
- 16- Do not modify the NIBP cuffs in any way. Modifications may affect patient safety, performance, and accuracy and void the warranty.
- 17- As with all medical equipment, carefully route hoses to reduce the possibility of patient entanglement or strangulation.
- 18- Before using the NIBP cuffs, read the instruction manual of the NIBP monitoring device or patient monitor for further information and warnings.

Manual Cleaning and Disinfection:

Manual Cleaning: The reusable NIBP cuffs must be cleaned and disinfected immediately after each use. Remove heavy debris from the NIBP cuff by wiping it with an appropriate, moist, lint-free cloth, sponge, or gauze pad. Prepare the detergent solution with reverse osmosis or distilled water according to the manufacturer's instructions and take the solution into the spray bottle. Spray the detergent solution liberally on the cuff until saturated, and at least allow a contact time of 2 minutes. Wipe smooth surfaces with a soft cloth. Use a soft plastic cleaning brush for adhered particles or dried secretions. If visible contamination is still present on the cuff at the end of the cleaning repeat the cleaning process. Moist a clean cloth, sponge, or gauze pad with tap water or reverse osmosis and thoroughly wipe the cuff to remove any debris and all traces of the cleaning solution. Dry manually by using a lint free cloth or absorbent paper.

Manual Disinfection: Only disinfect cleaned products that have dried slightly. Prepare the disinfectant solution with reverse osmosis or distilled water according to the manufacturer's instructions and take the solution into the spray bottle. Spray the disinfectant solution liberally on the cuff until saturated and at least allow a contact time of 6 minutes. Moist a clean cloth, sponge, or gauze pad with reverse osmosis or distilled water and thoroughly wipe the cuff to remove any residual debris and all traces of disinfectant solution. Dry manually by using a lint free cloth or absorbent paper.

Before cleaning or disinfecting the NIBP cuff, disconnect it from the NIBP hose or the equipment. Cuffs should never be immersed in cleaning, disinfecting, or rinsing solutions, or placed in an ultrasonic cleaning device. Care should be taken to prevent liquid from entering the cuff tubing. The liquid in the tubing may affect blood pressure determination accuracy and damage automatic or manual monitors. Avoid using excessive force when cleaning or disinfecting by hand which may reduce the product life. Please do not clean with hard or frictional materials, abrasive cleaners, or solvents. The NIBP cuffs should be completely dry before use, do not use wet cuffs. These instructions were validated using the below detergents and disinfectants.

Cleaning agent (manual cleaning): Enzo® / Cidezyme®, enzymatic detergent, Johnson&Johnson (add 8ml enzymatic detergent per 1 liter distilled water, for cable with dried-on organic matter use 16 ml enzymatic detergent per 1 liter distilled water/ warm water)

Disinfectant (manual disinfection): 10% solution of household bleach (5.25% sodium hypochlorite) diluted with distilled water.

Reverse osmosis/ distilled water.

Note: When cleaning and disinfecting, care should be taken to prevent liquid from entering the cuff tubing. For cuffs with the bladder, users can remove the inflation bag (bladder) from the cuff, and wash plugs (part number BP00) can be used for bladderless cuffs.

Caution: During the cleaning and disinfection of NIBP cuffs, disposable gloves, protective eyewear, and a filtration mask should be used to reduce the risk of transmission of infectious agents by splashing or creating aerosols.

Warning: FMT reusable NIBP cuffs are sold NON-STERILE. Clean and disinfect the cuffs before use.

Warning: Do not sterilize NIBP cuffs by autoclave, radiation, or steam. Do not press with a hot iron.

Warning: Avoid contact with strong, aromatic, chlorinated, ketone, ether, or ester solvents.

Warning: Prolonged use of these disinfectants at full strength may cause discoloration and degradation.

Environmental Requirements:

Operating Temperature : 5°C - 45°C (41°F - 113°F)

Storage Temperature : 0°C - 50°C (32°F - 122°F)

Relative Humidity : 20 - 80% (not condensing)

The product should be protected from sunlight during storage. It is recommended that it be stored in its original package until its first use.

Compatibility: In order to ensure compatibility, the NIBP cuff should only be used with the equipment for which they have been designed and is specified in the label on the product's primary packaging. For product models and compatible devices, please refer to the catalog at metkomedical.com/PDF/05NIBPCuffsandHoses.pdf

Safety:

Degree of protection from electric shocks: type BF

Degree of protection against the ingress of water: Keep dry

The FMT reusable NIBP cuffs are sold NON-STERILE.

The FMT reusable NIBP cuffs are not made with natural rubber latex and PVC.

The FMT reusable NIBP cuffs have not been evaluated for safety in the MR environment. Scanning a patient who has this device may result in patient injury.

For an explanation of symbols, refer to the separate leaflet "Description of Symbols" located within the product package.

Reporting of serious incidents: Any serious incident related to the use of this product should be reported to both the manufacturer and the health authority/competent authority where the product is installed. Either contact your local representative or report to: metko@metkomedical.com

Caution: Federal Law (U.S.A.) restricts this device to sale by or on the order of a physician.

Warranty: FMT reusable NIBP cuffs are under twelve (12) months warranty against material and workmanship defects from the date of original purchase. In the warranty period, METKO will be responsible for repairing or replacing the cuff free of charge if the defect is proven. This warranty does not cover any product that has been subjected to misuse, improper disinfection, neglect, or accident and that has been damaged by causes external to the product, or that has been used in violation of the operating instructions supplied with the product. The product life is 3 years of the production date.

FMT® is a registered trademark of Metko A.Ş.

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