# Reusable Incubator Skin Temperature Probe and Air Auxiliary Temperature Probe

English

## Indications:

The reusable skin temperature probe is indicated for the routine continuous monitoring of the skin temperature of the infant in the incubator, infant warmer, or phototherapy unit. The air auxiliary temperature probe is indicated for continuous measurement and control of the air temperature in the hood of the incubator. The reusable skin temperature probes and air auxiliary temperature probes are intended to be used with compatible incubators, infant warmers, or phototherapy units by trained and qualified medical personnel in hospitals and hospital-type facilities.

## Contraindications:

The use of a skin temperature probe is contraindicated for usage over traumatized and hypo-perfused areas.

## Adverse Reactions:

Adverse reactions reported during applications of skin temperature probes include skin sensitizations or irritations due to adhesive cover.

## Clinical Benefits:

A wired temperature sensor enables continuous measurement of patient skin temperature in a warmer or air temperature in the warmer.

## Product Description:

FMT reusable skin and air temperature probes are intended for use with compatible electronic thermometers and NICU devices. The use of the temperature probe with an incompatible electronic thermometer or device can result in degraded performance and accuracy. The operator is responsible for checking the compatibility of the device and the temperature probe before its use. Probes are accurate and interchangeable with any other probe in the same series at specified measurement temperatures. Within the measurement range of 25 - 45 °C, the accuracy for the probes is  $\pm 0.1^{\circ}\text{C}$ . The minimum measurement time for accurate temperature readings is 150 seconds. The temperature reaches a stabilized temperature in 5 minutes, depending on the body site. The system accuracy is the sum of the probe and instrument accuracies.

The reusable skin temperature probe is designed for placement on the skin surface. In order to fix the sensor and minimize the effect of the ambient temperatures on temperature measurement, the sensor must be covered with a heat-reflective adhesive cover. The heat-reflective adhesive cover is made of metallic-looking Mylar film, closed-cell foam, and backing paper to protect the adhesive surface. The adhesive on the foam surface holds the skin sensor in situ.

- FMT-INC/S** Skin temperature probe for Air Shields / Hillrom - Drager Incubators (OEM P/N 68 209 70 – MU06951)
- FMT-INC/SO** Skin temperature probe for Ohmeda / Ohio Incubators (OEM P/N 0208 0697 700 – 6600 0628 700)
- FMT-INC/SO2** Skin temperature probe for Ohmeda / Ohio Incubators (OEM P/N 0208 0161 300 – 0208 0754 700)
- FMT-INC/SIS** Skin temperature probe for Air Shields / Hillrom - Drager Isolette Incubators (OEM P/N 83 102 41 – MU12533)
- FMT-INC/SGR** Skin temperature probe for GE Ohmeda Giraffe Incubators (OEM P/N 6600 0875 700)
- FMT-INC/SA, FMT-INC/SA-S** Skin temperature probe for Atom Medical Incubators (OEM P/N 60882) (OEM P/N 60883)
- FMT-INC/ST, FMT-INC/ST-S** Skin temperature probe for Atom Medical Incubators (OEM P/N 60881)
- FMT-INC/SND** Skin temperature probe for Ningbo David Incubators (OEM P/N 02 025 00)
- FMT-INC/SND2** Skin temperature probe for Ningbo David Incubators (OEM P/N P02 025 00)
- FMT-INC/SND3** Skin temperature probe for Ningbo David Incubators (OEM P/N 23 07 064 00)
- FMT-INC/SFP** Skin temperature probe for Fisher & Paykel Incubators (OEM P/N 900IW001)
- FMT-INC/SAD** Skin temperature probe for Air Shields / Hillrom - Drager Incubators (OEM P/N 68 209 72)
- FMT-INC/SRS** Skin temperature probe for Air Shields / Hillrom - Drager Incubators (OEM P/N 81 300 05 – MU112205-3)
- FMT-INC/MDP** Skin temperature probe for Mediprema Satis Incubators (OEM P/N 3617)
- FMT-INC/MDP2** Skin temperature probe for Mediprema Nite Incubators (OEM P/N 5038)
- FMT-INC/A** Air Auxiliary probe (single sensor) for Air Shields / Hillrom - Drager Incubators (OEM P/N 68 209 80 – MU06955)
- FMT-INC/AA** Air Auxiliary probe (dual sensors) for Air Shields / Hillrom - Drager Incubators (OEM P/N 68 209 81– MU06956)

## Directions for Use:


The reusable temperature probe and its plug should be inspected for visual defects. Probes with a loose or damaged plug or contact pins or corrosion of metal contacts, damage in the insulation of cable (i.e. cuts, cracks, brittle or solid, permanent bends, burnt or melted), wires or metal parts that have come out of the insulation of the cable or plug, or not completely dried and ventilated probes should not be used.

### Application of skin temperature probe:

- 1- Before use, check that the product is intact and clean.
- 2- Clean and dry skin completely in the area intended for probe placement.
- 3- Make sure that the probe is compatible with the monitor to be connected and that the plug of the probe matches the socket on the monitor. To connect the probe to the temperature monitor, align the plug of the probe with the socket of the monitor or extension cable and push firmly to ensure full contact. Forced mating of connectors without proper alignment can result in damage to connectors and loss of electrical continuity.
- 4- There is a disc-shaped stainless steel part at the patient end of the skin temperature probe. The disc is covered on one side with epoxy. Always apply the probe with the metal side of the probe disc in contact with the skin.
- 5- Use a heat-reflective adhesive cover to fix the probe disc to the skin. Remove the backing paper of the reflective cover and place the adhesive side over the disk of the skin probe. Fix the skin probe to the selected skin site by gently pressing the edges of the probe cover. If you are using medical adhesive tape, choosing foam-type medical tapes provides better insulation against ambient temperature.
- 6- A piece of medical tape can be placed over the probe cable approximately 4 or 5 centimeters from the probe tip to stabilize the attached probe.
- 7- Wait to see the complete measurement reading after the response time is reached. Follow the monitoring device manufacturer's instructions of use and your hospital's protocols on temperature monitoring of the patients.
- 8- When temperature monitoring is completed, first remove the heat-reflective adhesive cover then carefully remove the skin probe tip from the patient's skin then disconnect the probe from the monitor. Hold the plastic plug firmly and pull. Do not pull on the cable.

### Application of air auxiliary temperature Probe:

Connect the probe's plug to the air auxiliary probe socket. Insert the probe through the hole at the top of the hood and press and insert the spacer plastic into the hole. Make sure the probe tip hangs freely in the incubator.

## Warnings:

- 1- All temperature probes are designed for use with specific patient monitors or electronic thermometers. Verify the compatibility of the monitor, extension cable, and probe and all functions including alarming perform properly before use, otherwise degraded performance, accuracy, and/or patient injury can result.
- 2- The application of temperature probes should be done by trained medical personnel following established procedures.
- 3- All wire-lead patient-connected transducer assemblies are subject to reading error, local heating, and possible damage from high-intensity sources of RF energy. Electrosurgical equipment represents one such source since capacitively-coupled currents may seek alternative paths to ground through probe cables and isolated instruments. Patient burns may result. If possible, remove the probe from patient contact before activating the surgical unit or other RF source. If probes must be used simultaneously with electrosurgical apparatus, the instruments to which the probes are connected should be checked for adequate isolation from electrical grounds at radio frequencies. Hazards can be reduced by selecting a temperature monitoring point located away from the expected RF current path between the active probe and the return pad.
- 4- It is not recommended to use the temperature probes within the range of Magnetic Resonance Imaging (MRI), Computed Tomography (CT) etc. applications. Conducted current may cause burns.
- 5- The temperature probe is in direct contact with the patient's skin and therefore allergic reactions may occur in susceptible individuals. In case of allergic or skin reactions discontinue use and consult a specialist.
- 6- As with all medical equipment, carefully route temperature probe cables to reduce the possibility of patient entanglement or strangulation.
- 7- Assure that probe cable is not positioned where they could become entangled around IV tubes and restrict IV flow.
- 8- The medical thermometer cannot measure the patient's temperature accurately if the probe is not correctly positioned or the probe tip is not firmly contacting the body site. Regularly check the probe remains correctly positioned on the patient and probe is always in direct contact with the patient's skin.
- 9- Make sure the skin temperature probe is free from pressure. Never place skin temperature probes underneath the patient.
- 10- Do not use a skin probe without a reflective probe cover. Do not remove the reflective foil from the reflective probe cover. Do not use any type of lotion or dressing between patient skin and adhesive probe cover. Replace the probe cover when repositioning the skin probe or when adhesive strength degrades.
- 11- Place skin temperature probe on soft muscle area. Avoid any bony area.
- 12- Do not apply the temperature probe to a skin surface that has been damaged by abrasion, laceration, incision, or burning.
- 13- When monitoring the infants, if the infant is on its back, the skin probe should be placed on the abdomen, halfway between the xyphoid and the umbilicus. When the infant is in the prone position, the skin probe should be on the infant's back. If the baby is placed in the lateral position, ensure the skin probe is located on the side of the infant and clear of the infant's arm.

- 14- To ensure the continued comfort of the infant, avoid pressing directly on the skin probe when attaching it to the infant. Always use caution when applying or removing a temperature probe from a patient.
- 15- Do not leave the patient unattended when using the incubator. Check the patient's temperature periodically to ensure the comfort and safety of the patient.
- 16- Do not place the incubator in direct sunlight or in proximity to any other source of radiant heat. External heat sources can cause elevated interior temperatures that may overheat the infant. Phototherapy units located too close to the incubator may affect the hood wall temperature, incubator temperature, and infant skin temperature. Patient injury could occur. Do not place phototherapy units too close to the incubator.
- 17- When using phototherapy lamps, the probe must be directly in the path of the radiant heat of the lamp, do not place the probe in an area shielded from the lamp's light. The phototherapy lamp may raise the infant's skin temperature.
- 18- Do not use an auxiliary or rectal temperature measurement to control the infant's skin temperature. It can result in overheating or under heating of the infant.
- 19- Mishandling of the probes could result in damage to internal wires and loss of electrical isolation or improper temperature readings. Don't tie the probes in knots or bunch tightly. Don't wrap probes around equipment to avoid damaging internal wires. Don't leave the probes dangling from the instrument.
- 20- If the Air Auxiliary Probe contacts the infant's skin erratic temperature control can result. Always ensure that the probe tip does not contact the infant. Doing so could result in personal injury.
- 21- Do not operate the incubator with Air Auxiliary Probe outside the incubator or other than as recommended. Overheating will occur, and personal injury could result.
- 22- Care must be taken to ensure that the probe cable is not subjected to mechanical stress, such as pulling, bending, and squeezing.
- 23- Portable and mobile RF communications equipment can affect medical electrical equipment and should be used no closer to any part of the equipment.
- 24- The probe is not protected against the effect of the discharge of a cardiac defibrillator. Remove the probe before cardiac defibrillation.
- 25- Never attempt to perform repairs. Do not modify the temperature probes in any way. Modifications may affect patient safety, and performance, and void the warranty.
- 26- Clean and disinfect the probe before attaching it to a new patient. Keep the probe tip clean and free from foreign matter, particularly adhesives.
- 27- Do not autoclave or boil the temperature probes. Do not immerse the probe in a liquid. Do not use wet probes.
- 28- Do not use damaged probes. Dispose of damaged probes according to local laws and regulations for medical waste.
- 29- For further information and warnings read the instructions of use accompanying the temperature monitoring device.

#### **Manual Cleaning and Disinfection:**

##### **Manual Cleaning:**

The temperature probes must be cleaned and disinfected immediately after each use. Remove heavy debris from the temperature probe by wiping with an appropriate, moist lint-free cloth, sponge, or gauze pad. Prepare the detergent solution with reverse osmosis or distilled water according to the manufacturer's instructions. Saturate a synthetic felt or sponge pad with the cleaning solution. Wipe all surfaces of the temperature probe cable until all visible contamination is removed. Wash the pad in the cleaning solution and repeat wiping the cable. Use a plastic cleaning brush for adhered particles or dried secretions. Wipe or brush all surfaces at least six strokes back and forth. At least allow a contact time of 2 minutes. If visible contamination is still present on the probe at the end of the cleaning repeat the cleaning process. Moist a clean cloth, sponge, or gauze pad with tap water or reverse osmosis and thoroughly wipe the temperature probe to remove any debris and all traces of the cleaning solution. Dry manually by using a lint-free cloth or absorbent paper and, in particular, for drying cavities, sterile compressed air.

##### **Manual Disinfection:**

Prepare the disinfectant solution according to the manufacturer's instructions. Only disinfect clean products that have dried slightly. Saturate a synthetic felt or sponge pad with the disinfectant solution. Wipe all surfaces of the temperature probe. Wash the pad in the disinfectant solution and repeat wiping the probe cable. Wipe or brush all surfaces at least six strokes back and forth. At least allow a contact time of 6 minutes. Moist a clean cloth, sponge, or gauze pad with reverse osmosis or distilled water and thoroughly wipe the temperature probe to remove any residual debris and all traces of disinfectant solution. Dry manually by using a lint-free cloth or absorbent paper and, in particular, for drying cavities, sterile compressed air.

Temperature probe, especially metallic parts of connector should never be immersed in cleaning, disinfecting, or rinsing solutions, or placed in an ultrasonic cleaning device. Soaking will reduce the life of the temperature probe by oxidation of metal parts and hardening of the cable jacket. Avoid using excessive force when cleaning or disinfecting by hand that may damage the internal wires of the probe and reduces the product life. Please do not clean with hard or frictional materials, abrasive cleaners, or solvents. Cleaning agents containing phenol or chloride are not suitable. The temperature probe should be completely dry before use, do not use wet temperature probes. These instructions were validated using the below detergents and disinfectants.

Cleaning agent (manual cleaning): Enzol® / Cidezyme®, enzymatic detergent, Johnson&Johnson (add 8ml enzymatic detergent per 1 liter distilled water, for cable with dried-on organic matter use 16 ml enzymatic detergent per 1 liter distilled water/ warm water)

Disinfectant (manual disinfection): Cidex OPA®, 0.55% ortho-phthalaldehyde solution, Johnson&Johnson (at a minimum temperature of 20°C/68°F).  
Reverse osmosis/ distilled water.

**Caution:** During the cleaning and disinfection of temperature probes disposable gloves, protective eyewear, and filtration mask should be used to reduce the risk of transmission of infectious agents by splashing or the creation of aerosols.

**Warning:** FMT temperature probes are sold NON-STERILE. Clean and disinfect the probes before first and every use.

**Warning:** Before cleaning or disinfecting the probe, disconnect it from the medical thermometer.

**Warning:** Temperature probes should be cleaned of overburden prior to disinfection to improve the effectiveness.

**Warning:** Do not sterilize temperature probes by autoclave, radiation, or steam.

**Warning:** Avoid contact with strong, aromatic, chlorinated, ketone, ether, or ester solvents. Do not use bleach on electrical contacts or connectors.

#### **Environmental Requirements:**

Operating Temperature : 5°C - 45°C (41°F - 113°F)

Storage Temperature : 0°C - 50°C (32°F - 122°F)

Relative Humidity : 20 - 80% (not condensing)

During the storage, the product should be protected from sunlight. It is recommended to store the product in its original package until the first use. Used probes should be loosely coiled in 10 to 15 cm diameter (4 to 6 inches) and must be stored individually in a protective container or packaging.

#### **Compatibility:**

In order to ensure compatibility and claimed accuracy, the temperature probe should only be used with the equipment for which they have been designed and is specified in the label on the product's primary packaging. For product models and compatible devices, please refer to the catalog at [metkomedical.com/PDF/04MedicalTemperatureProbes.pdf](http://metkomedical.com/PDF/04MedicalTemperatureProbes.pdf)

#### **Safety:**

Degree of protection from electric shocks: type BF

Degree of protection against the ingress of water: IPX2

The FMT reusable skin temperature probes and air auxiliary probes are sold NON-STERILE.

The FMT reusable temperature probes are not made with natural rubber latex and PVC.

The FMT reusable temperature probes have not been evaluated for safety in the MR environment. Scanning a patient who has this device may result in patient injury.

Refer to the separate leaflet "Description of Symbols" located within the product package for the explanation of symbols.

#### **Reporting of serious incidents:**

Any serious incident related to the use of this product should be reported to both the manufacturer and the health authority/competent authority where the product is installed. Either contact your local representative or report to: [metko@metkomedical.com](mailto:metko@metkomedical.com)

**Caution:** Federal Law (U.S.A.) restricts this device to sale by or on the order of a physician.

#### **Warranty:**

FMT reusable temperature probes are under twelve (12) months warranty against material and workmanship defects from the date of original purchase. In warranty period, METKO will be responsible for repairing the probe or change the probe free of charge if the defect is proven. This warranty does not extend to any product that has been subject to misuse, wrong disinfection, neglect or accident, or that has been damaged by causes external to the product, or that has been used in violation of the operating instructions supplied with the product. The product life is 3 years of the production date.

FMT® is a registered trademark of Metko A.Ş.

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