# Reusable Electrosurgical Dispersive Cable

English

**Indications:** FMT reusable electrosurgical dispersive cables are intended for connecting a non-corded disposable grounding plate (patient return electrode) to a compatible electrosurgical generator (high-frequency or electrosurgery unit). Reusable electrosurgical dispersive cables are intended to be used with appropriate electrosurgical equipment by trained and qualified medical personnel in hospitals and hospital-type facilities.

**Contraindications:** While not an absolute contraindication to the use of electrosurgery, implantable electrical devices (IEDs) such as cardiac and gastric pacemakers, implantable cardiac defibrillators (ICDs), cochlear implants, and deep brain, nerve, spinal cord, or bone stimulators create unique risks. HF current during electrosurgery can interfere with the functioning of these devices and damage them. Patients with these devices require a thorough pre-procedure evaluation and may require intraoperative monitoring and post-procedure device assessment.

**Adverse Reactions / Events:** Adverse events have been reported when using electrosurgical devices such as tissue damage and equipment damage at the wrong site caused by unintentional activation, fires involving surgical drapes and other flammable materials, alternating current paths that cause burns where the patient or operator is in contact with exposed metal, explosions caused by an electrosurgical spark in the flammable gas environment, organ perforation, sudden massive bleeding.

**Clinical Benefits:** Enables safe return for electrosurgical current during monopolar electrosurgery.

**Product Description:** FMT reusable electrosurgical dispersive cable is an accessory in an electrosurgical system. The reusable electrosurgical dispersive cable consists of a plastic or silicone-insulated cable, a metallic connector that fits into the electrosurgical generator, and a clamping mechanism that securely connects to a disposable grounding plate. Blue colored reusable dispersive cables are non-autoclavable. Autoclave sterilization is only applied to the black/dark grey colored autoclavable dispersive cables. Reusable dispersive cables have different sizes of disposable plate-side clamps and generator-side connectors for different electrosurgical generator (high-frequency unit) manufacturers. Always consider the information on the label of the product and/or its primary packaging for compatibility.

**RDC-X3Y** (3 m), **RDC-X5Y** (5 m), dispersive cable for standard disposable grounding plates with 23.5mm tab.

(X: A for Ø6.35 mm mono plug, B for ERBE Ø6.35 mm mono plug, C for VALLEYLAB REM plug, C for VALLEYLAB NON-REM plug, D for VALLEYLAB Ø4mm plug, F for BOVIE Ø8mm plug, G for ESCHMANN dual Ø4mm plug, H for ADJUSTABLE dual Ø4mm plug; I for PETAŞ REM plug; Y: empty for non-autoclavable, A for autoclavable)



**RDC-VX3Y** (3 m), **RDC-VX5Y** (5 m) Dispersive cable for Valleylab E7509 and E7509B disposable grounding plates with 32mm tab.

(X: A for Ø6.35 mm mono plug, B for ERBE Ø6.35 mm mono plug, L for VALLEYLAB REM plug, Y: empty for non-autoclavable, A for autoclavable)

Maximum accessory rated voltage (U<sub>max</sub>):  
Dispersive cables 0.5 kVp

**Directions for Use:** Reusable electrosurgical dispersive cables and their connector and clamp should be inspected for visual defects. Loose or damaged connector and clamp mating parts and corrosion of metal contacts, damage in the insulation of the cable (i.e. cuts, cracks, brittle or solid, permanent bends, burnt or melted), wires or metal parts that have come out of the insulation of the cable (or connector and clamp), lose or improperly attached connectors, and not completely dried and ventilated cables should not be used.

- 1- Dispersive cables are used to transmit the returned HF current to the ground. Therefore, a suitable disposable grounding plate should be connected to the clamp side of the cable. Connection to the high-frequency unit/generator is made with the plug of the dispersive cable.
- 2- Check whether the dispersive cable plug is suitable for the correct socket on the generator. Careless push or pushing with force may result in damage to the plug and/or the socket on the generator.
- 3- Ensure the clamp lever is in the open (up) position. Insert the tab with the foil contacts of the disposable plate into the slot on the clamp completely. Ensure that the entire foil contact area of the tab is covered by the clamp and inserted up to the disposable plate backing material.
- 4- To secure the disposable plate depress down the lever arm until it is flush with the body of the clamp so that the tab contact is held firmly by the clamp. Make sure the tab of the plate is centered in the jaws of the clamp.
- 5- Remove the back paper from the grounding plate and apply the disposable plate to the patient according to your hospital protocols and instructions for the grounding plate.
- 6- Verify the function of the generator's audible alarm circuit by turning on the generator before attaching the dispersive cable. An alarm should sound. If no alarm sounds, check the alarm volume adjustment. If not working properly, return the generator for service. When the alarm sounds, insert the dispersive cable plug into the generator's dispersive electrode receiver and the alarm should cease. If it does not, try the procedure again. If the alarm does not stop after repeating this procedure, exchange the generator and/or cable and try the test again.
- 7- Connect the dispersive cable plug to the generator. Ensure that the electrosurgical generator is in the OFF position. When connecting the dispersive cable, it is important to check whether good and continuous contact is obtained. This is obtained by pushing the cable plug into the mechanical borders of the socket on the generator covering all the metal contacts.
- 8- Once cable continuity is established, check the generator's CQM System if applicable. If the CQM System alarms at this time or during surgery, carefully check the plate-to-patient contact, cable, and clamp connections. If the alarm fault is not found, replace the dispersive electrode and/or cable. If this does not satisfy the alarm fault, replace the generator.
- 9- Before disconnecting the dispersive cables always switch the HF unit off. When connecting and disconnecting cables, hold only the plastic part of the plug. Otherwise, it may result in a burn or electric shock in patients or users. Don't pull on the cable. In preparation, transport, sterilization, and storage dispersive cables should not be bent, twisted tightly, or folded. Such movements reduce the life of the product. Do not load heavy or sharp objects on the cable. This may damage the cable insulation. Damage to the insulation may result in burns with the unwanted electric current passes.
- 10- At the conclusion of the surgical procedure, Lift the lever to release the clamp and disconnect the disposable plate from the cable then remove the disposable grounding plate from the patient. To avoid skin trauma when removing the disposable plate, peel it off slowly with one hand while supporting the underlying tissue with the other hand.

## Warnings:

- 1- The user of dispersive cables should have the necessary knowledge, expertise, and education to use the products. Read all instruction manuals for disposable grounding plate and generator thoroughly. An insufficient understanding of the dangers, warnings, cautions, and information in these manuals can result in death, serious injury, or equipment damage.
- 2- Before first and each use, an electrical continuity test must be performed and the insulation of the electrosurgical cable must be inspected for damage. If there is a break in the cable wire or the dispersive cable becomes otherwise electrically discontinuous, arcing may occur in the patient-return or active circuit and may burn the patient or user or create a fire. Visually inspect the entire cable, connector, and clamp. Visual inspection alone may not be sufficient to ensure that the insulation is intact. Inspection may be done visually under magnification or with a high-voltage insulation testing device. Do not use a cable with brittle or defective insulation or exhibiting pressure points.
- 3- Reusable dispersive cables should be inspected, cleaned, sterilized, and dried/aerated prior to use per this instruction manual prior to first and each subsequent use. The user is recommended to implement a program to track the number of uses of the cables. The serial number on the product allows the user to keep track of the dispersive cable and record reprocessing cycles or the service life of each cable.
- 4- Properly reprocess the dispersive cables before first and each subsequent use following the instructions in this manual. Improper and/or incomplete reprocessing can cause infection of the patient and/or medical personnel.
- 5- Do not perform electrosurgery in the presence of flammable anesthetics or oxidizing gases, other flammable gases, or in close proximity to volatile solvents, fluids, or objects, or in oxygen-enriched environments, as explosions and fire could result. Do not place electrosurgical instruments near or in contact with flammable materials (such as gauze or surgical drapes). Instruments that are activated or hot from use may cause a fire.
- 6- During the use of the high-frequency equipment in patients with cardiac pacemakers, implanted heart defibrillators, or pace electrodes there is a danger to ruin the function or the device itself. The pacemaker or electrodes should never be touched with high-frequency electrodes. Before use advice from cardiology departments must always be taken.
- 7- Dispersive cables should only be connected to the high-frequency generator when the unit is switched off. Make sure that the dispersive cable is not connected while preparing the dispersive grounding plate. Failure to do so may result in an injury or electrical shock to the patient or user.
- 8- Connect dispersive cable connectors to the return receptacle of the generator only. Improper connection of cables may result in inadvertent accessory activation or other potentially hazardous conditions.
- 9- Make sure that the patient is electrically isolated against electrical conductive elements and does not come into direct contact with other metal parts or electric conductive instruments which are grounded or which have an appreciable capacitance to ground. Metal devices can provide an alternate pathway for the electrical current and be responsible for patient burns. Avoid skin-to-skin contact of the patient.
- 10- Always provide a safe working environment. Avoid sharp edges and trip hazards. Avoid kinking or sharply bending the dispersive cables, or damaging the insulation. Failure to do so may result in damage to the insulation and result in burns with the unwanted electric current passes.
- 11- The dispersive cables should not touch directly the patient's skin or other cables and should be kept at minimum length and should be plugged in without any friction to any place. Do not lay high-frequency cables in loops.
- 12- To avoid interference with the camera cables, high-frequency cables should not be positioned in parallel with the camera cables.
- 13- Always plug and unplug the dispersive cables by holding on to the connector. Do not touch the metal parts of the connector. Never pull at the cable.
- 14- Dispersive cables whose connectors have not been thoroughly rinsed and dried/aerated may cause electrical shocks or burns to the patient or user.
- 15- Apparent low output or failure of the generator to function correctly at the normal operating settings may indicate faulty application of the dispersive (neutral) electrode or poor contact in its connections. In this case, prior to increasing the output power, check the adherence of the dispersive (neutral) electrode and its connections.
- 16- For patients with electrically conductive implants, a possible hazard exists due to the concentration or re-direction of HF currents. In case of doubt, qualified advice should be obtained.
- 17- The Contact Quality Monitoring System (CQMS, REM, ARM, NESSY, etc.) will not function if the incorrect dispersive cable is used. The type of generator will determine which cable should be used.
- 18- Aspirate fluid from the area before activating the HF instrument. Conductive fluids (e.g., blood or saline) in direct contact with or in close proximity to an active electrode may carry electrical current or heat away from target tissues, which may cause unintended burns to the patient. Flammable solutions pooling under the patient or in body depressions may cause burns to the patient.
- 19- Potentially hazardous conditions may exist when accessories are not appropriate for the type of generator output used. If unusual operating characteristics develop or are observed during procedures, the use of the product should be discontinued.
- 20- The dispersive cables should not be wrapped around metal objects or bundled together. This may induce currents that could lead to shocks, fires, or injury to the patient or user.
- 21- The patient monitoring electrodes, such as ECG electrodes, should be placed as far apart from the surgical site as feasible. Do not use ECG needle electrodes for monitoring. Needle monitoring electrodes provide an alternate pathway for the electrical current placing the patient at risk for burns. All ECG electrodes must be equipped with protective impedances or HF choke coils.

- 22- In operations where the high-frequency current is passing from the very minor parts of the body, it is preferred to use the bipolar technique instead of monopolar use to avoid unwanted coagulation.
- 23- In Monopolar electrosurgery, patient laying position, the use and the proper placement of a dispersive (neutral) electrode is a key element in the safe and effective use, particularly in the prevention of burns.
- 24- Check this instruction manual, the label, or the current product catalog for the maximum rated voltage of the HF cable. If anything is unclear, contact the manufacturer. In any combination with another electrosurgical accessory, the maximum rated voltage of the combination corresponds to the lowest rated voltage of the accessories used. Do not use bipolar or monopolar HF cables with a generator set to output voltage above the specified maximum accessory-rated voltage.
- 25- Reprocessing and mechanical stress damage the dispersive cable depending on the intensity of use. Damage caused by misuse, overuse, or normal wear may cause electrical shocks or burns to the patient and/or user.
- 26- If this product is used on a patient with or suspected of having CJD (Creutzfeldt-Jacob disease), vCJD (variant Creutzfeldt-Jacob disease), BSE (Bovine Spongiform Encephalopathy), and TSE (Transmissible spongiform Encephalopathy) the product cannot be reused and must be destroyed to eliminate the risk of cross-contamination. Such pathogens cannot be killed by standard disinfection and sterilization methods.
- 27- Never attempt to perform repairs. Do not modify the dispersive cables in any way. Modifications may affect patient safety, and performance, and void the warranty.
- 28- Do not use damaged electrosurgical cables. Dispose of damaged dispersive cables according to national and local laws and regulations for medical waste.
- 29- Refer to the manuals of the high-frequency generator and dispersive (neutral) electrode for additional cautions and warnings.

#### Manual Cleaning and Disinfection:

**Manual Cleaning:** The dispersive cables must be cleaned and disinfected immediately after each use. Remove heavy debris from the dispersive cable by wiping with an appropriate, moist lint-free cloth, sponge, or gauze pad. Prepare the detergent solution with deionized or distilled water according to the manufacturer's instructions. Saturate a synthetic felt or sponge pad with the cleaning solution. Wipe all surfaces of the dispersive cable until all visible contamination is removed. Wash the pad in the cleaning solution and repeat wiping the dispersive cable. Use a plastic cleaning brush for adhered particles or dried secretions. Wipe or brush all surfaces at least six strokes back and forth. At least allow a contact time of 2 minutes. If visible contamination is still present on the cable at the end of the cleaning repeat the cleaning process or safely dispose of the cable. Moist a clean cloth, sponge, or gauze pad with deionized or distilled water and thoroughly wipe the dispersive cable to remove any debris and all traces of the cleaning solution. Dry manually by using a lint-free cloth or absorbent paper and, in particular, for drying cavities, sterile compressed air.

**Manual Disinfection:** Prepare the disinfectant solution according to the manufacturer's instructions. Only disinfect clean products that have dried slightly. Saturate a synthetic felt or sponge pad with the disinfectant solution. Wipe all surfaces of the dispersive cable. Wash the pad in the disinfectant solution and repeat wiping the cable. Wipe or brush all surfaces at least six strokes back and forth. At least allow a contact time of 6 minutes. Moist a clean cloth, sponge, or gauze pad with deionized or distilled water and thoroughly wipe the dispersive cable to remove any residual debris and all traces of disinfectant solution. Dry manually by using a lint-free cloth or absorbent paper and, in particular, for drying cavities, sterile compressed air.

Dispersive cables should never be immersed in cleaning, disinfecting, or rinsing solutions, or placed in an ultrasonic cleaning device. Soaking will reduce the life of cable by oxidation of metal parts and hardening of the cable jacket. Avoid using excessive force when cleaning or disinfecting by hand. It is not recommended to use saline solutions as a final rinse as saline solutions may interfere with subsequent disinfection or sterilization steps. Please do not clean with hard or frictional materials, abrasive cleaners, or solvents. Cleaning agents containing phenol or chloride are not suitable. These instructions were validated using the below detergents and disinfectants.

Cleaning agent (manual cleaning): Enzo® / Cidezime®, enzymatic detergent, Johnson&Johnson (add 8ml enzymatic detergent per 1 liter distilled water, for cable with dried-on organic matter use 16 ml enzymatic detergent per 1 liter distilled water/ warm water)

Disinfectant (manual disinfection): Cidex OPA®, 0.55% ortho-phthalaldehyde solution, Johnson&Johnson (at a minimum temperature of 20°C/68°F).

Reverse osmosis/ distilled water

**Caution:** During the cleaning and disinfection of dispersive cables disposable gloves, protective eyewear, and filtration mask should be used to reduce the risk of transmission of infectious agents by splashing or the creation of aerosols.

**Warning:** FMT reusable dispersive cables are sold NON-STERILE. Clean and disinfect the cables before first and every use.

**Warning:** Before cleaning or disinfecting the cables, disconnect it from the high-frequency generator.

**Warning:** Dispersive cable should be cleaned of overburden prior to disinfection to improve the effectiveness.

**Warning:** Avoid contact with strong, aromatic, chlorinated, ketone, ether, or ester solvents. Do not use bleach on electrical contacts or connectors.

**Warning:** Avoid using flammable cleaning and disinfecting agents. If flammable agents are used, allow them to evaporate before HF surgery.

**Sterilization:** The dispersive cables must be cleaned and disinfected prior to sterilization. When preparing the dispersive cable for sterilization, do not coil the cable into a circle of less than 10cm (4") in diameter. Pack the cable in disposable sterilization packaging (single or double packaging) made from paper/polyethylene and/or put it in a sterilization container. Packaging must comply with the ISO 11607-1 and EN 868-5 standards for packaging sterilized instruments. Sterilization of the high-frequency cables should be done by steam sterilization using the fractionated pre-vacuumed method. The cables must be sterilized at a minimum of 134°C [273°F] and a maximum of 137°C [278°F] in saturated steam during a holding (exposure) time of at least 4 minutes to at most 20 minutes, then dried in a vacuum for at least 20 minutes. When sterilizing several products at the same time in a steam autoclave, ensure that the maximum load capacity of the autoclave specified by the manufacturer is not exceeded. Sterilization must be performed in accordance with EN ISO 17665 (Sterilization of Medical Devices in Moist Heat).

The above sterilization process has been validated to assure a Sterility Assurance Level of  $10^{-6}$  ( $SAL \leq 1:10^6$ ) and the proof of the effectiveness of steam sterilization has been provided by an independent, accredited test laboratory (test report number AYA-SVR100620-01 dated 10.06.2020). The use of other sterilization methods (e.g. ethylene oxide gas, formaldehyde, gamma radiation, or low-temperature plasma sterilization) or other autoclave sterilization cycles is outside the manufacturer's responsibility.

**Note:** Before use, the dispersive cables must be cooled to room temperature. It is the sole responsibility of the user to safeguard and maintain the sterile condition of the cables/ adapters after the sterilization process.

**Warning:** Only black/dark grey colored dispersive cables can be autoclaved. Do not autoclave blue colored non-autoclavable dispersive cables. Such an attempt will destroy the cable.

**Warning:** Do not sterilize with hot air. Do not use "Flash" autoclave procedures.

**Compatibility:** In order to ensure compatibility, the dispersive cable should only be used with the equipment and disposable grounding plate for which they have been designed and is specified in the label on the product's primary packaging. For product models and compatible devices, please refer to the catalog at [metkomedical.com/PDF/06ElectrosurgeryCablesandAccessories.pdf](http://metkomedical.com/PDF/06ElectrosurgeryCablesandAccessories.pdf)

#### Environmental Requirements:

Operating Temperature : 5°C - 45°C (41°F - 113°F)

Storage Temperature : 0°C - 50°C (32°F - 122°F)

Relative Humidity : 20 - 80% (not condensing)

Store the cables in a clean, cool, and dry place. During the storage, the products should be protected from sunlight. The storage of the product in its original packaging until their first use is recommended. Used cables should be loosely coiled in 10 to 15 cm diameter (4 to 6 inches) and must be stored individually in a protective container or sterile packaging.

#### Safety:

Degree of protection from electric shocks: type CF

Degree of protection against the ingress of water: Keep Dry

The FMT reusable dispersive cables are sold NON-STERILE. Clean, disinfect or sterilize prior to initial and each subsequent use.

The FMT reusable dispersive cables are not made with natural rubber latex.

Refer to the separate leaflet "Description of Symbols" located within the product package for the explanation of symbols.

**Reporting of serious incidents:** Any serious incident related to the use of this product should be reported to both the manufacturer and the health authority/competent authority where the product is installed. Either contact your local representative or report to: [metko@metkomedical.com](mailto:metko@metkomedical.com)

**Caution:** Federal Law (U.S.A.) restricts this device to sale by or on the order of a physician.

**Warranty:** FMT reusable dispersive cables are under twelve (12) months warranty against material and workmanship defects from the date of the purchase. Damage caused by misuse, overuse, mishandling, or normal wear is not covered under warranty. Autoclavable reusable dispersive cables are designed to withstand a minimum of 50 sterilization cycles when adequately cared for and sterilized according to instructions. After the first use, the product has a service life of 50 sterilization cycles or 12 months.

FMT® is a registered trademark of Metko A.Ş.

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