Indications: Reusable EKG clamp electrodes are used to receive electrocardiographic signals of a patient and transmit the signal to an electrocardiograph via an EKG cable. EKG clamp electrodes are intended to be used by qualified medical personnel in hospitals and hospital-type facilities. EKG clamps can be used in pediatric and adult patients depending on the selection of the appropriate size.


Contraindications: No known contraindications to the use of the product.

Adverse Reactions: No allergic reactions are observed at places where EKG clamp electrodes are in contact with the intact skin of a patient. However, such reaction cannot be completely excluded in people particularly susceptible to allergies, it should be noted that no such events have been reported to the manufacturer. Incorrect cleaning or disinfection, or a failure to perform them, as well as failure to observe hygiene procedures may cause the transfer of pathogens and infections.

Clinical Benefits: Enables ECG recording with diagnostic EKG equipment.

Product Description: Reusable EKG clamp electrodes with universal connection are used in diagnostic EKG recordings. The EKG clamp electrode consists of a Silver/Silver Chloride (Ag/AgCl) plated metal electrode placed on one of the two plastic arms tensioned by a stainless steel spring. The screw on the electrode and the press stud behind the screw allows Ø4.0mm banana, Ø3.2mm straight pin, snap, or grabber all patient terminations to be attached to the electrode.

- CEA11** EKG clamp electrode, set of 4 pieces, IEC color code, Adult
- CEA12** EKG clamp electrode, set of 4 pieces, IEC color code, Pediatric
- CEAA1** EKG clamp electrode, set of 4 pieces, AHA/AAMI color code, Adult
- CEAA2** EKG clamp electrode, set of 4 pieces, AHA/AAMI color code, Pediatric

Directions for Use:

- 1- Before use, check that the product is intact and clean.
- 2- Clean the skin surfaces with a cotton cloth soaked with alcohol. If there is hair on the application area, shave it before the procedure.
- 3- Insert the Ø4.0mm banana plug (or Ø3.2mm straight pin) of patient termination of EKG cable to the hole of the screw. Then tighten the screw till the clamp electrode is interlinked with the patient termination. For snap patient termination push the snap to the press stud at the rear of the screw. For the grabber termination clip the grabber latch to the press stud.
- 4- Spread conductive gel on the cleaned skin areas with an amount of enough to ensure complete adherence to the metallic plate.
- 5- Position the upper limbs about 10cm above the wrists on the internal forearm face and the lower limbs about 10cm above the ankle on the internal leg face.
- 6- Verify that the ECG signal appears on the electrocardiograph screen. Please refer to the electrocardiograph and EKG cable instructions for use for detailed information on diagnostic EKG recording. After the procedure, remove the clamp electrodes from the limbs.

Warnings:

- 1- Before each use, visually inspect the EKG clamp electrodes. The EKG clamp electrodes should be free of visible defects and damage. Do not use the clamp electrode if the clamp spring has lost tension, Silver/Silver Chloride (Ag/AgCl) plating on the contact surface is worn, scratched, or oxidized, or if the patient terminals cannot be securely attached to the clamp electrode. Dispose of damaged EKG clamp electrodes according to hospital disposal practices, local environmental and waste disposal laws, and regulations.
- 2- Do not use EKG clamp electrodes during Magnetic Resonance Imaging (MRI) scanning or in an MRI environment as it may result in physical harm.
- 3- Keep leadwires connected to the electrodes as far as possible from electromagnetic sources such as the power transformer of an electromedical device. Also, avoid placing the EKG cable alongside the power cord of an electromedical device. Strong electromagnetic fields may adversely interfere with the ECG examination.
- 4- To prevent shock hazards avoid any contact of EKG electrodes with the external current.
- 5- To avoid a shock hazard and interface from nearby electrical equipment, keep leadwires connected to electrodes away from grounded metal and other electrical equipment.
- 6- Do not modify the EKG clamp electrodes in any way. Modifications may affect patient safety, performance, and accuracy and void the warranty.
- 7- Do not apply electroconductive gel and electrodes on injured or irritated skin.
- 8- Mishandling of EKG clamp electrodes may cause the latch spring to lose tension, scratch and damage the Ag/AgCl layer on the electrode surface, or deform the clamp electrode.
- 9- Do not use clamp electrodes for monitoring during defibrillation or cardioversion.
- 10- Do not use clamp electrodes for medium or long-term monitoring (monitoring during surgical operations, ECG recording, or stress tests).
- 11- For further information and warnings read the instruction manual of the electrocardiograph and EKG cable before using the EKG clamp electrodes.

Manual Cleaning and Disinfection:

Manual Cleaning: EKG clamp electrodes must be cleaned and disinfected immediately after each use. Remove heavy debris from the EKG clamp electrodes by wiping with an appropriate, moist lint-free cloth, sponge, or gauze pad. Prepare the detergent solution with reverse osmosis or distilled water according to the manufacturer's instructions. Saturate a synthetic felt or sponge pad with the cleaning solution. Wipe all surfaces of the EKG clamp electrodes until all visible contamination is removed. Wash the pad in the cleaning solution and repeat wiping the clamp. Use a plastic cleaning brush for adhered particles or dried secretions. Wipe or brush all surfaces at least six strokes back and forth. At least allow a contact time of 2 minutes. If visible contamination is still present on the electrode at the end of the cleaning repeat the cleaning process. Moist a clean cloth, sponge, or gauze pad with tap water or reverse osmosis and thoroughly wipe the EKG clamp

electrodes to remove any debris and all traces of the cleaning solution. Dry manually by using a lint-free cloth or absorbent paper and, in particular, for drying cavities, sterile compressed air.

Manual Disinfection: Prepare the disinfectant solution according to the manufacturer's instructions. Only disinfect clean products that have dried slightly. Saturate a synthetic felt or sponge pad with the disinfectant solution. Wipe all surfaces of the EKG clamp electrodes. Wash the pad in the disinfectant solution and repeat wiping the clamp. Wipe or brush all surfaces at least six strokes back and forth. At least allow a contact time of 6 minutes. Moist a clean cloth, sponge, or gauze pad with reverse osmosis or distilled water and thoroughly wipe the EKG clamp electrodes to remove any residual debris and all traces of disinfectant solution. Dry manually by using a lint-free cloth or absorbent paper and, in particular, for drying cavities, sterile compressed air.

EKG clamp electrodes, especially metallic parts should never be immersed in cleaning, disinfecting, or rinsing solutions, or placed in an ultrasonic cleaning device. Soaking will reduce the life of EKG clamp electrodes by oxidation of metal parts. Please do not clean with hard or frictional materials, abrasive cleaners, or solvents. Cleaning agents containing phenol or chloride are not suitable. The EKG clamp electrodes should be completely dry before use, do not use wet EKG clamp electrodes. These instructions were validated using the below detergents and disinfectants.

Cleaning agent (manual cleaning): Enzol® / Cidezyme®, enzymatic detergent, Johnson&Johnson (add 8ml enzymatic detergent per 1 liter distilled water, for cable with dried-on organic matter use 16 ml enzymatic detergent per 1 liter distilled water/ warm water)

Disinfectant (manual disinfection): Cidex OPA®, 0.55% ortho-phthalaldehyde solution, Johnson&Johnson (at a minimum temperature of 20°C/68°F).

Reverse osmosis/ distilled water.

Caution: During the cleaning and disinfection of EKG clamp electrodes, disposable gloves, protective eyewear, and a filtration mask should be used to reduce the risk of transmission of infectious agents by splashing or creating aerosols.

Warning: FMT EKG clamp electrodes are sold NON-STERILE. Clean and disinfect the clamp electrodes before first and every use.

Warning: Before cleaning or disinfecting the clamp electrodes, disconnect them from the EKG cable and the electrocardiograph.

Warning: Improper cleaning can cause damage to the Ag/AgCl layer, which could result in signal interference and inaccurate records of ECG.

Warning: EKG clamp electrodes should be cleaned of overburden prior to disinfection to improve the effectiveness.

Warning: Do not sterilize EKG clamp electrodes by autoclave, radiation, or steam.

Warning: Avoid contact with strong, aromatic, chlorinated, ketone, ether, or ester solvents. Do not use bleach on electrical contacts.

Environmental Requirements:

Operating Temperature : 5°C - 45°C (41°F - 113°F)

Storage Temperature : 0°C - 50°C (32°F - 122°F)

Relative Humidity : 20 - 80% (not condensing)

The product should be protected from sunlight during storage. It is recommended that it be stored in its original package until its first use.

Compatibility: EKG clamp electrodes are fully compatible with diagnostic EKG cables with Ø4.0mm banana plug, Ø3.2mm straight pin, and standard snap, and grabber patient terminations.

Safety:

Degree of protection from electric shocks: type BF

Degree of protection against the ingress of water: Keep dry

FMT EKG clamp electrodes are sold NON-STERILE.

FMT EKG clamp electrodes are not made with natural rubber latex and PVC.

For an explanation of symbols, refer to the separate leaflet "Description of Symbols" located within the product package.

Reporting of Serious Incidents: Any serious incident related to the use of this product should be reported to both the manufacturer and the health authority/competent authority where the product is installed. Either contact your local representative or report to: metko@metkomedical.com

Caution: Federal Law (U.S.A.) restricts this device to sale by or on the order of a physician.

Warranty: FMT reusable EKG clamp electrodes are under three (3) months warranty against material and workmanship defects from the date of original purchase. In the warranty period, METKO will be responsible for repairing or changing the clamp electrode free of charge if the defect is proven. This warranty does not extend to any product that has been subject to misuse, wrong disinfection, neglect, or accident, or that has been damaged by causes external to the product, or that has been used in violation of the operating instructions supplied with the product. The product life is 3 years of the production date.

FMT® is a registered trademark of Metko A.Ş.

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