Indications: Reusable EKG adaptors are used to connect diagnostic EKG cables with 4.0 mm banana or 3.2 mm straight pin patient termination to pre-gelled disposable ECG electrodes with press stud or disposable ECG tab electrodes in diagnostic ECG recordings. Reusable EKG adaptors are intended to be used by qualified medical personnel in hospitals and hospital-type facilities. Reusable EKG adaptors can be used in pediatric and adult patients depending on selecting the appropriate disposable ECG electrode.

Contraindications: No known contraindications to the use of the product.

Adverse Reactions: No allergic reactions are observed at places where EKG adaptors are in contact with a patient's intact skin. However, such reactions cannot be completely excluded in people particularly susceptible to allergies. It should be noted that no such events have been reported to the manufacturer. Incorrect cleaning or disinfection, or a failure to perform them, as well as failure to observe hygiene procedures, may cause the transfer of pathogens and infections.

Clinical Benefits: Enables ECG recording with diagnostic EKG equipment.

Product Description:

BS3 converts a 3.2mm straight pin to snap, for use with a disposable ECG electrode with press stud, gray, set of 10 pieces

BS4 converts 4.0mm banana to snap, for use with disposable ECG electrode with press stud, white, set of 10 pieces

BG3I/10, BG3A/10 converts a 3.2mm straight pin to a grabber, for use with a disposable ECG electrode with press stud, gray, set of 10 pieces

BG4/10 converts a 4.0mm banana to a grabber, for use with a disposable ECG electrode with press stud, white, set of 10 pieces

BA3 converts a 3.2mm banana to an alligator, for use with a disposable ECG tab electrode, set of 10 pieces

BA4 converts a 4.0mm banana to an alligator, for use with a disposable ECG tab electrode, set of 10 pieces

Directions for Use:

- 1- Before use, check that the product is intact and clean.
- 2- Clean the skin surfaces with a cotton cloth soaked in alcohol. If hair is on the application area, shave it before the procedure.
- 3- Insert the Ø4.0mm banana plug or Ø3.2mm straight pin of the patient termination of the EKG cable into the EKG adaptor. Attach the disposable ECG electrodes to the EKG adaptors.
- 4- Peel the back paper of the disposable ECG electrodes and place them in the required position according to your hospital's protocols. Verify that the ECG signal appears on the electrocardiograph screen.
- 5- After the procedure, remove disposable ECG electrodes from the patient and then from the EKG adaptors.

Warnings:

- 1- Before each use, visually inspect the EKG adaptors. They should be free of visible defects and damage. Do not use the EKG adaptor if the metal plating on the contact surface is worn, scratched, or oxidized. Dispose of damaged EKG adaptors according to hospital disposal practices and local environmental and waste disposal laws and regulations.
- 2- Do not use EKG adaptors during Magnetic Resonance Imaging (MRI) scanning or in an MRI environment, as this may result in physical harm.
- 3- Keep leadwires connected to the adaptors as far as possible from electromagnetic sources such as an electromedical device's power transformer. Also, avoid placing the EKG cable alongside the power cord of an electromedical device. Strong electromagnetic fields may adversely interfere with the ECG examination.
- 4- To prevent shock hazards, avoid any contact with EKG adaptors with the external current.
- 5- To avoid a shock hazard and interface from nearby electrical equipment, keep leadwires connected to EKG adaptors away from grounded metal and other electrical equipment.
- 6- Do not modify the EKG adaptors in any way. Modifications may affect patient safety, performance, and accuracy and void the warranty.
- 7- Do not apply disposable ECG electrodes and EKG adaptors on injured or irritated skin.
- 8- Do not use EKG adaptors for monitoring during defibrillation or cardioversion.
- 9- Do not use EKG adaptors for medium or long-term monitoring (monitoring during surgical operations, ECG recording or stress tests).
- 10- Do not pull back the EKG adaptor to remove the disposable ECG electrode from the skin. Gently lift up the edge of the disposable electrode and peel it off slowly. If required, moisten the edges of the disposable electrode using sterile water and dampened cotton swap.
- 11- For further information and warnings read the instruction manual of the electrocardiograph and EKG cable before using the EKG adaptors.

Manual Cleaning and Disinfection:

Manual Cleaning: EKG adaptors must be cleaned and disinfected immediately after each use. Remove heavy debris from the EKG adaptors by wiping with an appropriate, moist, lint-free cloth, sponge, or gauze pad. Prepare the detergent solution with reverse osmosis or distilled water according to the manufacturer's instructions. Saturate a synthetic felt or sponge pad with the cleaning solution. Wipe all surfaces of the EKG adaptors until all visible contamination is removed. Wash the pad in the cleaning solution and repeat wiping the adaptor. Use a plastic cleaning brush for adhered particles or dried secretions. Wipe or brush all surfaces at least six strokes back and forth. At least allow a contact time of 2 minutes. If visible contamination is still present on the adaptor at the end of the cleaning repeat the cleaning process. Moist a clean cloth, sponge, or gauze pad with tap water or reverse osmosis and thoroughly wipe the EKG adaptors to remove any debris and all traces of the cleaning solution. Dry manually by using a lint-free cloth or absorbent paper and, in particular, for drying cavities, sterile compressed air.

Manual Disinfection: Prepare the disinfectant solution according to the manufacturer's instructions. Only disinfect clean products that have dried slightly. Saturate a synthetic felt or sponge pad with the disinfectant solution. Wipe all surfaces of the EKG adaptors. Wash the pad in the disinfectant solution and repeat wiping the adaptor. Wipe or brush all surfaces at least six strokes back and forth. At least allow a contact time of 6 minutes. Moist a clean cloth, sponge, or gauze pad with reverse osmosis or distilled water and thoroughly wipe the EKG adaptors to remove any residual debris and all traces of disinfectant solution. Dry manually by using a lint-free cloth or absorbent paper and, in particular, for drying cavities, sterile compressed air.

EKG adaptors, especially metallic parts, should never be immersed in cleaning, disinfecting, or rinsing solutions or placed in an ultrasonic cleaning device. Soaking will reduce the life of EKG adaptors by oxidation of metal parts. Please do not clean with hard or frictional materials, abrasive cleaners, or solvents. Cleaning agents containing phenol or chloride are not suitable. The EKG adaptors should be completely dry before use, do not use wet EKG adaptors. These instructions were validated using the below detergents and disinfectants.

Cleaning agent (manual cleaning): Enzol® / Cidezime®, enzymatic detergent, Johnson&Johnson (add 8ml enzymatic detergent per 1 liter distilled water, for cable with dried-on organic matter use 16 ml enzymatic detergent per 1 liter distilled water/ warm water)
Disinfectant (manual disinfection): Cidex OPA®, 0.55% ortho-phthalaldehyde solution, Johnson&Johnson (at a minimum temperature of 20°C/68°F).
Reverse osmosis/ distilled water.

Caution: During the cleaning and disinfection of EKG adaptors, disposable gloves, protective eyewear, and a filtration mask should be used to reduce the risk of transmission of infectious agents by splashing or creating aerosols.

Warning: FMT EKG adaptors are sold NON-STERILE. Clean and disinfect the EKG adaptors before the first and every use.

Warning: Before cleaning or disinfecting the EKG adaptors, disconnect them from the EKG cable and the electrocardiograph.

Warning: EKG adaptors should be cleaned of overburden prior to disinfection to improve the effectiveness.

Warning: Do not sterilize EKG adaptors by autoclave, radiation or steam.

Warning: Avoid contact with strong, aromatic, chlorinated, ketone, ether, or ester solvents. Do not use bleach on electrical contacts.

Environmental Requirements:

Operating Temperature : 5°C - 45°C (41°F - 113°F)

Storage Temperature : 0°C - 50°C (32°F - 122°F)

Relative Humidity : 20 - 80% (not condensing)

During the storage, the product should be protected from sunlight. It is recommended to store the product in its original package until the first use.

Compatibility: EKG adaptors are fully compatible with diagnostic EKG cables with Ø4.0mm banana plug or Ø3.2mm straight pin patient terminations.

Safety:

Degree of protection from electric shocks: type CF

Degree of protection against the ingress of water: Keep dry

FMT EKG adaptors are sold NON-STERILE.

FMT EKG adaptors are not made with natural rubber latex and PVC.

For an explanation of symbols, refer to the separate leaflet "Description of Symbols" located within the product package.

Reporting of serious incidents: Any serious incident related to the use of this product should be reported to both the manufacturer and the health authority/competent authority where the product is installed. Either contact your local representative or report to: metko@metkomedical.com

Caution: Federal Law (U.S.A.) restricts this device to sale by or on the order of a physician.

Warranty: FMT reusable EKG adaptors are under twelve (12) months warranty against material and workmanship defects from the date of original purchase. In the warranty period, METKO will be responsible for repairing or changing the EKG adaptor free of charge if the defect is proven. This warranty does not extend to any product that has been subject to misuse, wrong disinfection, neglect, or accident, or that has been damaged by causes external to the product or that has been used in violation of the operating instructions supplied with the product. The product life is 3 years from the production date.

FMT® is a registered trademark of Metko A.Ş.

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