

FMT

Medical Sensors, Cables and Accessories



SPECIAL CABLES & ACCESSORIES



Special Cables & Accessories

TENS/EMS Cables & Electrodes



• DJO Global /Chattanooga Intellect Lead Wires

Four female 1.0mm (0,04") sockets to dual 2.0mm (0,08") pins, 2.5 m

TEC-01/1 - Channel 1, OEM P/N 27312

TEC-01/2 - Channel 2, OEM P/N 27313

TEC-01/3 - Channel 3, OEM P/N 27314

TEC-01/4 - Channel 4, OEM P/N 27315

Compatible with Chattanooga Intellect Advanced Combo/Stim, Intellect Mobile Combo/Stim, Intellect Legend XT, Vectra Genisys, TransPort.



• Chattanooga CPS & Fizyomed Fizyotens Lead Wires

3.5 mm (1/8") Stereo Plug to dual 2.0mm (0,08") pins, 2.5 m

TEC-02/1 - Channel 1

TEC-02/2 - Channel 2

TEC-02/3 - Channel 3

TEC-02/4 - Channel 4

Compatible with Chattanooga CPS 200/400 Combo/Stim, Fizyotens 4000/6000



• ATI Lead Wires

3 pole Plug to dual 2.0mm (0,08") pins, 2.5 m

TEC-03/1 - Channel 1

TEC-03/2 - Channel 2

TEC-03/3 - Channel 3

TEC-03/4 - Channel 4

Compatible with ATI AT 001/ 002, GT-001/001 Plus / 002, GB-001, AT-001, ST-001 200/400



• TENS Lead Wire with 4.0mm Right Angle Female Socket

4.0 mm Right Angle Female Socket to dual 2.0mm (0,08") pins, 1.2 m (47")

TEC-04

Compatible with EMS 5000, 6500/ ProCura/ Promed 300 TENS/ Accurate / Protens/ Rice AR1-32/ TENS 2500, 3000, 6000/ Elle TENS 2, Elle TENS Plus / Trimbio Dual Channel / StimTec Prime, StimTec Target / TwinStimPlus, Twinstim Digital EMS & TENS Combo/ TENS 3000 / Axion



• Vacuum Electrode Cup

Ø60mm, 4mm Fitting

VEC-01

• Vacuum Electrode Sponge

Ø60mm

VS-01



• Vacuum Electrode Hose

1.85m, 4 mm Fitting

TVH-01

Special Cables & Accessories

Flow Sensor Adapter Cables

AEC-SLE

Flow sensor cable for use with SLE5000 and SLE4000 ventilators (2.0m/ 6,56 ft.)

OEM P/N: N6656

Compatible with N5302 single use flow sensors.



AEC-GEEC

Neonatal flow sensor cable for use with GE Engström Carestation ventilators (2.0m/ 6,56 ft.)

OEM P/N: 1500-5604-000

Compatible with 1505-3272-000 reusable neonatal flow sensors.



AEC-FAB

Flow Sensor Cable for use with Fabian Acutronic Ventilators (2.0m/ 6,56 ft.)

OEM P/N: 1016

Compatible with 1031 reusable and 7058.1 single use flow sensors.



Special Cables & Accessories



AEC-ARG
Aerogen Pro/Pro-X Controller cable for
Aerogen Solo Nebulizers, 1.8 m,
OEM P/N: AG-AP1085



AEC-ARG/GE
Aerogen Nebuliser Cable to GE/Datex-Ohmeda
ENGSTROM CARESTATION S/5 Ventilators, 1.8 m
OEM P/N: AG-AP1085-GE, 1505-5602-000



SAC 01
GE AM4/AM5 coiled patient cable. (4.6m/15ft.)
OEM P/N: 700044-201, 700044202, 700044-204

SAC 01-S
GE AM4/AM5 coiled patient cable. (1.3m/50 in.)
OEM P/N: 700044-203, 700044205, 700044-210

GE MAX-1 Exercise Testing System, GE CASE 12, CASE 15, CASE 16, MAC PC, MAC6, MAC8, MAC12, CENTRA Systems, MAX-Personal, MAC VU System with AM4 or AM5 Acquisition Module



SAC 02
GE CAM14 coiled patient cable. (4.6m/15ft.)
OEM P/N: 2016560-003, 2016560-002

SAC 02-S
GE CAM14 coiled patient cable. (1.3m/50")
OEM P/N: 2016560-001

GE MAC 5000 Resting ECG System, MAC 5000ST Exercise Testing System, CASE, CASE 8000 Exercise Testing Systems with CAM14 Acquisition Module



FMT EEG/S
EEG Elektrode Cable (1.55 m / 61 inch)
Silver Plated 12 pcs / set

FMT EEG/G
EEG Elektrode Cable (1.55 m / 61 inch)
Golden Plated 12 pcs / set

1.5mm DIN type Connector (DIN 42 802)

10 mm (0.39") elektrode diameter



FMT is registered trademark of **Metko Ltd.**



Headquarter & Factory

İvedik O.S.B. Ağaç İşleri Sanayi Sit.
1354. Cad. 1358.Sok. No:9
06378 Yenimahalle

ANKARA / TÜRKİYE

Tel (pbx) : +90 312 387 12 46
Fax : +90 312 387 12 51

Branch

Atatürk Mah. Sedef Cad. 38.Ada
Ata Plaza 3-5 Blokları Kat:2 No:57
34758 Ataşehir

İSTANBUL / TÜRKİYE

Tel (pbx) : +90 216 455 91 21
Fax : +90 216 455 48 24

www.metkold.com

Other products, logos and company names mentioned herein may be trademarks of their respective owners and are not to be taken as an endorsement of affiliation with FMT®. The information in this catalogue has been carefully checked and it is believed to be accurate. Images in this catalogue are representative of items offered. Actual product may have slight visual differences. In the interest of continued product development, Metko reserves the right to make changes and improvements to this catalogue and the product it described any time, without notice or obligation.